EC urine diluent

Cat. No.	Product Name	Packaging (Content)
REG00059	EC urine diluent	1 × 100 mL, Instruction for use

EN

INTENDED USE

Urine diluent is a buffered reagent used to dilute urine samples for the *in vitro* quantitative determination of sodium (Na⁺), potassium (K⁺) and chlorides (Cl⁻) on an EC 90 analyzer. The results are used to evaluate the balance management of the organism with ions. For professional use in clinical laboratory only.

PRINCIPLE

Urine must be diluted to overcome potential matrix effects and to permit accurate measurements.

REAGENT DESCRIPTION AND COMPOSITION

Active ingredients	Concentration
HEPES	<0.3 %
Magnesium acetate	<0.2 %
Triethanolamine	<0.3 %
Preservatives	<0.01 %

REAGENT PREPARATION

Ready to use.

MATERIAL REQUIRED BUT IS NOT PROVIDED WITH THE DEVICE

EC 90 electrolyte analyser, Cat. No. INS00079

STABILITY AND STORAGE

Store at 2–30 °C until expiration date on labels. Protect from freezing. Avoid contaminating the reagent. Keep the container tightly closed and store the reagent in the original container.

PROCEDURE

For specimen collection and preparation only use suitable tubes or collection containers. Dilution of urine samples to urine diluent is in 1:3 ratio. Mix the diluent and sample thoroughly before measurement. 210 µl of diluted urine is required for the electrolyte measurement. For further information on sample procedures please refer to the EC 90 analyzers manual.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and / or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Supplemental information:

EUH 208 Contains reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one. May produce an allergic reaction.

WASTE MANAGEMENT

Please refer to local legal requirements.

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