EC cartridge plus iCa

Cat. No.:	Pack Name:	Number of Tests:
REG00066	EC cartridge XS plus iCa	250
REG00061	EC cartridge S plus iCa	500
REG00062	EC cartridge M plus iCa	1000
REG00063	EC cartridge L plus iCa	3000

EN

INTENDED USE

EC cartridge is an IVD kit intended for calibration and quantification of sodium (Na⁺), potassium (K⁺), ionized calcium (iCa²⁺) and chloride (Cl⁻) on the EC 90 analyzer series. Intended only for IVD use by a professionally trained person.

CONTENT

1) EC cartridge

Calibration solution 1:

Na⁺ 150 mmol/l, K⁺ 5 mmol/l, iCa²⁺ 1,25 mmol/l, Cl⁻ 103 mmol/l, buffer, preservative, wetting agent

Calibration solution 2:

Na⁺ 105 mmol/l, K⁺ 2 mmol/l, iCa²⁺ 2,34 mmol/l, Cl⁻ 71,1 mmol/l, buffer, preservative, wetting agent

Reference solution:

K⁺ 1,2 mol/l, buffer, preservative, wetting agent

2) EC Biosensor

PREPARATION

Insert EC Biosensor into the cartridge according to the user manual.

DIRECTIONS FOR USE

- In the analyzer software, in the maintenance screen, click on CHANGE EC CARTRIDGE button. Follow instructions on the screen.
- Unpack the new cartridge, remove the cap and seal cap from the top and place it in the analyzer.
- In the analyzer software click on OK button. The software will detect new cartridge and perform priming automatically.
- Calibration will be performed automatically.

STABILITY AND STORAGE

Store at 2–30 °C. Note the expiration date on the labels. Protect the cartridge from frost. On board stability: 90 days after opening.

PERFORMANCE DATA

Performance is included in the user manual.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Calibration solution 1 and 2:

EUH 208 Contains reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one. May produce an allergic reaction.

Reference solution: The reagent is not classified as dangerous.

WASTE MANAGEMENT

Users should be exceedingly careful when working with this biological material and treat the waste as a potential infectious material. Please refer to local legal requirements.

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