Erba Protime

Protrombinový čas

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LITERATURA

- CAEN J., LARRIEU M.-J., SAMAMA M. : "L'hémostase. Errichetti, A.M., Holden, A., Ansell, J.: Management of Oral Anticoagulant Therapy: Experience with an Anticoagulation Clinic. Arch Inter Med. 144:1966-68, 1984.
- Hirsh, J., Dalen, J.E., Deykin, D., Poller, L.: Oral Anticoagulants: Mechanisms of Action, Clinical Effectiveness, and Optimal Therapeutic Range. Chest 102(Suppl): 312S-316S, 1992.
- NCCLS: Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays; Approved Guideline. NCCLS document H21-A4. NCCLS, Wayne, PA, 2003.
- Palmer, R.N., Gralnick, H.R.: Inhibition of the Cold Activation of Factor VII and the Prothrombin Time. Am J Clin Path. 61: 618-622, 1984.
- Young, D.S., Thomas, D.W., Friedman, R.B., et al: Effect of Drugs on Clinical Laboratory Tests. Clin Chem 18:1041, 1972.
- Dalen, J.E., Hirsh, J.: American College of Chest Physicians and the National Heart, Lung, and Blood Institute National Conference on Antithrombotic Therapy. Arch Inter Med. 146:462-472, 1986.
- Palaretti, G., Coccheri, S., Poggi, M., et al: Oral Anticoagulant Therapy Control: Evidence that the INR Expression Improves the Interlaboratory Comparability of Results. The Bologna Oral Anticoagulant ControlExercise. Thromb Haemostasis 58:905-910, 1987.

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INTENDED USE

Erba Protime is intended for the determination of the prothrombin time (PT).

CLINICAL SIGNIFICANCE

The PT is used as a screening tool and as a quantitative test for coagulation factors in the extrinsic and common pathways. This test will be prolonged in patients with acquired or congenital disorders that reduce the activity of factors I (fibrinogen), II (Prothrombin), V, VII, and X. The PT is also widely used to monitor oral anticoagulant therapy.^{1,2} Oral anticoagulants reduce the activity of vitamin-K dependent clotting factors (II, VII, IX, X, Protein C, and Protein S), and the PT is prolonged as a result.

PRINCIPLE

The one-stage PT measures the clotting time of plasma after adding a source of tissue factor (thromboplastin) and calcium. The recalcification of plasma in the presence of tissue factor generates activated Factor Xa (F.Xa). F.Xa in turn activates Prothrombin to thrombin, which converts fibrinogen to an insoluble fibrin clot. The Patient PT is compared to a normal standard.

COMPOSITION

Erba Protime: prepared from fresh rabbit cerebral tissues (< 2 %), 0.013 % of sodium azide, 5 % of buffers, salts, and stabilizers.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only. These reagents are to be used by certified medical laboratory personnel only.
- Do not ingest.
- Wear gloves when handling all kit components.
- Only use clean or single use laboratory equipment to avoid contaminations.
- The eventual rest of reagents should be disposed of in accordance with the internal regulations and in compliance with local and national regulations relating to the safe handling of waste.
- Reagents containing sodium azide should be discarded with care to prevent the formation of explosive metallic azides. If waste materials are dumped into sinks, use copious quantities of water to flush plumbing thoroughly.

WARNING – POTENTIAL BIOHAZARDOUS MATERIAL

Some reagents provided in these kits contain materials of human and/or animal origin.

Whenever human plasma is required for the preparation of these reagents, the plasmas are tested for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with regulatory safety precautions in the manipulation of these biological materials as if they were infectious.

WORKING REAGENT

Reconstitute each vial of Erba Protime with exactly 4 ml (EHL00001/EHL00002) or 10 ml (EHL00027) of distilled/deionized water. Then swirl the vial gently and allow the reconstituted material to stand at room temperature (18–25 °C) for 15 minutes. Mix gently before each use. Once reconstituted, the Erba Protime reagent presents turbidity but does not sediment during the course of a work shift, it is however recommended to swirl the vial gently before each shift. Avoid reagent contamination.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

Once reconstituted, the reagent remains stable for:

- 8 hours at 37 °C (open vial).
- 24 hours at 15–25 °C (open vial).
- 7 days at 2–8 °C (closed vial).

On-board the instrument, for optimal performance, it is advised to only take out of the vial the necessary volume to perform the daily tests, and keep the vial resealed in the refrigerator.

Do not freeze!

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout.

Blood (9 parts) should be collected into 3.2 % sodium citrate anticoagulant (1 part). Avoid hemolysis and contamination by tissue fluids. Samples that have less than 90 % of the expected fill volume should be rejected. Centrifuge blood for 15 minutes at 1500 x g. Test within 2 hours if samples are held at 22–24 °C. If testing is not completed within 24 hours, plasma should be frozen at -20 °C for up to two weeks or -70 °C for up to 6 months.

PROCEDURE

Manual method

- Pre-warm the Erba Protime to 37 °C.
- Add 100 µL of the patient plasma or control plasma to the cuvette and prewar to 37 °C.
- Add 200 µL of Erba Protime. Start simultaneously a timer.

The time required for clot formation is the PT.

Automated method

Refer to the instrument's operator's manual.

INTERPRETATION OF RESULTS

Note the clotting time of the patient's plasma and that of the reference normal plasma.

The PT if expressed in percentage activity compared to a normal can be calibrated using the calibrators Erba-PT – INR MultiCal (Cat. No.: EHL00013).

The INR is calculated using the ratio of the patient PT to the mean of the normal reference range (TPN), according to the following mathematical relationship:

$$INR = (Patient\ PT / TPN)^{ISI}$$

For the INR calculation, use the ISI value insert included in the box.

Ensure that the values obtained for the controls are within the ranges stated in the Assay Value provided in the control box. If the control values are outside the stated ranges, check all components of the test system to ensure that all are functioning correctly, i.e., assay conditions, reagents, integrity of the plasmas being tested, etc. If necessary, repeat the test-run.

REFERENCES VALUES

These values are intended as a guideline only. Each laboratory should establish a Normal Reference Range (NRR) using instrumentation, blood collection methods, and testing techniques used in that laboratory. The NRR should be reestablished or at least verified when changing lot numbers of the same reagent.^{3,4} A new NRR should be established with any change in reagents, instrumentation, blood collection techniques, or anticoagulant.

The clotting time of abnormal plasmas will depend on the ISI of the reagent lot in use.

QUALITY CONTROL

Each laboratory should establish a quality control program.

To ensure adequate quality, controls plasma are recommended. It is suggested to use two levels of control, one close to the normal patient values (Erba Control N, Cat. No.: EHL00014 or Erba Control N Plus, Cat. No.: EHL00016) and the second representative the pathologic values (Erba Control P, Cat. No.: EHL00015 or Erba Control P Plus, Cat. No.: EHL00017).

LIMITATIONS

Make sure that the control(s) and the patient's plasma are treated in the same way.

Sample

Reject any suspect sample. Microclots will induce an important shortening of the clotting time. On the contrary, a complete coagulation will prolong the clotting times because of consumption of factors and fibrinogen.

Do not store plasmas between 2–8 °C, as they may undergo cold activation resulting in significant shortening of the PT.⁴

Anticoagulant

Sodium oxalate, EDTA, and heparin are not suitable anticoagulants.

Respect the correct anticoagulant (Trisodium citrate)/blood sample volume ratio of 1 : 9. Plasma samples with hematocrits outside the range of 20–55 % may be improperly anticoagulated and should be adjusted appropriately.

Interferences

Results from lipemic, hemolyzed or icteric specimens should be interpreted with caution. Studies show the on ECL analyzer, levels of 5 g/l (500 mg/dl) of Triglycerides do not influence results of PT when the Fibrinogen of the sample is within normal. For low fibrinogen samples, however, the coagulation time can be shortened.

The PT may be prolonged by substances such as oral contraceptives, corticosteroids, EDTA, asparaginase, clofibrate, erythromycin, ethanol, tetracycline, and anticoagulants such as heparin and warfarin.

The PT may be shortened by substances including antihistamines, butabarbital, caffeine, oral contraceptives, phenobarbital, and vitamin K.

Equipment

Check the stability of the instrument's temperatures displayed for reaction and preheater at 37 °C (± 0.5 °C) during active distributions.

Ensure that the Erba Protime reagent is not contaminated by any trace of plasma.

Use only clean and disposable lab ware.

Heparins

A study shows that addition of Heparin up to 2 U/ml in normal samples has no significant effect on coagulation time when tested by Erba Protime. For AVK treated patients displaying INR of 2.5, an increase of 0.5 unit INR for each 1 U/ml of heparin was noted.

PERFORMANCES

These performances have been obtained using an analyzer ECL. Results may vary if a different instrument or manual procedure is used.

Precision:

	Intra-assay precision (N = 20)		Inter-assay precision (N = 20)	
Mean (s)	12.9	46.4	13.1	36.6
CV (%)	0.63	2.78	2.27	2.84

REFERENCE

- CAEN J., LARRIEU M.-J., SAMAMA M.: "L'hémostase. Errichetti, A.M., Holden, A., Ansell, J.: Management of Oral Anticoagulant Therapy: Experience with an Anticoagulation Clinic. Arch Inter Med. 144:1966-68, 1984.
- Hirsh, J., Dalen, J.E., Deykin, D., Poller, L.: Oral Anticoagulants: Mechanisms of Action, Clinical Effectiveness, and Optimal Therapeutic Range. Chest 102(Suppl): 312S-316S, 1992.
- NCCLS: Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays; Approved Guideline. NCCLS document H21-A4. NCCLS, Wayne, PA, 2003.
- Palmer, R.N., Gralnick, H.R.: Inhibition of the Cold Activation of Factor VII and the Prothrombin Time. Am J Clin Path. 81: 618-622, 1984.
- Young, D.S., Thomas, D.W., Friedman, R.B., et al: Effect of Drugs on Clinical Laboratory Tests. Clin Chem 18:1041, 1972.
- Dalen, J.E., Hirsh, J.: American College of Chest Physicians and the National Heart, Lung, and Blood Institute National Conference on Antithrombotic Therapy. Arch Inter Med. 146:462-472, 1986.
- Palaueri, G., Coccheri, S., Poggi, M., et al: Oral Anticoagulant Therapy Control: Evidence that the INR Expression Improves the Interlaboratory Comparability of Results. The Bologna Oral Anticoagulant Control Exercise. Thromb Haemostasis 58:905-910, 1987.

USED SYMBOLS

REF	Catalogue Number	IVD	In vitro Diagnostics	I	See Instruction for Use
LOT	Lot Number	M	Manufacturer	CONT	Content
Exp	Expiry Date	T	Storage Temperature		

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erbamannheim.com, www.erbamannheim.com

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- CAEN J., LARRIEU M.-J., SAMAMA M. : "L'hémostase. Errichetti, A.M., Holden, A., Ansell, J.: Management of Oral Anticoagulant Therapy: Experience with an Anticoagulation Clinic. Arch Inter Med. 144:1966-68, 1984.
- Hirsh, J., Dalen, J.E., Deykin, D., Poller, L.: Oral Anticoagulants: Mechanisms of Action, Clinical Effectiveness, and Optimal Therapeutic Range. Chest 102(Suppl): 312S-316S, 1992.
- NCCLS: Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays; Approved Guideline. NCCLS document H21-A4. NCCLS, Wayne, PA, 2003.
- Palmer, R.N., Gralnick, H.R.: Inhibition of the Cold Activation of Factor VII and the Prothrombin Time. Am J Clin Path. 81: 618-622, 1984.
- Young, D.S., Thomas, D.W., Friedman, R.B., et al: Effect of Drugs on Clinical Laboratory Tests. Clin Chem 18:1041, 1972.
- Dalen, J.E., Hirsh, J.: American College of Chest Physicians and the National Heart, Lung, and Blood Institute National Conference on Antithrombotic Therapy. Arch Inter Med. 146:462-472, 1986.
- Palaereti, G., Coccheri, S., Poggi, M., et al: Oral Anticoagulant Therapy Control: Evidence that the INR Expression Improves the Interlaboratory Comparability of Results. The Bologna Oral Anticoagulant ControlExercise. Thromb Haemostasis 58:905-910, 1987.

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REFERENCIAS

- CAEN J., LARRIEU M.-J., SAMAMA M.: "L'hémostase. Errichetti, A.M., Holden, A., Ansell, J.: Management of Oral Anticoagulant Therapy: Experience with an Anticoagulation Clinic. Arch Inter Med. 144:1966-68, 1984.
- Hirsh, J., Dalen, J.E., Deykin, D., Poller, L.: Oral Anticoagulants: Mechanisms of Action, Clinical Effectiveness, and Optimal Therapeutic Range. Chest 102(Suppl): 312S-316S, 1992.
- NCCLS: Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays; Approved Guideline. NCCLS document H21-A4. NCCLS, Wayne, PA, 2003.
- Palmer, R.N., Gralnick, H.R.: Inhibition of the Cold Activation of Factor VII and the Prothrombin Time. Am J Clin Path. 81 618-622, 1984.
- Young, D.S., Thomas, D.W., Friedman, R.B., et al: Effect of Drugs on Clinical Laboratory Tests. Clin Chem 18:1041, 1972.
- Dalen, J.E., Hirsh, J.: American College of Chest Physicians and the National Heart, Lung, and Blood Institute National Conference on Antithrombotic Therapy. Arch Inter Med. 146:462-472, 1986.
- Palaereti, G., Coccheri, S., Poggi, M., et al: Oral Anticoagulant Therapy Control: Evidence that the INR Expression Improves the Interlaboratory Comparability of Results. The Bologna Oral Anticoagulant Control Exercise. Thromb Haemostasis 58:905-910, 1987.

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