Erba Actime

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LITERATURA

- Langdell R, Wagner R, Brinkhous K (1953) Effect of Antihemophilic Factor on One Stage Clotting Tests, J. Lab. Clin. Med, 41: 637
- Proctor R, Rapaport S (1961) The Partial Thromboplastin Time with Kaolin, Am. J. Clin. Path, 36: 212
- Brandt JT and Triplett DA (1981) Laboratory Monitoring of Heparin. Effect of Reagents and Instruments on the Activated Partial Thromboplastin Time, AJCP, 76: 530-537
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Young DS et al. Effects of Drugs on Clinical Laboratory Tests, 3rd ed., AACCC Press, Washington, D.C., 1990
- Kelton JG and Hirsh J (1980) Bleeding Associated with Antithrombotic Therapy, Semin Hematol, 17: 259-91
- Estes J and Poulin P (1974) Pharmacokinetics of Heparin, Thromb Diath Haemorr, 33: 26-35
- Solomon G, Hilgartner M, Kutt H (1972) Coagulation defects caused by diphenylhydantoin, Neurology, 22: 1165-71
- Ambrus JL, Schimert G, Lajos TZ, Ambrus CM, Mink IB, Lassman HB, Moore RH, Melzer J (1971) Effect of antifibrinolytic agents and estrogens on blood loss and blood coagulation factors during open heart surgery, J Med, 2: 65-81
- Crowell EB Jr, Clatanoff DV, Kiekhof W (1971) The effect of oral contraceptives on factor VIII levels, Lab. Clin. Med, 77: 551-7

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10020364
10020407
10020802
10020803

Erba Actime

Activated Partial Thromboplastin Time

| Cat. No.: | Pack name: | Packaging (Content): |
|-----------|----------------|----------------------|
| EHL00003 | Erba Actime 30 | 6 × 5 ml |
| EHL00004 | Erba Actime 50 | 10 × 5 ml |

INTENDED USE

Erba Actime is an *in vitro* diagnostic assay for quantitative determination of activated partial thromboplastin time (APTT) in plasma, using a manual or automated method. Intended for screening of intrinsic coagulation pathway, detection of deficiencies in Factors VIII, IX, XI, XII, and Prekallikrein and high molecular weight kininogen. The APTT is also commonly used to monitor heparin therapy since APTT prolongation is directly proportional to increasing amounts of heparin. For professional use in clinical laboratory only.

CLINICAL SIGNIFICANCE

The activated partial thromboplastin time (APTT) is useful as a screening tool for the intrinsic coagulation pathway. It is sensitive to deficiencies of all plasma clotting factors except Factor VII. However, it is mainly used to detect deficiencies in Factors VIII, IX, XI, XII, and Prekallikrein and high molecular weight kininogen.

The APTT is also commonly used to monitor heparin therapy since APTT prolongation is directly proportional to increasing amounts of heparin.³

The APTT assay is not recommended for monitoring oral anticoagulant therapy, for which is best monitored by Prothrombin Time (PT) test.

PRINCIPLE

The APTT test is performed by adding reagent containing a plasma activator and phospholipid to the test specimen. This mixture is incubated for 3 minutes at +37 °C for optimum activation. Calcium chloride is added and clot formation is timed. Clot detection can be by mechanical, manual (tilt-tube), or photo-optical measurement.

COMPOSITION

Erba Actime contains a near-colloidal particle activator (magnesium-aluminium-silicate) for optimum sensitivity to factor deficiencies and to heparin. The reagent also contains phospholipids with buffer and stabilisers.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only. These reagents are to be used by certified medical laboratory personnel only.
- Do not ingest.
- Wear gloves when handling all kit components.
- Only use clean or single use laboratory equipment to avoid contaminations.
- The eventual rest of reagents should be disposed of in accordance with the internal regulations and in compliance with local and national regulations relating to the safe handling of waste. Reagent of the kit are not classified like dangerous.

WARNING – POTENTIAL BIOHAZARDOUS MATERIAL

Some reagents provided in these kits contain materials of human and/or animal origin. Therefore, users of reagents of these types must exercise extreme care in full compliance with regulatory safety precautions in the manipulation of these biological materials as if they were infectious.

WORKING REAGENT

Reagent is ready to use.

Bring to room temperature and mix well by swirling or inversions prior to use.

Avoid reagent contamination.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

The reagent is stable for 30 days after opening.

Do not freeze.

REAGENTS REQUIRED BUT NOT PROVIDED

Erba Calcium Chloride – 0.025M CaCl₂ (Cat.No.: EHL00020)

SAMPLE COLLECTION AND PREPARATION

It is recommended that specimen collection and storage be carried out in accordance with the CLSI guideline H21-A5.⁴

Plastic or siliconised glass should be used throughout.

Blood (9 parts) should be collected into 3.2 % or 3.8 % sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes.

Plasma should be kept between 2–8 °C or 18–24 °C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20 °C for 2 weeks or -70 °C for 6 months. Thaw quickly at 37 °C prior to testing. Do not keep at 37 °C for more than 5 minutes. This will minimize the neutralization of the lupus inhibitor.⁴

Erroneous results may be caused by contamination with tissue fluids or stasis. Avoid agitation, air bubbles or foaming. For the effects of commonly administered drugs, refer to Young, et al.⁵

PROCEDURE

Manual method

- Pre-warm well mixed Erba Actime and Erba Calcium Chloride to 37 °C.
- Add 100 µl of the patient or control plasma into a reaction tube and incubate at 37 °C for 2 minutes.
- Add 100 µl of Erba Actime, mix and incubate at 37 °C for exactly 5 minutes.
- Add 100 µl of prewarmed Erba Calcium Chloride and mix. Start simultaneously a timer.

The time required for clot formation is the APTT (s).

If using the manual-tilt method, gently rock the reaction tube back and forth in the 37 °C water bath, continually looking for clot formation. Perform all tests in duplicate.

Automated method

Refer to the instrument's operator's manual.

REFERENCES VALUES

References values may vary depending on local conditions (type of population ...), the techniques and systems in use. Therefore, it is necessary that each laboratory establish its own normal ranges and acceptable control values for their particular local patient population. In general, values are considered normal if they fall within the range of mean ± 2 standard deviations.

QUALITY CONTROL

Each laboratory should establish a quality control program.

To ensure adequate quality, control plasmas are recommended. It is suggested to use two levels of control, one close to the normal patient values (Erba Control N, Cat. No.: EHL00014 or Erba Control N Plus, Cat. No.: EHL00016) and the second representative the pathologic values (Erba Control P, Cat. No.: EHL00015 or Erba Control P Plus, Cat. No.: EHL00017).

LIMITATIONS

Expected values for the APTT test will vary from one laboratory to another, depending on the technique used. The method of clot detection, temperature, pH, collection technique, type of anticoagulant and time and method of specimen storage are all very important. Plasma sample collection and storage conditions should be standardized and carefully controlled. Unexpected results should be confirmed by additional tests. Platelet fragments present in a specimen may cause the release of phospholipids, and thus the neutralization of any lupus inhibitor present in the specimen. The use of specimens with small plasma volumes should be avoided due to possible physiological pH changes. Testing could be affected by several drugs.⁵ An increase in the APTT results may be caused by the administration of diphenylhydantoin, heparin, warfarin and radiographic agents.^{5,8} Decreased APTT values may be seen during the use of oral contraceptives, or male estrogen therapy.^{9,10} Thus, laboratories should establish their own expected values for patients and well defined performance standards for the control.

ANALYTICAL PERFORMANCE

Results may vary if a different instrument or manual procedure is used. Each laboratory should establish its own performance data.

Precision

| Intra assay | ECL 105/412 | | ECL 760 | |
|----------------|--------------------|--------|--------------------|--------|
| | Clot formation (s) | CV (%) | Clot formation (s) | CV (%) |
| Erba Control N | 29.29 | 2.39 | 28.20 | 1.96 |
| Erba Control P | 54.61 | 2.43 | 53.61 | 3.36 |
| High Control | 97.26 | 2.59 | 80.03 | 1.77 |

| Inter assay | ECL 105/412 | | ECL 760 | |
|----------------|--------------------|--------|--------------------|--------|
| | Clot formation (s) | CV (%) | Clot formation (s) | CV (%) |
| Erba Control N | 28.57 | 3.56 | 29.23 | 4.55 |
| Erba Control P | 52.11 | 4.42 | 51.79 | 5.13 |

Factor sensitivity

| % Factor | Factor VIII (s) | Factor IX (s) | Factor XI (s) |
|----------|-----------------|---------------|---------------|
| <1 | 85.7 | 70.9 | 92.4 |
| 10 | 45.9 | 44.4 | 50.7 |
| 40 | 34.2 | 33.9 | 34.4 |
| 100 | 28.8 | 28.8 | 28.8 |

Heparin sensitivity

| Heparin (IU/mL) | Clot formation (s) |
|-----------------|--------------------|
| 0 | 29.9 |
| 0.2 | 70.0 |
| 0.4 | 174.0 |

INTERFERENCES

APPT is not affected by substances in concentrations up to:

| Interfering substance | ECL 105/412 | ECL 760 |
|-----------------------|----------------|----------------|
| Hemoglobin | min. 7 g/L | min. 2 g/L |
| Triglycerides | min. 535 mg/dL | min. 320 mg/dL |
| Bilirubin | min. 24 mg/dL | min. 24 mg/dL |

REFERENCES

- Langdell R, Wagner R, Brinkhous K (1953) Effect of Antihemophilic Factor on One Stage Clotting Tests, J. Lab. Clin. Med, 41: 637
- Proctor R, Rapaport S (1961) The Partial Thromboplastin Time with Kaolin, Am. J. Clin. Path, 36: 212
- Brandt JT and Triplett DA (1981) Laboratory Monitoring of Heparin. Effect of Reagents and Instruments on the Activated Partial Thromboplastin Time, AJCP, 76: 530-537
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Young DS et al. Effects of Drugs on Clinical Laboratory Tests, 3rd ed., AACC Press, Washington, D.C., 1990
- Kelton JG and Hirsh J (1980) Bleeding Associated with Antithrombotic Therapy, Semin Hematol, 17: 259-91
- Estes J and Poulin P (1974) Pharmacokinetics of Heparin, Thromb Diath Haemorr, 33: 28-35
- Solomon G, Hilgartner M, Kutt H (1972) Coagulation defects caused by diphenylhydantoin, Neurology, 22: 1165-71
- Ambrus JL, Schimert G, Lajos TZ, Ambrus CM, Mink IB, Lassman HB, Moore RH, Melzer J (1971) Effect of antifibrinolytic agents and estrogens on blood loss and blood coagulation factors during open heart surgery, J Med, 2: 65-81
- Crowell EB Jr, Clatanoff DV, Kiekhof W (1971) The effect of oral contraceptives on factor VIII levels, Lab. Clin. Med, 77: 551-7

USED SYMBOLS

| | | | | |
|-------------|------------------|---------------------|----------------------|-------------------------|
| LOT | Lot Number | IVD | In vitro Diagnostics | See Instruction for Use |
| REF | Catalogue Number | Manufacturer | CONT | Content |
| Expiry Date | | Storage Temperature | | |

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REFERENCIAS

- Langdell R, Wagner R, Brinkhous K (1953) Effect of Antihemophilic Factor on One Stage Clotting Tests, J. Lab. Clin. Med, 41: 637
- Proctor R, Rapaport S (1961) The Partial Thromboplastin Time with Kaolin, Am. J. Clin. Path, 36: 212
- Brandt JT and Triplett DA (1981) Laboratory Monitoring of Heparin. Effect of Reagents and Instruments on the Activated Partial Thromboplastin Time, AJCP, 76: 530-537
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Young DS et al. Effects of Drugs on Clinical Laboratory Tests, 3rd ed., AACC Press, Washington, D.C., 1990
- Kelton JG and Hirsh J (1980) Bleeding Associated with Antithrombotic Therapy, Semin Hematol, 17: 259-91
- Estes J and Poulin P (1974) Pharmacokinetics of Heparin, Thromb Diath Haemorr, 33: 26-35
- Solomon G, Hilgartner M, Kutt H (1972) Coagulation defects caused by diphenylhydantoin, Neurology, 22: 1165-71
- Ambrus JL, Schimert G, Lajos TZ, Ambrus CM, Mink IB, Lassman HB, Moore RH, Melzer J (1971) Effect of antifibrinolytic agents and estrogens on blood loss and blood coagulation factors during open heart surgery, J Med, 2: 65-81
- Crowell EB Jr, Clatanoff DV, Kiekhof W (1971) The effect of oral contraceptives on factor VIII levels, Lab. Clin. Med, 77: 551-7

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- Langdell R, Wagner R, Brinkhous K (1953) Effect of Antihemophilic Factor on One Stage Clotting Tests, J. Lab. Clin. Med, 41: 637
- Proctor R, Rapaport S (1961) The Partial Thromboplastin Time with Kaolin, Am. J. Clin. Path, 36: 212
- Brandt JT and Triplett DA (1981) Laboratory Monitoring of Heparin. Effect of Reagents and Instruments on the Activated Partial Thromboplastin Time, AJCP, 76: 530-537
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Young DS et al. Effects of Drugs on Clinical Laboratory Tests, 3rd ed., AACC Press, Washington, D.C., 1990
- Kelton JG and Hirsh J (1980) Bleeding Associated with Antithrombotic Therapy, Semin Hematol, 17: 259-91
- Estes J and Poulin P (1974) Pharmacokinetics of Heparin, Thromb Diath Haemorr, 33: 28-35
- Solomon G, Hilgartner M, Kutt H (1972) Coagulation defects caused by diphenylhydanto-in, Neurology, 22: 1165-71
- Ambrus JL, Schimert G, Lajos TZ, Ambrus CM, Mink IB, Lassman HB, Moore RH, Melzer J (1971) Effect of antifibrinolytic agents and estrogens on blood loss and red blood coagulation factors during open heart surgery, J Med, 2: 65-81
- Crowell EB Jr, Clatanoff DV, Kiehofer W (1971) The effect of oral contraceptives on factor VIII levels, Lab. Clin. Med, 77: 551-7

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