Erba Thrombin Reagent

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Erba Thrombin Reagent

For Fibrinogen Determination

Cat. No.	Pack name	Packaging (Content)
EHL00005	Erba Thrombin Reagent 10	5 × 2 mL

EN

INTENDED USE

Erba Thrombin Reagent is intended for use in the quantitative determination of fibrinogen in human plasma using the Clauss Method.¹

CLINICAL SIGNIFICANCE

Fibrinogen (also called Factor I) is a glycoprotein synthesised in the liver hepatocytes and megakaryocytes. Fibrinogen is composed of two sets of three different chains (α , β , and γ), linked to each other by disulfide bonds. On the alpha and beta chains, there is small peptide sequence (called a fibrinopeptide). Thrombin (Factor IIa) is responsible for converting fibrinogen into fibrin, by splitting out fibrinopeptides monomers. The monomers aggregate to form Fibrin, which is then cross linked by Factor XIII a to form an insoluble fibrin clot.

Fibrinogen is an acute phase reactant protein. Therefore its blood level can be increased as a response of inflammatory states, with infections, during pregnancy, and after trauma.² A decrease of the fibrinogen level is found in congenital hypofibrinogenemia and afibrinogenemia. It can also be decreased in cases of Disseminated Intra-vascular Coagulation (DIC), systemic fibrinolysis, pancreatic and severe hepatic dysfunction.

PRINCIPLE

Clauss¹ developed a simple method for the quantitative determination of fibrinogen by measuring the clotting time of dilute plasma after the addition of Thrombin (> 30 NIH units/ml). This clot time is proportional to the fibrinogen concentration.

COMPOSITION

Erba Thrombin Reagent; approximately 100 NIH units/ml bovine thrombin with stabilisers.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only. These reagents are to be used by certified medical laboratory personnel only.
- Do not ingest.
- Wear gloves when handling all kit components.
- Only use clean or single use laboratory equipment to avoid contaminations.
- The eventual rest of reagents should be disposed of in accordance with the internal regulations and in compliance with local and national regulations relating to the safe handling of waste.

Reagent contains <30 % bovine thrombin.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Warning

Hazard statement:

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

Precautionary statement:

P302 + P352 IF ON SKIN: Wash with plenty of water.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

WARNING – POTENTIAL BIOHAZARDOUS MATERIAL

The Erba Thrombin Reagent contains bovine thrombin. Plasma products have been screened and founded negative (unless otherwise stated on the box or vial) for the presence of Hepatitis B Antigen (HBsAg), HIV 1 and 2

antibody and HCV antibody. Therefore, users of reagents of these types must exercise extreme care in full compliance with regulatory safety precautions in the manipulation of these biological materials as if they were infectious.

WORKING REAGENT

Reconstitute each vial of Erba Thrombin Reagent with exactly 2 mL of purified water. Take care when pipetting to avoid contamination. Swirl gently and allow to stand for 15 minutes. Mix well immediately before use. Do not shake.

Avoid reagent contamination

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

The reagent in the original closed vial after reconstitution is stable:

- 8 hours at 15–30 °C
- 1 week at 2–8 °C
- 1 month at -20 °C.

REAGENTS REQUIRED BUT NOT PROVIDE

Erba Standard Plasma (Cat. No.: EHL00012)

Erba Owrens Veronal Buffer (Cat. No.: EHL00021)

SAMPLE COLLECTION AND PREPARATION

It is recommended that specimen collection and storage be carried out in accordance with the CLSI guideline H21-A5.³ Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2 % sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept between 2–8 °C or between 18–24 °C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20 °C for 2 weeks or -70 °C for 6 months. Thaw quickly at 37 °C prior to testing. Do not keep at 37 °C for more than 5 minutes.³

PROCEDURE

Manual method

1. Standard Curve Preparation

A new standard curve should be used when reagent lot numbers change or if expected values in quality control testing change. Prepare the following dilutions of Erba Standard Plasma (Cat.No.: EHL00012) in Erba Owrens Veronal Buffer (Cat. No.: EHL00021)

Tube	ECL 760 Recommended dilution	Dilution (Manual, ECL 105, ECL 412)	Erba Standard Plasma (mL)	Owren's buffer (mL)
1	2/1	1 + 4	0.2	0.8
2	1/1	1 + 9	0.1	0.9
3	1/2	1 + 19	0.1	1.9
4	1/3	1 + 29	0.1	2.9
5	1/4	1 + 39	0.1	3.9

2. Patient Sample Preparation:

- Prepare 1 + 9 dilutions of the patient plasma or control plasma in Owrens Buffer.
- Mix without shaking.

3. Testing

- Pipette 0.2 mL of standard, patient or control dilution into a reaction tube and incubate at 37 °C for 2 minutes.
- Add 0.1 mL of Erba Thrombin Reagent (15–30 °C) and start simultaneously a timer.
- Note the time for clot formation to the nearest 0.1 second.
- Plot means Standard Clot Time (Y-axis) versus Fibrinogen Level (X-axis) on the fibrinogen graph paper. A straight line should be obtained. Assign the Erba Standard Plasma (Cat. No.: EHL00012) reference value to the 1 + 9 dilution to allow direct interpolation of patient and control values from the standard curve.

utomated method

Refer to the instrument's operator's manual.

REFERENCES VALUES

Reference limits^{4,5} for Fibrinogen (g/L) in plasma:

0–1 year: 1.50–3.40

1–6 years: 1.70–4.00

6–11 years: 1.55–4.00

11–16 years: 1.55– 4.50

16–18 years: 1.60– 4.20

≥18 years: 1.80– 4.20

Reference intervals vary from laboratory to laboratory depending on the population, the technique and reagent lot. Therefore, each laboratory must establish its own reference intervals or verify them whenever one or more of the aforementioned variables are changed. For more information on establishing reference intervals see CLSI document C28-A3.⁶

QUALITY CONTROL

Each laboratory should establish a quality control program.

To ensure adequate quality, control plasmas are recommended. It is suggested to use two levels of control, one close to the normal patient values (Erba Control N, Cat. No.: EHL00014 or Erba Control N Plus, Cat. No.: EHL00016) and the second representative the pathologic values (Erba Control P, Cat. No.: EHL00015 or Erba Control P Plus, Cat. No.: EHL00017).

LIMITATIONS

Heparin levels >0.6 units/ml and fibrinolytic degradation product >100 mg/ml may cause falsely low fibrinogen quantitation. If values fall outside the standard curve values for the patient samples, re-assay using an appropriate dilution to bring values into the standard range.

PERFORMANCES

These performances have been obtained using an analyzers ECL. Results may vary if a different instrument or manual procedure is used.

Precision

Intra assay	ECL 105/412		ECL 760	
	Clot formation (g/L)	CV (%)	Clot formation (g/L)	CV (%)
Erba Control N	2.59	6.68	3.30	3.38
Erba Control P	2.66	5.05	2.93	5.43
High Control	0.91	7.03	0.99	5.18

Inter assay	ECL 105/412		ECL 760	
	Clot formation (g/L)	CV (%)	Clot formation (g/L)	CV (%)
Erba Control N	2.46	7.31	3.51	7.58
Erba Control P	2.37	5.86	3.11	9.67

Limit of quantification	
ECL 105/412	ECL 760
0.368 g/L	0.517 g/L

Measuring range	
ECL 105	0.80–7.50 (g/L)
ECL 412	0.80–6.81 (g/L)
ECL 760	1.39–5.56 (g/L)

Interfering substance	ECL 105/412	ECL 760
Hemoglobin	min. 1 g/L	min. 3 g/L
Triglycerides	min. 320 mg/dL	min. 680 mg/dL
Bilirubin	min. 24 mg/dL	min. 28 mg/dL

REFERENCES

- Clauss A (1957) Gerinnungsphysiologische Schnell-methode zur bestimmung des Fibrinogens, Acta Haematol, 17:237-246
- Shaw TS (1977) Assays for Fibrinogen and its Derivatives CRC, Crit.Rev.Clin.Lab.Sci, 8:145-192
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI:H21-A5
- Hematology of Infancy and Childhood, Nathan and Orkin, Saunders 2003, P. Monagle: Developmental haemostasis: Impact for clinical haemostasis laboratories, Thromb Haemost 2006; 95: 362-72
- Czech haematological Society CzMA JEP
- Clinical and Laboratory Standards Institute. How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline. CLSI document C28-A3

USED SYMBOLS

LOT	Lot Number	IVD	In vitro diagnostic medical device	i	See Instructions for Use
REF	Catalogue Number	Manufacturer	Manufacturer	CONT	Content
Expiry	Expiry Date	Storage	Storage Temperature	Biological	Biological Risks

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erbamannheim.com, www.erbamannheim.com

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REFERENCIAS

- Clauss A (1957) Gerinnungsphysiologische Schnell-methode zur-bestimmung des Fibrinogens, Acta Haematol, 17:237-246
- Shaw TS (1977) Assays for Fibrinogen and its Derivatives CRC, Crit.Rev.Clin.Lab.Sci, 8:145-192
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI:H21-A5
- Hematology of Infancy and Childhood, Nathan and Orkin, Saunders 2003, P. Monagle: Developmental haemostasis: Impact for clinical haemostasis laboratories, Thromb Haemost 2006; 95: 362-72
- Czech haematological Society CzMA JEP
- Clinical and Laboratory Standards Institute. How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline. CLSI document C28-A3

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