# Erba Thrombin Time

## Thrombin Time

Cat. No.	Pack name	Packaging (Content)
EHL00007	Erba Thrombin Time	10 x 2 ml
EHL00026	Erba Thrombin Time	10 x 5 ml

EN

### INTENDED USE

Erba Thrombin Time is intended for the determination of the thrombin time.

### CLINICAL SIGNIFICANCE

The Erba Thrombin Time is commonly applied to detect various sources of interference with normal blood coagulation.

Prolongation of the thrombin clotting time can be taken as a qualitative indication of abnormal fibrinogen levels (high or low), or the presence of interfering substances such as Fibrin Degradation Products (FDPs) or heparin.

Quantitative evaluation of the possible causes of prolonged thrombin clotting time should be performed as follow-up studies, such as aPTT or chromogenic assay for heparin, Clauss fibrinogen, FDP determinations, heparin neutralisation by protamine sulphate<sup>2</sup> or polybrene<sup>3</sup> or reptilase assay<sup>4</sup> to distinguish between hypofibrinogenemia and FDP effects.

### PRINCIPLE

Clotting test.

### COMPOSITION

Erba Thrombin Time: lyophilised preparation of bovine thrombin with buffers and stabilisers. The reconstituted reagent contains approximately 10 NIH units/ml of thrombin. The reagent should be a white lyophilised plug.

### WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only. These reagents are to be used by certified medical laboratory personnel only.
  - Do not ingest.
  - Wear gloves when handling all kit components.
  - Only use clean or single use laboratory equipment to avoid contaminations.
  - The eventual rest of reagents should be disposed of in accordance with the internal regulations and in compliance with local and national regulations relating to the safe handling of waste.
- Reagent contains < 1 % thrombin.

### Hazards identification in accordance with Regulation (EC) No 1272/2008

#### Supplemental information:

EUH208 Contains Thrombin. May produce an allergic reaction.

EUH210 Safety data sheet available on request.

### WARNING – POTENTIAL BIOHAZARDOUS MATERIAL

The Erba Thrombin Time reagent contains bovine thrombin. Bovine thrombin contains bovine material. All donor animals were sourced from BSE-free herds. The cattle received ante- and post mortem health inspection by a veterinarian, and they were apparently free from infectious and contagious material. However, the material should be treated as potentially infectious.

### WORKING REAGENT

Reconstitute each vial with exactly 2 ml (EHL00007) or 5 ml (EHL00026) of purified water. Allow to stand for 5 minutes then mix gently by inversion and transfer to a plastic tube. Reconstituted reagent should be a clear, colourless solution. Avoid reagent contamination.

### STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

Once reconstituted, the reagent is stable for:

- 14 days at 2–8 °C
- 1 month at -20 °C

### SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout.

Blood (9 parts) should be collected into 3.2 % or 3.8 % sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at 15–25 °C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20 °C for 2 weeks or -70 °C for 6 months. Thaw quickly at 37 °C prior to testing. Do not keep at 37 °C for more than 5 minutes.<sup>4</sup>

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### LITERATURA

- Laposata et al. The clinical haemostasis handbook, yearbook medical publisher Inc., p219, 1989
- Thompson AR and Harker LA. Manual of haemostasis and thrombosis. 3rd Ed., F.A. Davis Co., p62, 1983
- DeMott WR. Laboratory Test Handbook, 2nd Ed., Jacobs D.S. et al Eds., Lexi-Comp Inc., p432–433, 1990
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for testing plasma based coagulation assays and molecular haemostasis assays: approved guideline, 5th edn. CLSI: H21-A5
- Gastineau DA et al. (1991) inhibitor of the thrombin time in systemic amyloidosis: a common coagulation abnormality. Blood, 77: 2637–40

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- Laposata et al. The clinical haemostasis handbook, yearbook medical publisher Inc., p219, 1989
- Thompson AR and Harker LA. Manual of haemostasis and thrombosis. 3rd Ed., F.A.Davis Co., p62, 1983
- DeMott WR. Laboratory Test Handbook, 2nd Ed., Jacobs D.S. et al Eds., Lexi-Comp Inc., p432-433, 1990
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for testing plasma based coagulation assays and molecular haemostasis assays: approved guideline, 5th edn. CLSI: H21-A5
- Gastineau DA et al. (1991) inhibitor of the thrombin time in systemic amyloidosis: a common coagulation abnormality, Blood, 77: 2637-40

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