# Erba Standard Plasma

## Haemostasis Calibration Plasma

Cat. No.:	Pack name:	Packaging (Content):
EHL00012	Erba Standard Plasma	5 × 1 ml

EN

### INTENDED USE

Erba Standard Plasma may be used as reference plasma for the following tests:

- Prothrombin Time (PT)
- Fibrinogen<sup>3</sup>
- Anti-Thrombin (AT III)
- Protein C
- Factors

### COMPOSITION

Erba Standard Plasma: prepared from a frozen pool of citrated plasma from healthy donors and is buffered and lyophilised to ensure stability of all plasma constituents.<sup>5</sup>

### WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only. These reagents are to be used by certified medical laboratory personnel only.
- Do not ingest.
- Wear gloves when handling all kit components.
- Only use clean or single use laboratory equipment to avoid contaminations.
- The eventual rest of reagents should be disposed of in accordance with the internal regulations and in compliance with local and national regulations relating to the safe handling of waste.

### WARNING – POTENTIAL BIOHAZARDOUS MATERIAL

Some reagents provided in these kits contain materials of human and/or animal origin. Whenever human plasma is required for the preparation of these reagents, the plasmas are tested for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with regulatory safety precautions in the manipulation of these biological materials as if they were infectious.

### WORKING REAGENT

Reconstitute each vial of Erba Standard Plasma with exactly 1 ml of distilled water. Swirl gently. Allow to stand for 20 minutes for complete dissolution. Avoid calibrator contamination

### STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. Once reconstituted, the control remains stable: • 4 hours at 2–8 °C.

### PROCEDURE

The calibrator should be treated in the same manner as the patients' sample in accordance with the instructions outlined in each particular assay procedure.

### EXPECTED VALUES

The calibrator values for each of the parameters may vary from one lot to another. Refer to the Assay Value insert provided in the box.

### QUALITY CONTROL

To ensure adequate quality, control plasmas are recommended. It is suggested to use two levels of control, one close to the normal patient values (Erba Control N, Cat. No.: EHL00014, Erba Control N Plus, Cat. No.: EHL00016) and the second representative the pathologic values (Erba Control P, Cat. No.: EHL00015, Erba Control P Plus, Cat. No.: EHL00017).

### LIMITATIONS

The results obtained with Erba Standard Plasma depend on several factors strongly associated with instrumentation, types of reagents, deficient substrates and laboratory to laboratory variations.<sup>7,8,9</sup> Each laboratory should establish an expected range for the particular instrument-reagent system.

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### LITERATURA

1. Babson AL and Flanagan ML (1975) Quantitative One Stage Assays for Factors V and X, AJCP, 64: 817-819.
2. Hardisty RM et al. (1962) A One Stage Factor VIII Assay and Its Use on Venous and Capillary Plasma. Thrombosis et Diathesis Haemorrhagica, 7:215-229.
3. Morse EE et al. (1971) Automated Fibrinogen Determination, AJCP, 55:671-676.
4. Elodi S et al (1978) Some Sources of Error in the One-Stage Assay of Factor VIII, Haemostasis, 7:1-9.
5. Thelin M (1968) Preparation and Standardization of a Stable AHF Plasma, Thrombosis et Diathesis Haemorrhagica, 19:423.
6. Kirkwood TBL et al. (1977) Identification of Sources of Variation in Factor VIII Assay, British Journal of Haematology, 37:559-568.
7. Goldenfarb MD (1971) Reproducibility in Coagulation Assays, AJCP, 55:561-564.
8. Palkuti HA and Longberry JR (1973) A Precision Study of Coagulation Factor Assay Techniques, AJCP, 59:231-235

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ  
e-mail: diagnostics@erbamannheim.com, www.erbamannheim.com

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- ЛІТЕРАТУРА
1. Babson AL and Flanagan ML (1975) Quantitative One Stage Assays for Factors V and X, AJCP, 64: 817-819.
2. Hardisty RM et al. (1962) A One Stage Factor VIII Assay and Its Use on Venous and Capillary Plasma. Thrombosis et Diathesis Haemorrhagica, 7: 215-229.
3. Morse EE et al. (1971) Automated Fibrinogen Determination, AJCP, 55: 671-676.
4. Elodi S et al (1978) Some Sources of Error in the One-Stage Assay of Factor VIII, Haemostasis, 7: 1-9.
5. Thelin M (1968) Preparation and Standardization of a Stable AHF Plasma, Thrombosis et Diathesis Haemorrhagica, 19: 423.
6. Kirkwood TBL et al. (1977) Identification of Sources of Variation in Factor VIII Assay, British Journal of Haematology, 37: 559-568.
7. Goldenfarb MD (1971) Reproducibility in Coagulation Assays, AJCP, 55: 6 61-664.
8. Palkuti HA and Longberry JR (1973) A Precision Study of Coagulation Factor Assay Techniques, AJCP, 59: 231-235

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