# Erba PT-INR Multical PT % / INR Calibrator

Cat. No.:	Pack name:	Packaging (Content):
EHL00013	Erba PT-INR MultiCal	4 x 1 mL (4 levels)

## INTENDED USE

Erba PT-INR MultiCal is a set of citrated plasma intended for establishing calibration curves required in testing for prothrombin time (PT) and for aiding in the verification of the ISI value of the Erba Protine and Protine LS (Cat. No.: EHL00001, EHL00002, EHL00024, EHL00046, EHL00047) used.

## COMPOSITION

Calibrators 1, 2, 3 and 4: human citrated plasma, lyophilised. Calibrator 1 is a lyophilised normal pool of plasma produced from the selected citrated plasma, of at least 100 blood donors. Calibrators 2, 3 and 4 are produced from the plasma of donors on long term oral anticoagulant therapy. As a result the Prothrombin complex factors (II, VII, IX, X, Proteins C and S) are present at the level normally expected in plasma from patients undergoing the same treatments, together with the proteins induced by Vitamin K absence (PIVKA inhibitors). No preservatives are added to the lyophilised plasmas.

## WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only. These reagents are to be used by certified medical laboratory personnel only.
- Do not ingest.
- Wear gloves when handling all kit components.
- Only use clean or single use laboratory equipment to avoid contaminations.
- The eventual rest of reagents should be disposed of in accordance with the internal regulations and in compliance with local and national regulations relating to the safe handling of waste.

## Hazards Identification in accordance with Regulation (EC) No 1272/2008

Reagent is not classified as dangerous.

## WARNING – POTENTIAL BIOHAZARDOUS MATERIAL

Some reagents provided in these kits contain materials of human and/or animal origin. Whenever human plasma is required for the preparation of these reagents, the plasmas are tested for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with regulatory safety precautions in the manipulation of these biological materials as if they were infectious.

## WORKING REAGENT

Reconstitute each vial with exactly 1 mL of distilled water. Allow the reconstituted material to stand at room temperature (between 18 and 25 °C) for 15 minutes. Then, swirl the vial gently before use. Avoid calibrators contamination.

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

Once reconstituted, they remain stable for:

- 4 hours at 15–20 °C
- 6 hours at 15 °C or on board ECL analyzer
- 8 hours at 2–8 °C
- 1 month frozen at -20 °C

Do not freeze / thaw more than once.

## PROCEDURE

Each calibrator should be treated in the same manner as the patients' sample in accordance with the instructions outlined in each particular test protocol. Refer to the instrument's operator's manual / application manual / reagents' technical insert.

## For the ISI verification:

Run the Calibrators 1, 2, 3 and 4, and write down the clotting times. Using a spreadsheet, plot a graph using for the X axis the Log (sec of the 4 calibrators), and for the Y axis the Log (INR). Add a regression line and select to visualize the Equation. The equation should be:  $a \times b$  format Where  $a = 1 / \text{ISI}$ . To verify the correct ISI, calculate  $\text{ISI} = 1 / a$ .

## EXPECTED VALUES

The prothrombin levels in Erba PT-INR MultiCal may vary from one lot to another. Refer to the Assay Value insert provided in the box.

## QUALITY CONTROL

To ensure adequate quality, controls plasma are recommended. It is suggested to use two levels of control, one close to the normal patient values (Erba Control N, Cat. No.: EHL00014 or Erba Control N Plus, Cat. No.: EHL00016) and the second representative the pathologic values (Erba Control P, Cat. No.: EHL00015 or Erba Control P Plus, Cat. No.: EHL00017).

## LIMITATIONS

A new calibration is required for each lot of reagents and for each instrument.

The results obtained with Erba PT-INR MultiCal depend on several factors strongly associated with instrumentation, types of reagents and laboratory to laboratory variations.

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