# Erba DDimer Control N + P Low & High D-Dimer Controls

Cat. No.:	Pack name:	Packaging (Content):
EHL00019	Erba DDimer Control N + P	R1: 5 × 1 ml Control N R2: 5 × 1 ml Control P

## INTENDED USE

Erba DDimer Control N + P may be used as a low and a high controls respectively, for the D-Dimer tests.

## COMPOSITION

Erba DDimer Control N + P: lyophilised D-Dimer enriched human plasma.

## WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only. These reagents are to be used by certified medical laboratory personnel only.
- Do not ingest.
- Wear gloves when handling all kit components.
- Only use clean or single use laboratory equipment to avoid contaminations.
- The eventual rest of reagents should be disposed of in accordance with the internal regulations and in compliance with local and national regulations relating to the safe handling of waste.

## WARNING – POTENTIAL BIOHAZARDOUS MATERIAL

Some reagents provided in these kits contain materials of human and/or animal origin. Whenever human plasma is required for the preparation of these reagents, the plasmas are tested for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with regulatory safety precautions in the manipulation of these biological materials as if they were infectious.

## WORKING REAGENT

Reconstitute each vial of Erba DDimer Control N+P with exactly 1 ml of distilled water / saline solution. Allow to stand for 5 minutes and mix well before use.  
Do not shake.

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.  
Once reconstituted, the control remains stable for:  
• 5 days at 2–8 °C.  
• 3 month at -20 °C.  
Do not freeze / thaw more than once.

## PROCEDURE

Each control should be treated in the same manner as the patients' sample in accordance with the instructions outlined in each particular assay procedure.

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10020391  
10020439

## EXPECTED VALUES

The control values for each of the parameters may vary from one lot to another. Refer to assay sheet section of this instruction. For greater control sensitivity each laboratory should establish its own mean and acceptable range and periodically re-evaluate the mean.

## QUALITY CONTROL

When properly used as recommended, the values of the Erba DDimer Control N+P should be found to be within the acceptable ranges, as indicated in the Assay Value Insert. If the stated control values for the lot being used cannot be reproduced, check that all the components of the test system are functioning correctly. If necessary, repeat the tests.

## LIMITATIONS

The values listed are for use with Erba D-Dimer procedures only. The use of other D-Dimer procedures and products may lead to erroneous results.

## ASSAY VALUE

LOT	XXXXXXX
EXP	MM/RRRR

Analyzer/Test		Unit	Control N	
			Value	Range
ECL 412	DDU <sup>1</sup>	ng/ml	000	000–000
ECL 105	FEU <sup>2</sup>	ng/ml	000	000–000
ECL 760	DDU <sup>1</sup>	ng/ml	000	000–000
	FEU <sup>2</sup>	ng/ml	000	000–000
Analyzer/Test		Unit	Control P	
			Value	Range
ECL 412	DDU <sup>1</sup>	ng/ml	000	000–000
ECL 105	FEU <sup>2</sup>	ng/ml	000	000–000
ECL 760	DDU <sup>1</sup>	ng/ml	000	000–000
	FEU <sup>2</sup>	ng/ml	000	000–000

<sup>1</sup> D-dimer Units (DDU) are assigned against D-dimer in-house reference material. Units of ng/ml = µg/l = 10<sup>-3</sup> µg/ml.

The range is specified in relation to the assigned concentration, in the each respective unit.

<sup>2</sup> Fibrinogen Equivalent Units (FEU) are calculated from DDU by multiplication with the correlation factor 2.5. Units of ng/ml = µg/l = 10<sup>-3</sup> µg/ml.

## REFERENCES

- Elms MJ et al (1986) Rapid Detection of Cross-Linked Fibrin Degradation Products in Plasma Using Monoclonal Antibody-Coated Latex Particles, J. Clin. Pathol, 85:360-64.
- Declerck PV et al (1987) Fibrinolytic Response and Fibrin Fragment D-dimer in Patients With Deep Vein Thrombosis, Thromb. Haemost, 58:1025-9.
- Ballegeer V et al (1987) Fibrinolytic Response to Venous Occlusion and Fibrin Fragment D-dimer Levels in Normal and Complicated Pregnancy, Thromb. Haemost, 58:1030-2.
- Hansson PO et al (1994) Can Laboratory Testing Improve Screening Strategies for Deep Vein Thrombosis at an Emergency Unit, J. Intern. Med, 235:143-51.

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- Elms MJ et al (1986) Rapid Detection of Cross-Linked Fibrin Degradation Products in Plasma Using Monoclonal Antibody-Coated Latex Particles, J. Clin. Pathol, 85:360-64.
- Declerck PV et al (1987) Fibrinolytic Response and Fibrin Fragment D-dimer in Patients With Deep Vein Thrombosis, Thromb. Haemost, 58:1025-9.
- Ballegeer V et al (1987) Fibrinolytic Response to Venous Occlusion and Fibrin Fragment D-dimer Levels in Normal and Complicated Pregnancy, Thromb. Haemost, 58:1030-2.
- Hansson PO et al (1994) Can Laboratory Testing Improve Screening Strategies for Deep Vein Thrombosis at an Emergency Unit?, J. Intern Med, 235:143-51.

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## REFERENCES

- Elms MJ et al (1986) Rapid Detection of Cross-Linked Fibrin Degradation Products in Plasma Using Monoclonal Antibody-Coated Latex Particles, J. Clin. Pathol, 85:360-64.
- Declerck PV et al (1987) Fibrinolytic Response and Fibrin Fragment D-dimer in Patients With Deep Vein Thrombosis, Thromb. Haemost, 58:1025-9.
- Ballegeer V et al (1987) Fibrinolytic Response to Venous Occlusion and Fibrin Fragment D-dimer Levels in Normal and Complicated Pregnancy, Thromb. Haemost, 58:1030-2.
- Hansson PO et al (1994) Can Laboratory Testing Improve Screening

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