## Erba Calcium Chloride

### Calcium Chloride Solution 0.025 M

Cat. No.:	Pack name:	Packaging (Content):
EHL00020	Erba Calcium Chloride	10 x 10 ml

#### INTENDED USE

Erba Calcium Chloride – 0.025 M  $\text{CaCl}_2$  is intended for use with Erba Actime reagent for the determination of the activated partial thromboplastin time (APTT) or for the assays of intrinsic pathway factors.

#### COMPOSITION

Erba Calcium Chloride: 0.025 M solution of calcium chloride.

#### WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only. These reagents are to be used by certified medical laboratory personnel only.
- Do not ingest.
- Wear gloves when handling all kit components.
- Only use clean or single use laboratory equipment to avoid contaminations.
- The eventual rest of reagents should be disposed of in accordance with the internal regulations and in compliance with local and national regulations relating to the safe handling of waste.

#### WORKING REAGENT

Ready to use.  
Avoid reagent contamination.

#### STABILITY AND STORAGE

Store between 2 and 8 °C. Stable until the expiry date stated on the label. A cloudy solution may be indicative of product degradation. Do not freeze.

#### SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2 % or 3.8 % sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept between 2–8 °C or between 18–24 °C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20 °C for 2 weeks or -70 °C for 6 months. Thaw quickly at 37 °C prior to testing. Do not keep at 37 °C for more than 5 minutes. This will minimize the neutralization of the lupus inhibitor.<sup>1</sup> Erroneous results may be caused by contamination with tissue fluids or stasis. Avoid agitation, air bubbles or foaming. For the effects of commonly administered drugs, refer to Young, et al.<sup>2</sup>

#### PROCEDURE

Refer to the instrument's operator's manual.

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
### LITERATURE / LITERATURA / REFERENCIAS / ЛІТЕРАТУРА

- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Young DS *et al.* Effects of Drugs on clinical Laboratory tests, 3rd ed., AACC Press, Washington, D.C., 1990
- Kelton JG and Hirsh J (1980) Bleeding Associated with Antithrombotic Therapy, *Semin Hematol*, 17: 259-91
- Estes J and Pouplin P (1974) Pharmacokinetics of Heparin, *Thromb Diath Haemorr*, 33:26-35
- Solomon G, Hilgartner M, Kutt H (1972) Coagulation defects caused by diphenylhydantoin, *Neurology*, 22: 1165-71
- Ambrus JL, Schimert G, Lajos TZ, Ambrus CM, Mink IB, Lassman HB, Moore RH, Melzer J (1971) Effect of antifibrinolytic agents and estrogens on blood coagulation factors during open heart surgery, *J Med*, 2: 65-81
- Crowell EB Jr, Clatanoff DV, Kiekhof W (1971) The effect of oral contraceptives on factors VIII levels, *Lab Clin. Med*, 77: 551-7
- Langdell R, Wagner R, Brinkhous K (1953) Effect of Antihemophilic Factors on One Stage Clotting Tests, *J. Lab. Clin. Med*, 41: 637
- Proctor R, Rapaport S (1961) The Partial Thromboplastin Time with Kaolin, *Am. J. Clin. Path*, 36: 212
- Brandt JT and Triplett DA (1981) Laboratory Monitoring of Heparin. Effect of Reagents and Instruments on the Activated Partial Thromboplastin Time, *AJCP*, 76: 530-537