Erba Protime LS

Protrombinový čas

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LITERATURA

1. Quick AJ (1935) A Study of the Coagulation Defect in Hemophilia and Jaundice, Am. J. Med. Sci, 190: 501.
2. Biggs R (1976) Human Blood Coagulation, Haemostasis and Thrombosis, 2nd Edition, Blackwell Scientific Publications, London.
3. Hirsh J, Poller L, Deykin D, Levine J, Dalen JE (1989) Optimal Therapeutic Range for Oral Anticoagulants, Chest, 95: 5S-11S.
4. Poller L (1986) Laboratory Control of Anticoagulant Therapy, Sem. Thromb. Haemostasis, 12: 13-19.
5. World Health Organisation (1984) Expert Committee on Biological Standards, Technical Series, 700: 19.
6. Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
7. Poller L., Triplett DA, Hirsh J, Carroll J, Clarke K (1995) The value of plasma calibrants in correcting coagulometer effects on International Normalised Ratios (INR): An international multicentre study, Amer. J. Clin. Pathol, 103: 358-365.
8. Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) A comparison of lyophilised artificially depleted plasmas and lyophilised plasmas from warfarin treated patients in correcting for coagulometer effects on International Normalised Ratios, Amer. J. Clin. Pathol, 103: 366-371.
9. Keeling D (2011) Guidelines on Oral Anticoagulation with warfarin: Forth Edition, British Journal of Haematology, 154(3): 311-324.

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INTENDED USE

Erba Protime LS is a liquid ready to use reagent, that contains rabbit brain thromboplastin intended for the determination of the prothrombin time (PT).

CLINICAL SIGNIFICANCE

Erba Protime LS thromboplastin is a rabbit brain extract, which contains tissue factor, lipids and calcium ions. The first standardized one-stage prothrombin time test was developed by Dr. Armand Quick in 1935.¹ It has now become the basic coagulation screening test for the diagnosis of congenital and acquired deficiencies of clotting factors from the extrinsic pathway (factors II, V, VII and X).² It is also used for the induction and monitoring of oral anticoagulant therapy³ and can be used to assess the protein synthesis capability of the liver in chronic or acute hepatic disorders.

Therefore, the PT by Erba Protime LS is optimally used for presurgical screening and monitoring for oral anticoagulant therapy (OAT)⁴ also. Erba Protime LS reagent with the corresponding deficient plasmas is also suitable for determination of activity of extrinsic coagulation pathway (factor II, V, VII and X).

PRINCIPLE

The one-stage PT measures the clotting time of plasma after adding a source of tissue factor (thromboplastin) and calcium. The recalcification of plasma in the presence of tissue factor generates activated factor Xa. Factor Xa in turn activates prothrombin to thrombin, which converts fibrinogen to an insoluble fibrin clot.

The time of this clotting process is measurable manually or with optical or mechanical coagulation analysers.

COMPOSITION

Erba Protime LS is a tissue thromboplastin from rabbit brain, which contains calcium ions and sodium azide (< 0,01 %) as preservative.

WARNINGS AND PRECAUTIONS

For *in vitro* diagnostic use only. These reagents are to be used by certified medical laboratory personnel only.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1

Contains: 2-methylisothiazol-3(2H)-one

Erba Protime LS 10,50 UFI: VCJV-EWQG-1J5S-73D0

Erba Protime LS 100 UFI: YEJV-XWDV-CJ58-VEY2

Warning

Hazard statement(s):

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

H412 Harmful to aquatic life with long lasting effects.

Precautionary statement(s):

P264 Wash hands and exposed parts of the body thoroughly after handling.

P280 Wear protective gloves/protective clothing/eye protection.

P314 Get medical advice/attention if you feel unwell.

WARNING – POTENTIAL BIOHAZARDOUS MATERIAL

Some reagents provided in these kits contain materials of human and/or animal origin.

Whenever human plasma is required for the preparation of these reagents, the plasmas are

tested for the antibodies to HIV 1. HIV 2 and HCV, and for hepatitis B surface antigen. and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with regulatory safety precautions in the manipulation of these biological materials as if they were infectious.⁵

WORKING REAGENT

Erba Protime LS is ready to use.

The contents of the vial should be mixed well before use (5 minutes).

Using of stirring bar is necessary!

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

Opened vials are stable:

- 2 months at 2–8 °C

- 5 days at 15 °C

- 6 hours at 37 °C

Do not freeze.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout.

Blood (9 parts) should be collected into 3.2 % or 3.8 % sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at 15–25 °C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20 °C for 2 weeks or -70 °C for 6 months. Thaw quickly at 37 °C prior to testing. Do not keep at 37 °C for more than 5 minutes.⁶

PROCEDURE

Erba Protime LS is a one-stage PT test which can be used with semi-automated coagulation analysers according to the protocol detailed below.

Manual Method

- Mix sufficient Erba Protime LS reagent to complete the anticipated testing for the day and incubate reagent at 37 °C no more than 4 hours.

- Add 50 µl of the patient plasma or control plasma into a reaction tube and incubate at 37 °C for 2 minutes.

- Add 100 µl of freshly mixed reagent and start simultaneously a timer.

- Note the time for clot formation to the nearest 0.1 seconds.

Automated Method

Refer to the instrument's operator's manual.

INTERPRETATION OF RESULTS

Erba Protime LS test results can be reported in the following units, lot specific sheet in the box will help in the calculation:

1. Seconds, which means the observed clotting time.

2. Percentage, which means the proportional part of the normal PT activity, which is calculable from the calibration curve. Method dependent master curve in the issued sheet can be used for the calculation.

3. International Normalized Ratio (INR), calculated as $INR = (PT / MNPT)^{ISI}$ Where PT means the clotting time of the sample, MNPT means Mean Normal Prothrombin Time and ISI means International Sensitivity Index. The ISI method dependent values in the lot specific values sheet can be used for the calculation.

- The INR is the only officially recognized dimension of the result at vitamin K antagonists treated patients. The normal range expressed in INR is 0,8–1,2.

- Every laboratory should determine its own MNPT value and reference range. MNPT value in the issued sheet is only for information, because it depends on the measuring circumstances and population.

- Accurate and general conversion of percentage into INR (or back) is not possible!

REFERENCES VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges. This is particularly important for local ISI calibration.^{7,8,9} Using an optical coagulation analyzer, normal values ranging from 11.5–14.6 seconds; 0.93–1.16 INR; 79.1–112.8 % PT are typical.

QUALITY CONTROL

Each laboratory should establish a quality control program.

To ensure adequate quality, control plasmas are recommended. It is suggested to use two levels of control, one close to the normal patient values (Erba Control N, Cat. No.: EHL00014 or

Erba Control N Plus, Cat. No.: EHL00016) and the second representative the pathologic values (Erba Control P, Cat. No.: EHL00015 or Erba Control P Plus, Cat. No.: EHL00017).

LIMITATIONS

Erba Protime LS is insensitive to the following substances:

Heparin: up to 2 IU/mL

Bilirubin: up to 0.5 g/l

Heamoglobin: up to 10 g/l

Triglycerides: up to 3.75 g/l

PERFORMANCES

These performances have been obtained using ECL coagulation analyzer. Results may vary if a different instrument or manual procedure is used.

Precision:

	Intra-assay precision (N = 20)		Inter-assay precision (N = 10)	
Mean (s)	12.4	25.6	12.7	24.8
CV (%)	0.49	1.46	2.18	1.95

REFERENCES

1. Quick AJ (1935) A Study of the Coagulation Defect in Hemophilia and Jaundice, Am. J. Med. Sci, 190: 501.
2. Biggs R (1976) Human Blood Coagulation, Haemostasis and Thrombosis, 2nd Edition, Blackwell Scientific Publications, London.
3. Hirsh J, Poller L, Deykin D, Levine J, Dalen JE (1989) Optimal Therapeutic Range for Oral Anticoagulants, Chest, 95: 5S-11S.
4. Poller L (1986) Laboratory Control of Anticoagulant Therapy, Sem. Thromb. Haemostasis, 12: 13-19.
5. World Health Organisation (1984) Expert Committee on Biological Standards, Technical Series, 700: 19.
6. Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
7. Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) The value of plasma calibrants in correcting coagulometer effects on International Normalised Ratios (INR): An international multicentre study, Amer. J. Clin. Pathol, 103: 358-365.
8. Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) A comparison of lyophilised artificially depleted plasmas and lyophilised plasmas from warfarin treated patients in correcting for coagulometer effects on International Normalised Ratios, Amer. J. Clin. Pathol, 103: 366-371.
9. Keeling D (2011) Guidelines on Oral Anticoagulation with warfarin: Forth Edition, British Journal of Haematology, 154(3): 311-324.

USED SYMBOLS

LOT	Lot Number	IVD	In vitro Diagnostics	I	See Instruction for Use
REF	Catalogue Number	M	Manufacturer	CONT	Content
E	Expiry Date	S	Storage Temperature		

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REFERENCIAS

1. Quick AJ (1935) A Study of the Coagulation Defect in Hemophilia and Jaundice, Am. J. Med. Sci, 190: 501.
2. Biggs R (1976) Human Blood Coagulation, Haemostasis and Thrombosis, 2nd Edition, Blackwell Scientific Publications, London.
3. Hirsh J, Poller L, Deykin D, Levine J, Dalen JE (1989) Optimal Therapeutic Range for Oral Anticoagulants, Chest, 95: 5S-11S.
4. Poller L (1986) Laboratory Control of Anticoagulant Therapy, Sem. Thromb. Haemostasis, 12: 13-19.
5. World Health Organisation (1984) Expert Committee on Biological Standards, Technical Series, 700: 19.
6. Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
7. Poller L., Triplett DA, Hirsh J, Carroll J, Clarke K (1995) The value of plasma calibrants in correcting coagulometer effects on International Normalised Ratios (INR): An international multicentre study, Amer. J. Clin. Pathol, 103: 358-365.
8. Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) A comparison of lyophilised artificially depleted plasmas and lyophilised plasmas from warfarin treated patients in correcting for coagulometer effects on International Normalised Ratios, Amer. J. Clin. Pathol, 103: 366-371
9. Keeling D (2011) Guidelines on Oral Anticoagulation with warfarin: Forth Edition, British Journal of Haematology, 154(3): 311-324

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