Erba Factor XII Deficient Plasma

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LITERATURA

- Penner JA (1979) The University of Michigan Medical School Blood Coagulation Laboratory Manual, 14th Ed., University Publications, Ann Arbor, 72-78.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
- Kirkwood TBL et al. (1977) Identification of Sources of Variation in Factor VIII Assay, British Journal of Haematology, 37:559-568.
- Goldenfarb MD (1971) Reproducibility in Coagulation Assays, AJCP, 55:561-564.
- Palkuti HA and Longberry JR (1973) A Precision Study of Coagulation Factor Assay Techniques, AJCP, 59:231-235.
- Triplett DA, Harms CS (1981) Procedures for the Coagulation Laboratory, Am. Society for Clin. Path., Chicago, 36

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INTENDED USE

Erba Factor XII Deficient Plasma is intended for the quantitative determination of the Factor XII in patients suspected of having a congenital or acquired deficiency of this coagulation protein.

PRINCIPLE

Quantitative measurement of individual coagulation factors by the one-stage method requires substrate plasma lacking the factor to be measured. A dilution of the test plasma is mixed with the factor deficient plasma and the clot time of the mixture determined. The degree of clot time correction with the patient plasma is compared to the correction with a reference material, allowing the % activity of the patient plasma to be determined.¹ The Erba Factor XII Deficient Plasma can be used on any instrument capable of performing APTT-based factor assay testing.

COMPOSITION

Erba Factor XII Deficient Plasma is derived from human plasma and contain less than 1 % residual factor activity.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only. These reagents are to be used by certified medical laboratory personnel only.
- Do not ingest.
- Wear gloves when handling all kit components.
- Only use clean or single use laboratory equipment to avoid contaminations.
- The eventual rest of reagents should be disposed of in accordance with the internal regulations and in compliance with local and national regulations relating to the safe handling of waste.

WARNING – POTENTIAL BIOHAZARDOUS MATERIAL

Some reagents provided in this kit contain materials of human and/or animal origin. Whenever human plasma is required for the preparation of these reagents, the plasmas are tested for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen. and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with regulatory safety precautions in the manipulation of these biological materials as if they were infectious.

WORKING REAGENT

Reconstitute each vial of Erba Factor XII Deficient Plasma with 1 ml of purified water. Swirl gently and allow to stand for 15 minutes at room temperature (18–25 °C). Mix well before use (do not shake).

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

Reconstituted reagent is stable:

- 8 hours at 2–8 °C.

The lyophilised product should appear as a dry, straw coloured plug or pieces.

REAGENTS REQUIRED BUT NOT PROVIDED

Erba Standard Plasma (Cat. No.: EHL00012)
Erba Owrens Veronal Buffer (Cat. No.: EHL00021)
Erba Actime (Cat. No.: EHL00003, EHL00004)
Erba Calcium Chloride (Cat. No.: EHL00020)

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2 % or 3.8 % sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes.

Plasma should be kept at 2–8 °C or 18–24 °C. Testing should be completed within 4 hours of sample collection.

Plasma can be storage frozen at -20 °C for 2 weeks or -70 °C for 6 months. Thaw quickly at 37 °C prior to testing. Do not keep at 37 °C for more than 5 minutes.²

PROCEDURE

Manual Method

- Prepare all reagents as instructed with each pack.
- Pre-warm Erba Calcium Chloride solution to 37 °C.

1. Standard Curve Preparation:

- Prepare the following dilutions in Erba Owrens Veronal Buffer:

Tube	Erba Standard Plasma (ml)	Owren's Buffer (ml)	Activity (%)
1	0.1	0.4	100
2	0.1	0.9	50
3	0.1	1.9	25
4	0.1	3.9	12.5

- Mix without shaking.

2. Patient Sample Preparation:

- Prepare a 1 + 4 dilution of the patient plasma or control plasma in Owrens Buffer.
- Mix without shaking.

3. Testing:

- Pipette, in duplicate, 0.05 ml of Erba Factor XII Deficient Plasma into a reaction tube.
- Add 0.05 ml of standard, patient or control plasma dilution and incubate at 37 °C for 2 minutes.
- Add 0.1 ml of APTT reagent and incubate for 5 minutes at 37 °C.
- Add 0.1 ml of 0.025 M calcium chloride solution while simultaneously starting a stopwatch.
- Determine the clot time for each of the standard, control or patient dilutions.
- Plot % Activity (X-axis) versus Mean Clot Time (Y-axis) for the standards on 2 cycle log-log graph paper.

Automated Method

Refer to the instrument's operator's manual.

INTERPRETATION OF RESULTS

Interpolate patient values from the curve. Calculate exact patient or control values by correcting for differences in calibration plasma values as follows:

Correct Activity (%) =

(Calibration plasma Reference Value / 100) x Interpolated Patient or Control Value

REFERENCES VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own normal range. Expected values for factor activity are 50–150 %.⁶

QUALITY CONTROL

Each laboratory should establish a quality control program.

To ensure adequate quality, control plasmas are recommended. It is suggested to use two levels of control, one close to the normal patient values (Erba Control N Plus, Cat. No.: EHL00016) and the second representative the pathologic values (Erba Control P Plus, Cat. No.: EHL00017).

Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

LIMITATIONS

The results obtained with Erba Factor XII Deficient Plasma depend on several factors strongly associated with instrumentation, types of reagents, deficient substrates and laboratory to laboratory variations.^{3, 4, 5} Each laboratory should establish an expected range for the particular instrument-reagent system.


PERFORMANCES

Assumed linearity of standard curve is from 10–150 %. Within run and between run precisions are expected to be < 5 %, using a range of automated instruments.

REFERENCES

- Penner JA (1979) The University of Michigan Medical School Blood Coagulation Laboratory Manual, 14th Ed., University Publications, Ann Arbor, 72-78.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
- Kirkwood TBL et al. (1977) Identification of Sources of Variation in Factor VIII Assay, British Journal of Haematology, 37:559-568.
- Goldenfarb MD (1971) Reproducibility in Coagulation Assays, AJCP, 55:561-564.
- Palkuti HA and Longberry JR (1973) A Precision Study of Coagulation Factor Assay Techniques, AJCP, 59:231-235.
- Triplett DA, Harms CS (1981) Procedures for the Coagulation Laboratory, Am. Society for Clin. Path, Chicago, 36

USED SYMBOLS

 REF	Catalogue Number	 IVD	In vitro Diagnostics	 I	See Instruction for Use
 LOT	Lot Number		Manufacturer	 CONT	Content
	Expiry Date		Storage Temperature		

 Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erbamannheim.com, www.erbamannheim.com

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