URINORM XL

Cat. No.: REG00060

Store at (+2 to +8) °C

| Lot No. | Expiry Date |
|---------|-------------|
| xxxxxxx | YYYY-MM |

INTENDED USE

URINORM XL control urines are liquid synthetic material with concentrations in the normal (URINORM XL N) and pathological (URINORM XL P) ranges. URINORM XL control urines are designed to verify the precision and accuracy of PHAN AUTO diagnostic strips and automatic urine analyzer LAURA XL intended for *in vitro* diagnostics. This product is intended also for evaluation of sediment. For professional use only.

PRODUCT DESCRIPTION:

URINORM XL control urines are human urine-based containing constituents of human and animal origin as well as purified chemicals. URINORM XL control urines contain preservatives added to maintain product stability and avoid the microbial growth. URINORM XL urine control is a ready-to-use.

REAGENTS:

URINORM XL N – normal level 120 ml

URINORM XL P – pathological level 120 ml

STORAGE AND STABILITY:

URINORM XL control urines are stable until the expiration date unopened stored at 2–8 °C and protected from light and microbial contamination.

Once opened URINORM XL control urines are stable for 30 days at 20–25 °C and 18 months at 2–8 °C. Do not store the product in freezers.

INSTRUCTIONS:

Control urines URINORM XL should be treated the same as patient specimens.

1. Allow the control urines bottles to reach the room temperature.
2. Gently invert the bottles to homogenize the sample.
3. Open the cap and pour 1–1.5 mL into the Eppendorf tube.
4. Perform the QC measurement on LAURA XL urine analyzer – device automatically evaluates the results and compare them with declared values.
5. Store appropriately.

DECLARED VALUES:

The following table summarizes declared values for each analyte for automatic urine analyser LAURA XL.

| CHEMISTRY | | | |
|-------------------|---------------|-----------------------|----------------------|
| SCALE 1 | | | |
| Analyte | URINORM XL N | URINORM XL P | Units |
| Specific gravity | X – X | X – X | |
| Nitrites | X | X | |
| pH | X – X | X – X | |
| Protein | X | X – X / X – X / X – X | g/L / mg/dL / ARB |
| Ascorbic acid | X | X – X / X – X / X – X | mmol/L / mg/dL / ARB |
| Glucose | X | X – X / X – X / X – X | mmol/L / mg/dL / ARB |
| Ketones | X | X – X / X – X / X – X | mmol/L / mg/dL / ARB |
| Urobilinogen | X | X – X / X – X / X – X | μmol/L / mg/dL / ARB |
| Bilirubin | X | X – X / X – X / X – X | μmol/L / mg/dL / ARB |
| Leucocytes | X | X – X / X – X / X – X | Leu / μL / ARB |
| Blood | X | X – X / X – X / X – X | Ery / μL / ARB |
| SEDIMENT | | | |
| Red blood cells | X – X / X – X | X – X / X – X | N/hpf / N/μL |
| White blood cells | X – X / X – X | X – X / X – X | N/hpf / N/μL |

| CHEMISTRY | | | |
|-------------------|---------------|-----------------------|----------------------|
| SCALE 2 | | | |
| Analyte | URINORM XL N | URINORM XL P | Units |
| Specific gravity | X – X | X – X | |
| Nitrites | X | X | |
| pH | X – X | X – X | |
| Protein | X | X – X / X – X / X – X | g/L / mg/dL / ARB |
| Ascorbic acid | X | X – X / X – X / X – X | mmol/L / mg/dL / ARB |
| Glucose | X | X – X / X – X / X – X | mmol/L / mg/dL / ARB |
| Ketones | X | X – X / X – X / X – X | mmol/L / mg/dL / ARB |
| Urobilinogen | X | X – X / X – X / X – X | μmol/L / mg/dL / ARB |
| Bilirubin | X | X – X / X – X / X – X | μmol/L / mg/dL / ARB |
| Leucocytes | X | X – X / X – X / X – X | Leu / μL / ARB |
| Blood | X | X – X / X – X / X – X | Ery / μL / ARB |
| SEDIMENT | | | |
| Red blood cells | X – X / X – X | X – X / X – X | N/hpf / N/μL |
| White blood cells | X – X / X – X | X – X / X – X | N/hpf / N/μL |

LIMITATIONS:

Control urines should not be used past the expiration date on the label. Do not use turbid control urines where there is an evidence microbial contamination.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled solely by an authorized, trained personnel. Human source material. Treat as potentially infectious.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagents of the kit are not classified as dangerous.

FIRST AID:

In case of an accidental ingestion wash up the mouth and drink approximately 0.5 L of water. In case of eye contact rinse the eye quickly and thoroughly using a stream of clean water. Contaminated skin is to be washed using warm water and soap. Consult a physician in all serious cases of health damage.

WASTE DISPOSAL:

Residues of control urines are to be disposed of in accordance with the internal regulations and by-laws. Paper and other packaging materials are to be discarded in the form of sorted waste (paper, glass, plastic).

Date of revision: 20. 2. 2024

USED SYMBOLS

| | | | | | |
|-----|------------------|---------------------|----------------------|--------------|--------------------------|
| REF | Catalogue Number | IVD | In Vitro Diagnostics | Manufacturer | See Instructions for Use |
| LOT | Lot Number | Storage Temperature | Expiry Date | | |

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