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Σ 50

DIAPHAN[®] FOR SELF-TESTING

Carefully follow the instructions for use!
Not intended for internal use.

REF 10006983

Introduction:
The diagnostic test strips DIAPHAN[®] For Self-testing are intended for the semi-quantitative determination of glucose and ketone concentrations in urine. They are intended primarily for patients with diabetes who are not treated with insulin. The frequency of measurement is individual based on consultation with a physician.
In general, regular checking 1-2x per week is recommended (fasting). In an initial phase of usage or in case of decompensation it is recommend to check daily in the morning (fasting) and then 1-2 hours after meals.
For diabetes on insulin who have a glucometer, DIAPHAN[®] For Self-testing can be recommended as an additional self-monitoring tool or on days when they do not measure glucose, i.e. morning urine checking before the application of insulin.

Reagents: 1 test contains
Glucose: glucose oxidase 1.3 %; peroxidase 1.3 %; tetramethylbenzidine 21 %
Ketones: sodium nitroprusside 4.9 %

Principle of test:
Glucose - The test is based on specific glucose oxidase/peroxidase reaction and is specific for D-glucose. The reagent pad does not react with other sugars. The reaction is independent of pH and the presence of ketone bodies. The presence of reducing substances that may be present in the urine, such as ascorbic acid, does not cause a false negative reaction.
Ketones - The test is based on the principle of Legal's test and is more susceptible to acetoacetic acid than to acetone. The test does not react with β-hydroxybutyric acid. Drugs and diagnostics on the basis of phenolphthalein or sulphophthalein may cause the reagent pad to turn red to purple due to an alkaline reaction. The colour scale is calibrated for acetoacetic acid.

Product composition:
Diagnostic test strips in primary protective packaging (tagged tube) and secondary packaging (box) containing instructions for use.

Storage:
Keep diagnostic test strips in tightly closed original tubes at +2 to +30°C. The strips must be protected from frost, high humidity (see caution to close the tube carefully after removing the strip) and elevated temperature. Do not store the product in the refrigerator or near heat sources. Also, protect from direct sunlight to avoid fading of the colour scale on the label. The product is usable until the last day of the month of the year of use indicated on the label and the box.

Notice:
One diagnostic strip is designed for a single usage. Do not use the same strip repeatedly. For the proper functioning of the product, the storage conditions, the urine sampling procedure, and the test procedure must be strictly followed. Do not use the product after the expiry date mentioned on the package and do not touch the strip at the test zone.
The manufacturer does not expect the use of diagnostic strips by children and people with mental disabilities.
Consideration should be given to the possibility of misinterpretation of results in the case of individual color perception, color sensitivity disturbance or color blindness.

Equipment:
Clean and dry collection container (glass, plastic), without trace of detergents and disinfectants. Stop-watch or watch with a second hand.

Urine sampling procedure:
Wash your hands with soap and water.
Take the urine collection container to the toilet, do not touch the inside of the collection container with your fingers.

Test procedure:
Immediately before use, remove the strip from the tube. (Pic. 1)
Tightly close the tube with the original cap containing a drying agent.
Do not touch the reagent zone of the strips.
Place the strip vertically in the urine for 1-2 seconds so that the zone soaks the urine. (Pic. 2)
Wipe the strip on the rim of the collecting container to remove excess urine. (Pic. 3)
Place the strip on a horizontal, clean, dry, and non-absorbent pad.
After 60 seconds, evaluate the colour of the reagent zone by aligning it with the colour scale on the label, in daylight if possible. (Pic. 4)

Result and its meaning:
Glucose
A negative result is represented by the zone at the beginning of the colour comparison scale on the label marked „neg“. Glucose-free urine will provide a colour response corresponding to the shade that is labelled as „neg“ on the label.
Positive results are represented by the additional 4 colour zones of the comparison scale on the label with the indication from 2.8 to 55 mmol/l (or from 50 to 1000 mg/dl) respectively.

Ketones
A negative result is represented by the zone at the beginning of the colour comparison scale on the label marked „neg“. Ketone-free urine provides a very light grey colouring of the diagnostic zone, which is labelled as „neg“ on the comparison scale on the label.
Positive results are represented by the additional 3 colour zones of the comparison scale on the label with the indication from 1.5 to 15 mmol/l (or from 16 to 156 mg/dl) respectively. In presence of ketones in urine, the colour of the diagnostic zone changes significantly to pink to pinkish red. The colour intensity depends on the concentration of ketones in urine.

Inform your doctor about any positive result of glucose or ketones in urine.
Insulin doses must not be adjusted based on the results of urine testing!

Limitations and possible errors:
If the urine is too diluted, the sample may not contain enough glucose.
For this reason, it is desirable to collect and test the first-morning urine, where the concentration of the individual analytes is usually the highest and steady.
The influence of drugs or their metabolites on the glucose test has not yet been fully elucidated. Always consult the controversial results with your doctor. If necessary, the doctor may recommend retesting urine after discontinuing of medication. The sensitivity of the test is influenced by the variability of urine composition.
As with all diagnostic tests, a definitive diagnosis should not be determined on the basis of a single test, but should only be determined by a physician after all clinical and laboratory findings have been evaluated.

Interference:
Glucose Zone:
False negativity - ascorbic acid (above 200 mg/l), gentisic acid, dihydroxy-phenylalanine, acetic acid (above 5 g/l), bilirubin (60 m/l/l).
In women we avoid examination during menstruation.

Health protection:
The product does not contain toxic substances.
The principles of personal hygiene should be taken into account when assessing the test.

Disposal:
Pour urine into a toilet and flush. Throw empty packaging and used strips into municipal waste.
Performance characteristics: Available on request from the manufacturer.

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