# PHAN® AUTO

  DIAGNOSTIC TEST STRIPS FOR URINALYSIS FOR ANALYSERS LAURA XL / ДИАГНОСТИЧЕСКИЕ ПОЛОСКИ ДЛЯ ИССЛЕДОВАНИЯ МОЧИ АНАЛИЗАТОРАМИ LAURA XL / DIAGNOSTICKÉ PŘOUŽKY K VYŠETŘENÍ MOČE ANALÝZÁTORY LAURA XL / DIAGNOSTICKÉ PRŮŽKY NA VYŠETRENIE MOČU ANALÝZÁTORMÍ LAURA XL / PASKI DIAGNOSTYCZNE DO BADANIA MOCZU ZA POMOCĄ ANALIZATORÓW LAURA XL

	REF		SG	LEU	NIT	pH	ASC	PRO	GLU	KET	UBG	BIL	BLD	CP
DEKAPHAN® AUTO	URPH0030	100	■	■	■	■		■	■	■	■	■	■	■
UNDEKAPHAN® AUTO	URPH0031	100	■	■	■	■	■	■	■	■	■	■	■	■

RFID

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 **Utilisation:**  
The diagnostic test strips PHAN® AUTO are intended for the semiquantitative analysis of the urine. The diagnostic test strips PHAN® AUTO are intended for *in vitro* diagnostic for professional use only.  
**Instruction:**  
For the objective evaluation of the diagnostic test strips, please, use the user manual for reader LAURAXL®.  
Use only a freshly voided, well-mixed, uncentrifuged specimen without preservatives collected in a clean container. The urine should not be more than 4 hours old when tested.  
**Notes:** For visual evaluation compare the tests pads to the corresponding colour scale on the label after approx. 60 sec. and for leucocytes after 120 sec. According to the different spectral sensitivity of the human eye and the optical system it is not possible to guarantee the exact correspondence between visual reading and reading by analyser.  
1. Remove only as many test strips as are required, and reseal the tube immediately after use.  
2. Do not touch test pads of the strips.  
3. Completely immerse all reagent pads in specimen (no longer than 1–2 sec.).  
4. Run edge of the strip against rim of urine container to remove excess urine. Hold the strip in horizontal position.

**Working reagent concentrations:**  
**Specific gravity:** poly(methylvinylether/maleic acid) 32%; bromthymol blue 5.1% / **Leucocytes:** indoxyl ester 0.43%; diazonium salt 0.05%  
**Nitrite:** sulphanilamide 5.1%; tetrahydrobenzo-[h]-quinoline 5.8% / **pH:** methyl red 0.71%; bromthymol blue 12.1%; **Ascorbic Acid:** phosphomolybdic acid 26%  
**Protein:** tetra bromphenolphthalein ester 0.21%; tetra bromphenol blue 0.35%; **Glucose:** glucose oxidase 1.3%; peroxidase 1.3%; tetramethylbenzidine 21.0%  
**Ketones:** sodium nitroprusside 4.9% / **Urobilinogen:** diazonium salt 2.3% / **Bilirubin:** diazonium salt 0.75% / **Blood:** tetramethylbenzidine 1.5%; cumene hydroperoxide 15.2%  
**Principle:**  
**Specific gravity** – The test is based on the principle of ion exchange, which runs between polyelectrolyte and ions present in the urine. Its result is a colour change of an acid-base indicator from the blue-green colour in the urine with low concentration of ions, through green and yellow-green in the urine with increased concentration of ions, to umber yellow colour. Using this test it is possible to determine the specific gravity of the urine in the range of 1.000 up to 1.030.  
**Leucocytes** – The test is based on the enzymatic reaction. The test pad contains an indoxyl ester, that is cleaved by the granulocyte esterases. The released indoxyl reacts with a diazonium salt and violet coloration is formed. The colour intensity is proportional to the leucocytes amount in a sample of tested urine and is evaluated after 120 seconds.  
**Nitrites** – The test is based on conversion of nitrate to nitrite by the action of certain species of bacteria contained in the urine. The colour test is based on the principle of the Griess test. Some degree of pink coloration should be interpreted as a positive nitrite test suggestive of 10<sup>5</sup> or more organisms in 1 ml of the urine specimen.  
**pH** – The test is based on the double indicator principle and gives a range of colours from orange through yellow and green to blue and permits differentiation to within 0.5 pH unit in the range pH 5 to 9.

**Ascorbic Acid** – The test is based on the reaction of phosphomolybdic acid which is reduced by ascorbic acid to molybdenum blue. The test is not specific for ascorbic acid because the green to greyish-blue colour of the test pad is exhibited also by other strongly reducing substances present in urine, such as gentisic acid and other acetylsalicylic acid metabolites. We recommend to carry out determination of ascorbic acid in urine especially in cases, in which ascorbic acid may disturb the tests for other urine constituents, such as glucose, blood and nitrite.  
**Proteins** – The test is based on colour change of acid-base indicator, which is caused by presence of proteins. It is particularly sensitive to albumin, but is much less sensitive to globulin, mucoprotein, haemoglobin and Bence-Jones protein.  
**Glucose** – The test is based on the specific glucose oxidase/peroxidase reaction and is specific for D-glucose. The reagent pad does not react with other sugars, it reacts with presence of D-glucose by green to the dark green coloration.  
**Ketones** – The test is based on the principle of Legal' s test and is more susceptible to the acetoacetic acid than to acetone. Test does not react with the β<sub>2</sub>-hydroxybutyric acid. The colour scale is calibrated for the acetoacetic acid.  
**Urobilinogen** – The test is based on the coupling of urobilinogen with a stabilized reagent. The test is specific for urobilinogen and stercobilinogen and is not susceptible to the interfering factors known at Ehrlich's test.

**Bilirubin** – The test is based on the coupling of bilirubin with stabilized reagent. The reaction is not affected by pH of the urine.  
**Blood** – The test is based on the peroxidase activity of haemoglobin which catalyzes the oxidation of the indicator due to the presence of the organic hydroperoxide contained in the diagnostic pad. The label contains two colour scales; for detection of intact erythrocytes and free haemoglobin. The test is highly sensitive to free haemoglobin and may detect its presence from concentrations corresponding approx. to 5 Ery/μl.  
**Compensation field** – Pad, which isn't impregnated with any reagents. The compensation field is used for suppressing of the dark colour of the urine sample, because the dark colour could have the effect on the evaluation of reagent pads.

**Limitations:**  
**Specific gravity** – This reaction is affected by pH values of urine over 6.5, this can cause shift in colour response towards lower values of specific gravity.  
**Leucocytes** – In case when the urine sample is more markedly coloured (e.g. increased bilirubin content), the resulting colour could be affected by a sample coloration. The intensity of the colour reaction is increased by alkaline pH and higher urine density.  
**Nitrites** – Before testing the patient should intake vegetable-rich meals and discontinue antibiotic therapy for 3 days prior to the test. Sensitivity of this test declines with the high specific gravity of the urine. Increased diuresis can cause the false negative results. Limited fluid intake prior testing can prevent from the excessive dilution of the urine. The test can be applied only at the fresh urine. Inaccurate results may occur at the stale urines, in which nitrite can be formed by contamination of the specimen.  
**Proteins** – In strongly alkaline urines (pH>8) from patients on medication with quinine or quinine containing drugs false positive reading may be obtained. False positive results may be found when the urine collection vessel contains traces of disinfectants with quarternary ammonium groups. On the other hand, in the presence of non-ionic or anionic detergents, false-negative results may occur. Do not take the colour of the dry pad into consideration.  
**Glucose** – The reaction is independent on pH and presence of ketone bodies.  
**Ketones** – Drugs and diagnostics on the basis of phenolphthalein or sulphophthalein may turn red to purple because of alkaline reaction of the pad.  
**Urobilinogen** – The reaction is not affected by pH of the urine. The presence of bilirubin gives yellow colour. This colour which turns slowly to greenish-blue does not interfere with urobilinogen determination provided the reading is made 1 minute after wetting. The urine specimen should not be exposed to the direct sunlight as this promotes the oxidation of urobilinogen and thus leads to artificially low or false negative readings.

**Bilirubin** – The urine specimen should not be exposed to direct sunlight as this promotes the oxidation of bilirubin and thus leads to artificially low or false negative readings. High concentrations of urobilinogen (above 100 μmol/l) interfere with the test. Also the red substances or the substances, which are turning red in low pH may interfere (e.g. Phenazopyridine).  
**Blood** – Microbial peroxidase associated with urinary tract infection may cause a false positive reaction. The sensitivity of the test is influenced by specific gravity or by inhibitors of pharmaceutical origin.  
All diagnostic pads do not interfere with the common concentration of the ascorbic acid. More information are stated at web site [www.erbalachema.com](http://www.erbalachema.com)  
**Reference values:** The reference values are written in the Table I. The reference values is only approximate, it is recommended that each laboratory verify the reference values for their particular examined population.  
**The results:** The results for reader LAURA XL® are written in the Table I. The colour scale for visual reading is on the label.

**Please Note:** Knowledge of the effects of drugs or their metabolites upon the individual tests is not yet complete. In doubtful cases, it is advisable to repeat the test after dis-continuing a drug. For the objective evaluation by using the analyser LAURA XL® the sensitivity of tests can be depended upon the variability of urines. The diagnostic test strips PHAN® AUTO can only be read by analyser LAURA XL® (see the scheme of the diagnostic strips above), they are not suitable for use with other instrument for analyse of urine. The analyser LAURA XL® use other method for evaluation of urine, than is the visual evaluation, therefore the concentration scale measured with analysers LAURA XL® are different from the colour scale on the label on the tube. This colour scale is determined for the orientation visual evaluation only. The semiquantitative analysis is not sufficient for the completing of diagnosis.  
**Quality control:**  
For quality control of precision and accuracy URINORM XL control urines (Cat. No. REG00060) are recommended. Obtained results are compared with the declared values in the instruction for use of URINORM XL Bio-Rad Liquechik Urinalysis Control can be also used as third-party control material. It is recommended to perform QC measurements according to local laboratory guidelines. More information about quality control can be found in LAURA XL user manual.  
**Storage:** Keep the diagnostic test strips in tightly closed original tubes in a dry and dark place at (+2 to +30) °C. The strips must be kept away from moisture, direct sunlight, elevated temperature and chemical fumes in the laboratory. When stored under these conditions, test strips are stable to the expiry date given on the pack.  
**Waste disposal:** Used strip should be treated as potentially infectious and should be liquidated in accordance with local and national regulations relating to the safe handling of such materials. Let waste recycle or put it to municipal waste.

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