Erba H360 Lyse

Haematology Analyser Lyse

Cat. No.	Product Name	Packaging (Content)
HEM00029	Erba H360 Lyse	500 mL, Instruction for use, RFID card

INTENDED USE

Erba H360 Lyse is an *in vitro* diagnostic reagent used to lyse red blood cells for quantitative counting of white blood cells, 3-part differentiation of white blood cells and measuring of haemoglobin in human blood on automated haematology analysers H360. It is designed to report on the haematological status of patient samples, supporting the monitoring and diagnosis of pathological conditions. The Erba H360 Lyse should be used by trained professionals in the haematology laboratory.

PRINCIPLE

The product is added to lyse the red blood cells (RBC) and react with the released haemoglobin (HGB) for the measurement of HGB, white blood cell (WBC) count and differential.

DESCRIPTION AND COMPOSITION

Surfactant; Buffering Agent; Antifungal and antibacterial agents.

PREPARATION

This reagent is ready for use and can be applied straight from the container; no special reagent preparation is necessary. Avoid reagent contaminations. Before using, leave the reagent at room temperature for at least 24 hours. For reagent connection to the analyser, refer to the H360 User manual.

MATERIAL REQUIRED BUT IS NOT PROVIDED WITH THE DEVICE

Haematology analyser H360, Cat. No. INS00077

Reagent Erba H360 Dil, Cat. No. HEM00028

SPECIFICATION

The product is a clear liquid without particles, precipitate or floccule. Background result: WBC $\leq 0.2 \times 10^9/L$, RBC $\leq 0.02 \times 10^{12}/L$, HGB $\leq 1 g/L$, PLT $\leq 10 \times 10^9/L$.

STABILITY AND STORAGE

The product is stable till the expiry date mentioned on the label when stored at (+2 to +30) °C. Discard opened container after 90 days.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and / or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Product is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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