Činidlo pro test ACETOIN

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Reagent for ACETOINtest

Cat. No.	Product Name	Packaging (content)
MLT00016	Reagent for ACETOINtest	Reagent for VPT I – weight: 1×1.8 mg, Reagent for VPT I – solvent: 1×18 mL, Reagent for VPT I – dropping bottle, Reagent for VPT II: 1×18 mL, Instruction for use

INTENDED USE

Reagent for ACETOIN test is an *in vitro* diagnostic reagent used as supplementary colour forming reagent for some of MIKROLATEST® kits with pure bacterial cultures obtained from various clinical materials in the general population. The kits are dedicated for visual evaluation or for determination on automatic system ErbaExpert. The kit is intended as assistance in establishing a diagnosis or pathological state which can be caused by the tested pathogen. For professional use in clinical laboratory only.

PRINCIPLE

Reagent for ACETOIN test is colour forming reagent for the test for the formation of acetoin (the Voges-Proskauer reaction).

DESCRIPTION AND COMPOSITION

Reagent for VPT I – weight
α-naphthol
Reagent for VPT I – solution
Ethanol anhydrous, denatured 5 % methanol
Reagent for VPT II
Potassium hydroxide

WORKING SOLUTION PREPARATION

Reagent for VPT I – Pour the content of the bottle labelled Reagent for VPT I – solvent into the bottle labelled Reagent for VPT I – weight. After the content has dissolved, pour carefully the solution into the dropping bottle labelled Reagent for VPT I. It is necessary to record immediately the DATE OF PREPARATION and EXPIRY (see STABILITY AND STORAGE) in the bottle label.
Reagent for VPT II – ready to use. Avoid to contamination.

STABILITY AND STORAGE

Reagent for the ACETOIN test should be stored in the dark at the temperature of (+2 to +25) °C. Expiration is marked at each package.
The dropping bottle with the reagent for VPT I must be stored in the dark at the temperature of (+2 to +8) °C, the prepared solution is usable up to 4 months since the date of preparation – not to exceed the expiry date of the kit.
The dropping bottle with reagent for VPT II can be stored at the temperature of (+2 to +25) °C.

ASSAY PROCEDURE

The Reagent for ACETOIN test should be used in accordance with the instructions for use of the specific kit.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagent for VPT I – weight

Contains: 1-naphtol (Index: 604-029-00-5; CAS: 90-15-3)

Hazard statement:

H315 Causes skin irritation.
H318 Causes serious eye damage.
H335 May cause respiratory irritation.
H302+H312 Harmful if swallowed or in contact with skin.

Precautionary statement:

P261 Avoid breathing dust.
P264 Wash hands and exposed parts of the body thoroughly after handling.
P280 Wear protective gloves/protective clothing/eye protection.
P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P308+P313 IF exposed or concerned: Get medical advice/attention.

Reagent for VPT I – solvent

UFI: KSCS-DWAR-EJ5N-U6DX

Contains: methanol

Hazard statement:

H225 Highly flammable liquid and vapour.
H302 Harmful if swallowed.
H319 Causes serious eye irritation.
H371 May cause damage to organs.

Precautionary statement:

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P261 Avoid breathing vapours.
P280 Wear protective gloves/protective clothing/eye protection.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P308+P313 IF exposed or concerned: Get medical advice/attention.
P370+P378 In case of fire: Use powder extinguisher/sand/carbon dioxide to extinguish.
Dropping bottle for Reagent VPT I:

Contains: 1-naphtol, methanol
Hazard statement:
H225 Highly flammable liquid and vapour.
H302 Harmful if swallowed.
H315 Causes skin irritation.
H318 Causes serious eye damage.
H371 May cause damage to organs.

Precautionary statement:

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P261 Avoid breathing vapours.
P280 Wear protective gloves/protective clothing/eye protection.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P308+P313 IF exposed or concerned: Get medical advice/attention.
P370+P378 In case of fire: Use powder extinguisher/sand/carbon dioxide to extinguish.

Reagent for VPT II:

Contains: potassium hydroxide, (Index: 604-029-00-5; CAS: 90-15-3)

Hazard statement:

H290 May be corrosive to metals.
H302 Harmful if swallowed.
H314 Causes severe skin burns and eye damage.

Precautionary statement:

P280 Wear protective gloves/eye protection/face protection.
P301+P330+P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P308+P313 IF exposed or concerned: Get medical advice/attention.

WASTE MANAGEMENT

Please refer to local legal requirements.

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