Reagent for HIPPURATetest

Cat. No.	Product name	Packaging (Content)
MLT00019	Reagent for HIPPURATetest	Reagent for HIP – weight: 2 × 0.9 g, Reagent for HIP – solvent: 2 × 18 mL, Reagent for HIP – dropping bottle, Instructions for use

INTENDED USE

Reagent for HIPPURATetest is an *in vitro* diagnostic reagent used as supplementary colour forming reagent for the diagnostic strips HIPPURATetest with pure bacterial cultures obtained from various clinical materials in the general population. The kits are dedicated for visual evaluation. The kit is intended as assistance in establishing a diagnosis or pathological state which can be caused by the tested pathogen. For professional use in clinical laboratory only.

PRINCIPLE

Reagent for HIPPURATetest is colour forming reagent for the test detecting the sodium hippurate hydrolysis and serves for colour expression of this test.

DESCRIPTION AND COMPOSITION

Reagent for HIP – weight

Ninhydrin

Reagent for HIP – solution

Aceton : Butanol

WORKING SOLUTION PREPARATION

Pour the content of the bottle labelled Reagent for HIP – solvent into the bottle labelled Reagent for HIP – weight. After the content has dissolved, pour carefully the solution arisen into the dropping bottle labelled Reagent for HIP. It is necessary to record immediately the DATE OF PREPARATION AND EXPIRY (see STABILITY AND STORAGE) in the bottle label. After the content of the dropping bottle Reagent for HIP has been used up, prepare again the solution in the same way, using the contents of remaining two flasks. Again record immediately the DATE OF PREPARATION – 2 and EXPIRY – 2 (see STABILITY AND STORAGE) in the bottle label.

STABILITY AND STORAGE

Reagent for the HIPPURATetest should be stored in the dark at the temperature of (+2 to +25) °C. Expiration is marked at each package.

The dropping bottle with the reagent for HIP must be stored in the dark at the temperature of (+2 to +8) °C, the prepared solution is usable up to 4 months since the date of preparation – not to exceed the expiry date of the kit.

ASSAY PROCEDURE

The Reagent for HIPPURATetest should be used in accordance with the instructions for use of the specific kit.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagent for HIP – weight

Warning

Hazard statement:

H302 Harmful if swallowed.

H315 Causes skin irritation.

H319 Causes serious eye irritation.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Reagent for HIP – solvent

UFI: DHX7-WJ9Q-NE7W-EKKG

Danger

Contains: acetone, butan-1-ol

Hazard statement:

H225 Highly flammable liquid and vapour.

H302 Harmful if swallowed.

H315 Causes skin irritation.

H318 Causes serious eye damage.

H335 May cause respiratory irritation.

H336 May cause drowsiness or dizziness.

Precautionary statement:

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P281 Avoid breathing dust/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection.


P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P308+P311 IF exposed or concerned: Call a POISON CENTER or doctor.

EUH066 Repeated exposure may cause skin dryness or cracking.

Dropping bottle for Reagent HIP

Danger

Contains: acetone, butan-1-ol, ninhydrine

Hazard statement:

H225 Highly flammable liquid and vapour.

H302 Harmful if swallowed.

H315 Causes skin irritation.

H318 Causes serious eye damage.

H335 May cause respiratory irritation.

H336 May cause drowsiness or dizziness.

Precautionary statement:

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P281 Avoid breathing dust/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P308+P311 IF exposed or concerned: Call a POISON CENTER or doctor.

EUH066 Repeated exposure may cause skin dryness or cracking.

WASTE MANAGEMENT

Please refer to local legal requirements.

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