Činidlo pro test HIPPURÁT

Kat. č.: MLT00019

Pro mikrobiologii

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Cat. No.: MLT00019

For microbiology

The set of the Reagent for HIPPURATtest is an auxiliary colour-producing reagent for the diagnostic strips MIKROLATEST® HIPPURATtest and for MIKROLATEST® kit STREPTtest 16, intended for a test, detecting the sodium hippurate (HIP) hydrolysis. The directions for use of the reagent are included in the operating instructions for the corresponding preparation of the MIKROLATEST®.

Preparation for use:

Pour the content of the bottle labelled Reagent for HIP-Solvent into the bottle labelled Reagent for HIP-weight. After the content has dissolved, pour carefully the solution arisen into the dropping bottle labelled Reagent for HIP. It is necessary to record immediately the DATE OF PREPARATION – 1 and EXPIRY – 1 on the bottle label (see Storage, expiration). After the content of the dropping bottle Reagent for HIP has been used up, prepare again the solution in the same way, using the contents of remaining two flasks. Again record the DATE OF PREPARATION – 2 and EXPIRY – 2 (see Storage, expiration).

Reagent for HIPPURATtest contains:

Reagent for HIP – weight	2 x 0,9 g
Reagent for HIP – solvent	2 x 18 ml
Dropping bottle of Reagent for HIP	1 pc

Information leaflet

Storage, expiration:

Reagents for the HIPPURATE test are stored in dark at the temperature of (+2 to +25)°C. Expiration is marked at each package.

The dropping bottle with the reagent for HIP is stored in dark at the temperature of (+2 to +8)°C, the prepared solution is usable up to 4 months since the date of preparation – not to exceed the expiry date of the set.

Warning:

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagents of the set contain highly flammable aceton (50 %, HIP – solvent) and flammable, harmful butan-1-ol (50 %, HIP – solvent) and harmful ninhydrin (100 %, HIP – weight).

Reagent for HIP – weight

Hazard statement:

H302 Harmful if swallowed.

H315 Causes skin irritation.

H319 Causes serious eye irritation.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Reagent for HIP – solvent

Danger

Hazard statement:

H225 Highly flammable liquid and vapour.

H302 Harmful if swallowed.

H315 Causes skin irritation.

H318 Causes serious eye damage.

H335 May cause respiratory irritation.

H336 May cause drowsiness or dizziness.

EUH066 Repeated exposure may cause skin dryness or cracking.

Precautionary statement:

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P261 Avoid breathing dust/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P308+P311 IF exposed or concerned: Call a POISON CENTER or doctor.

Reagent for HIP – classification see Reagent for HIP – solvent.

When preparing the solution of Reagent for HIP and when working with it, proceed with caution. Staining your skin with the reagent results in coloured blemishes, which are hardly to remove!

Waste disposal:

Contingent rest of the reagent should be liquidated in accordance with local and national regulations relating to the safe handling of such materials.

Put packaging paper waste and containers to recycling.

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