Reagent for NITRATEtest

Cat. No.	Product Name	Packaging (Content)
MLT00021	Reagent for NITRATEtest	Reagent for NIT: 1 × 18 mL, Reagent for NIT: Zinc powder: 1 × 10 g, Instruction for use

INTENDED USE

Reagent for NITRATE test is an *in vitro* diagnostic reagent used as supplementary colour forming reagent for some of MIKROLATEST® kits with pure bacterial cultures obtained from various clinical materials in the general population. The kits are dedicated for visual evaluation or for determination on automatic system ErbaExpert. The kit is intended as assistance in establishing a diagnosis or pathological state which can be caused by the tested pathogen. For professional use in clinical laboratory only.

PRINCIPLE

Reagent for NITRATE test is colour forming reagent for test for the demonstration of the proof of reduction of nitrate or nitrite (NIT, Nitrate, NO₃ or Nitrite, NO₂) and serves for colour expression of this test.

DESCRIPTION AND COMPOSITION

Reagent for NIT

Sulphanilic acid
N,N-dimethyl-1-naftylamin
HCl
Distilled water
Zinc powder

WORKING SOLUTION PREPARATION

Ready to use. Avoid to contamination.

STABILITY AND STORAGE

The unopened Reagent for NITRATE test is stable till the expiry date stated on the bottle and kit label when stored at (+2 to +8) °C. Store away from direct sunlight.

ASSAY PROCEDURE

The Reagent for NITRATE test should be used in accordance with the instructions for use of the specific kit.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagent

UFI: 3EK8-AJET-QE7Q-WSSF

Danger

Contains: acetic acid, sulphanilic acid

Hazard statement:

H314 Causes severe skin burns and eye damage.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.
P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
EUH208 Contains sulphanilic acid. May produce an allergic reaction.

Zinc powder

Warning

Hazard statement:

H410 Very toxic to aquatic life with long lasting effects.

Precautionary statement:

P260 Do not breathe dust.
P273 Avoid release to the environment.

WASTE MANAGEMENT

Please refer to local legal requirements.

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