Reagent for OXIDASEtest

Cat. No.	Product Name	Packaging (Content)
MLT00022	Reagent for OXIDASEtest	Reagent for OXI: 1 x 18 mL, Instruction for use

INTENDED USE

Reagent for OXIDASE test is an *in vitro* diagnostic reagent used as supplementary colour forming reagent for the diagnostic strips OXItest with pure bacterial cultures obtained from various clinical materials in the general population. The kits are dedicated for visual evaluation. The kit is intended as assistance in establishing a diagnosis or pathological state which can be caused by the tested pathogen. For professional use in clinical laboratory only.

PRINCIPLE

Reagent for OXIDASE test is colour forming reagent for test for the detection of bacterial cytochrome oxidase and serves for increasing of the strips OXItest sensitivity.

DESCRIPTION AND COMPOSITION

NaCl
Boric acid
Sodium tetraborate
Purified water

WORKING SOLUTION PREPARATION

Ready to use. Avoid to contamination.

STABILITY AND STORAGE

The unopened Reagent for the OXIDASE test is stable till the expiry date stated on the bottle and kit label when stored at (+2 to +25) °C. Store away from direct sunlight.

ASSAY PROCEDURE

The Reagent for OXIDASE test should be used in accordance with the instructions for use of the specific kit.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

UFI: VMK8-AJTM-AE7Q-7FXK

Danger

Contains: disodium tetraborate decahydrate

Hazard statement:

H360FD May damage fertility. May damage the unborn child.

Precautionary statement:

P263 Avoid contact during pregnancy and while nursing.
P280 Wear protective gloves.
P308 + P313 IF exposed or concerned: Get medical advice/attention.

WASTE MANAGEMENT

Please refer to local legal requirements.

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