URE-HPtest

Kat. č.: MLT00033 Pro mikrobiologii

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URE-HPtest

Cat. No.: MLT00033 For microbiology

URE-HPtest is intended for rapid detection of Helicobacter pylori on the basis of urease activity of tested sample. URE-HPtest is placed in divided microwell plate. One plate, containing 12 rows with 8 wells, enables up to 96 determinations.

Composition of the URE-HPtest kit:

- Divided microtitration plates with dessicant, 3 pieces
- Lid
- Storage bag
- 12 report sheets
- Instruction

Storage, expiration:

The URE-HPtest should be stored in a refrigerator at (+2 to +8) °C until the expiry date indicated on the package.

OPERATING PROCEDURE

Caution:

For professional use only.

Required material:

- Scalpel
- Sterile saline
- Paraffin oil, sterilized, Cat. No. MLT00042 – 750 determinations
- Desinfection

Observe the principles for working with infectious material!

Operating procedure:

- Open the aluminium sachet close to the weld and take out the plate
- Cut off required number of strips from plate
- Remove the adhesive tape from individual strips and insert them into prepared frame. In case, that you work with MIKROLATEST® kit for the first time and an empty frame is not available, use the frame of the first plate. The unused strips of the first plate put into the storage bag freely.
- Put the rest of the plate with dessicant to the storage bag and store it in a refrigerator for further use; protect it from humidity.
- Dispense 0.2 ml of sterile saline into prepared wells and cover the frame with a lid.
- Put remaining strips of microwells into refrigerator for further use.
- After about 10 min. it is prepared for testing.
- Leave one well as a negative control.
- Mark the position of microwell (e.g. well No. A1) and time when you inserted the sample into well.
- After inserting of sample overlay the well with two drops of paraffin oil; this prevents the influence of volatile ammonia released during urea hydrolysis in adjacent wells.
- Cover the frame with a lid and incubate at (+15 to +25) °C.
- Record the reaction after 5–10 min., 30–60 min. and 3–4 hours of incubation.
- Compare the reaction with a negative control and with a table "Interpretation of reaction".
- Record the results into report sheet.

Note:

Any uneven distribution of substrate in the well does not affect the functionality of the test.

Interpretation of reaction:

Reaction	Colour
positive	red, orange
negative	yellow

Performance:

The kit was tested on a set of 244 samples. The determination of 87 % was identical with the histological examination, 13 % of the determinations were different.

Disposal of used material:

- After use, all strips must be autoclaved or incinerated.
- Put paper packaging waste to recycling.

Health protection:

Components of the kit do not contain any dangerous substances.

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