Hb-Vario HbA1c Calibrator Set CM

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Lot: XXXXXXX
Expiry: MM/YYYY

ATTESTED VALUES

Reference system	HbA1c Level 1	HbA1c Level 2
NGSP % (National Glycohemoglobin Standardization Program)	0.0	00.0
IFCC % (International Federation of Clinical Chemistry)	0.0	00.0
IFCC mmol/mol (International Federation of Clinical Chemistry)	0.0	00.0

INTENDED USE

Hb-Vario HbA1c Calibrator Set CM is intended for calibration of glyca-
ted haemoglobin on Hb-Vario Analyser.

For *in vitro* diagnostic use only. These calibrators are intended for use
with instruments for the measurement of HbA1c and are designed to
allow calibration to either the American NGSP system or the Internati-
onal IFCC system.

The kit is intended for use with Hb-Vario Kit, Cat. No. REG00038 with
the same LOT.

PRINCIPLE

HbA1c measurement is used clinically to assess glycaemic control in dia-
betes. Calibrators should be used whenever a column or reagent kit is
replaced. The values assigned to these calibrators have been carefully
determined with reference to NGSP and IFCC reference materials. The-
se materials should be used in accordance with Hb-Vario user manual.

REAGENT COMPOSITION

The kit contains 2 chromatographic columns and 6 vials - 3 × 2 levels
of calibrators (0.5 ml). Each vial of calibrator contains a preparation of
freeze dried human whole blood.

STABILITY / STORAGE

Stored unopened between 2 and 8°C, these products are stable until
the expiration date stated on the product label (given that the column is
kept in its humid environment). After reconstitution, the calibrators are
stable for 1 day between 2 and 8°C.

REAGENT PREPARATION

The column is ready to use, just ensure that the humidifier is still wet
upon opening the container.

For calibrators: remove the stopper carefully. Using a volumetric pipette,
measure 0.5 ml of Rinse solution and add to the vial to reconstitute the
calibrator. Replace the stopper firmly and mix gently.

Allow to stand for 5 minutes. Visually inspect that all the material has
dissolved before use. Treat the reconstituted calibrator as whole blood
and use for calibration of the Hb-Vario.

Do not use time-expired material. Do not use incompletely dissolved
material. Do not use material that has not been stored in accordance
with the recommendations above.

WARNING AND PRECAUTIONS

All the individual blood donor lots used to prepare this material were
tested by FDA approved methods and found to be non-reactive for he-
patitis B surface antigen (HbsAG), hepatitis C (HCV), HIV-1 and HIV-2.
No test methods can offer complete assurance that these materials are
free from infective agents. Handle, use and dispose of these calibrators
as if they were human specimens and potentially infective.

Hazards identification in accordance with Regulation (EC) No 1272/2008 Level 1, 2

The reagents are not classified as dangerous.




FIRST AID

In case of an accidental ingestion, wash up the mouth and drink about 0.5
l of water. On eye contact rinse the eye quickly and thoroughly with the jet
of tap of water. Contaminated skin should be washed with warm water and
soap. In all serious cases of health damage consult a physician.

WASTE MANAGEMENT

Please refer to local legal requirements.

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