Hb-Vario HbA1c Calibrator Set

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Lot: XXXXXXXX
Expiry: MM/YYYY

ATTESTED VALUES

Reference system	HbA1c Level 1	HbA1c Level 2
NGSP % (National Glycohemoglobin Standardization Program)	0.0	00.0
IFCC % (International Federation of Clinical Chemistry)	0.0	00.0
IFCC mmol/mol (International Federation of Clinical Chemistry)	0.0	00.0

INTENDED USE

Hb-Vario HbA1c Calibrator Set is intended for calibration of glycated haemoglobin on Hb-Vario Analyser.
For *in vitro* diagnostic use only. These calibrators are intended for use with instruments for the measurement of HbA1c and are designed to allow calibration to either the American NGSP system or the International IFCC system. The kit is intended for use with Hb-Vario Kit, Cat. No. REG00038.

PRINCIPLE

HbA1c measurement is used clinically to assess glycaemic control in diabetes. Calibrators should be used whenever a column or reagent kit is replaced. The values assigned to these calibrators have been carefully determined with reference to NGSP and IFCC reference materials. These materials should be used in accordance with Hb-Vario user manual.

REAGENT COMPOSITION

The kit contains 4 vials -2 × 2 levels of calibrators (0.5 ml). Each vial of calibrator contains a preparation of freeze dried human whole blood.

STABILITY / STORAGE

Stored unopened between 2 and 8°C, these products are stable until the expiration date stated on the product label. After reconstitution, the calibrators are stable for 1 day between 2 and 8°C.

REAGENT PREPARATION

Remove the stopper carefully. Using a volumetric pipette, measure 0.5 ml of Rinse solution and add to the vial to reconstitute the calibrator. Replace the stopper firmly and mix gently. Allow to stand for 5 minutes. Visually inspect that all the material has dissolved before use. Treat the reconstituted calibrator as whole blood and use for calibration of the Hb-Vario. Do not use time-expired material. Do not use incompletely dissolved material. Do not use material that has not been stored in accordance with the recommendations above.

WARNING AND PRECAUTIONS

All the individual blood donor lots used to prepare this material were tested by FDA approved methods and found to be nonreactive for hepatitis B surface antigen (HbsAG), hepatitis C (HCV), HIV-1 and HIV-2. No test methods can offer complete assurance that these materials are free from infective agents. Handle, use and dispose of these calibrators as if they were human specimens and potentially infective.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Level 1, 2

The reagents are not classified as dangerous.

FIRST AID

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

WASTE MANAGEMENT

Please refer to local legal requirements.

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