# Hb-Vario HbA1c Control Set

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## INTENDED USE

For *in vitro* diagnostic use only. These controls are intended for use in a quality control procedure for determining imprecision and accuracy of glycated haemoglobin measurements. The values and ranges included herein are aligned to the American NGSP system and the International IFCC system.

The routine use of Control Material is recommended for the long term monitoring of imprecision of analytical methods and is a requirement of Good Laboratory Practice.

The kit is intended for use with Hb-Vario Kit, REG00038.

## REAGENT COMPOSITION

The kit contains 4 vials – 2 × 2 levels of control (0.5 ml). Each vial of control contains a preparation of freeze dried human whole blood.

## STABILITY / STORAGE

Stored unopened between 2 and 8 °C this product is stable until the expiration date. After reconstitution, the material is stable for 3 days at 2–8 °C and 3 weeks at -20 °C when stored as whole blood aliquots (100 µl aliquots for the whole blood mode or 25 µL aliquots for the lysate mode). Avoid multiple freeze-thaw cycles. In case of usage of the controls in the lysate mode, fresh lysate samples should be prepared daily from the whole blood aliquot for optimal performance.

## REAGENT PREPARATION

Remove the stopper carefully. Using a volumetric pipette, measure 0.5 ml of Rinse solution and add to the vial to reconstitute the control. Replace the stopper firmly and mix gently. Allow to stand for 5 minutes. Visually inspect that all the material has dissolved before use. Treat the reconstituted control as whole blood.

Do not use time-expired material. Do not use incompletely dissolved material. Do not use material that has not been stored in accordance with the recommendations above.

## WARNING AND PRECAUTIONS

All the individual blood donor lots used to prepare this material were tested by FDA approved methods and found to be non-reactive for hepatitis B surface antigen (HbsAG), hepatitis C (HCV), HIV-1 and HIV-2. No test methods can offer complete assurance that these materials are free from infective agents. Handle, use and dispose of these calibrators as if they were human specimens and potentially infective.

## Hazards identification in accordance with Regulation (EC) No 1272/2008 Level 1, 2

The reagents are not classified as dangerous.

## FIRST AID

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

## WASTE MANAGEMENT

Please refer to local legal requirements.

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