# ALBUMIN

Cat. No.	Pack Name	Packaging (Content)
BLT00001	ALB 250	R1: 5 × 50 ml, R2 standard: 1 × 3 ml
BLT00002	ALB 500	R1: 5 × 100 ml

## INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Albumin in human serum and plasma.

## CLINICAL SIGNIFICANCE

Albumin, a major plasma protein, is synthesised in the liver from amino acids which are absorbed from the ileum. It's functions include regulation of distribution of extra-cellular fluid, transportation of various hormones, vitamins and trace metals.

### Increased levels are observed in

Dehydration due to reduced plasma water content.

Stasis during venipuncture which causes fluid to escape into the extravascular compartment.

### Decreased levels are observed in

Excessive protein loss (mainly albumin) – from kidney, skin or intestine.

Decreased synthesis due to dietary, hepatic disease or malabsorption.

Increased catabolism in fever, untreated diabetes mellitus and hypertension.

## PRINCIPLE

Albumin binds with Bromo Cresol Green (BCG) at pH 4.2 causing a shift in absorbance of the yellow BCG dye. The blue-green colour formed is proportional to the concentration of albumin, when measured photometrically between 540–630 nm with maximum absorbance at 625 nm.

## REAGENT COMPOSITION

### R1

Bromocresol green 0.21 mmol/l

Succinate Buffer 100 mmol/l

Sodium Azide 0.5 g/l

**R2 standard** See bottle label

## REAGENT PREPARATION

Reagents are liquid, ready to use.

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

## SPECIMEN COLLECTION AND HANDLING

Use unheamolytic serum or plasma (EDTA, heparin)

It is recommended to follow NCCLS procedures (or similar standardized conditions).

**Stability in serum:** 1 month at 2–8°C  
1 week at 15–25°C  
at least 3 months at -20°C

Discard contaminated specimens.

## CALIBRATION

Calibration with the standard included in the kit or the calibrator XL MULTICAL 4x3 (Cat.No. XSYS0034) or XL MULTICAL 10x3 (Cat.No. XSYS0122) is recommended.

## QUALITY CONTROL

For quality control ERBA NORM 4x5 (Cat.No. BLT00080) or ERBA NORM 10x5 (Cat.No. XSYS0123) and ERBA PATH 4x5 (Cat.No. BLT00081) or ERBA PATH 10x5 (Cat.No. XSYS0124) are recommended.

## UNIT CONVERSION

g/dl × 10 = g/l

## EXPECTED VALUES 2

Serum:

0–4 d 2.8–4.4 g/dl

4 d – 14 y 3.8–5.4 g/dl

14–18 y 3.2–4.5 g/dl

Adult

20–60 y 3.5–5.2 g/dl

60–90 y 3.2–4.6 g/dl

**It is recommended that each laboratory verify this range or derives reference interval for the population it serves.**

## PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

**Limit of quantification:** 0.1 g/dl

**Linearity:** 7.2 g/dl

**Measuring range:** 0.1–7.2 g/dl

## PRECISION

Intra-assay precision Within run (n=20)	Mean (g/dl)	SD (g/dl)	CV (%)
Sample 1	3.303	0.017	0.51
Sample 2	4.926	0.044	0.88

Inter-assay precision Run to run (n=20)	Mean (g/dl)	SD (g/dl)	CV (%)
Sample 1	3.345	0.043	1.29
Sample 2	4.612	0.043	0.93

## COMPARISON

A comparison between XL-Systems Albumin (y) and a commercially available test (x) using 40 samples gave following results:

y = 1.000 × + 0.065 g/dl

r = 0.997

## INTERFERENCES

Following substances do not interfere:

haemoglobin up to 10 g/l, bilirubin up to 40 mg/dl, triglycerides up to 2000 mg/dl.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

### Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagent 1 contains < 1.2 % succinic acid.

R1 UFI: Y8RN-RJU5-0E7S-53CM

**Warning**

### Hazard statement:

H319 Causes serious eye irritation.

### Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rising.  
P337+P313 If eye irritation persists: Get medical advice/attention.

## WASTE MANAGEMENT

Please refer to local legal requirements.

## ASSAY PROCEDURE

**Wavelength** 578 (540–630) nm  
**Cuvette** 1 cm

	Reagent blank	Standard (Cal.)	Sample
Reagent 1	1.00 ml	1.00 ml	1.00 ml
Sample	–	–	0.01 ml
Standard (Cal.)	–	0.01 ml	–
Distilled water	0.01 ml	–	–

Mix and incubate 1–5 min. at 37 °C. Measure absorbance of the sample  $A_{sam}$  and standard  $A_{st}$  against reagent blank.

## CALCULATION

Albumin (g/dl) =  $\frac{\Delta A_{sam}}{\Delta A_{st}} \times C_{st}$   $C_{st}$  = standard (calibrator) concentration

## ASSAY PARAMETERS FOR PHOTOMETERS

Mode	End point
Wavelength 1 (nm)	600
Sample Volume (μl)	5/10
Reagent Volume (μl)	500/1000
Incubation time (min.)	1
Incubation temp. (°C)	37
Normal Low (g/dl)	3.5
Normal High (g/dl)	5.2
Linearity Low (g/dl)	0
Linearity High (g/dl)	7.2
Concentration of Standard	See bottle label
Blank with	Reagent
Absorbance limit (max.)	0.4
Units	g/dl

Applications for automatic analysers are available on request.

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