# ALKALINE PHOSPHATASE MEG

Cat. No.	Pack Name	Packing (Content)
BLT00005	ALP MEG 500	R1: 4 x 100 ml, R2: 1 x 100 ml

EN

## INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of ALP in human serum and plasma.

## CLINICAL SIGNIFICANCE

Human ALP consists of a group of enzymes which hydrolyse phosphates at an alkaline pH. ALP is found in practically all tissues of the body but in high concentrations in the osteoblasts

of bone, liver, placenta, kidney, intestinal wall and lactating mammary glands. In adults the ALP normally found circulating in the serum is largely derived from the liver. In children or in adolescents going through pubertal growth spurts, there is an additional contribution from bone and this accounts for the higher reference interval for these groups. Pregnancy also raises the normal values of ALP.

Raised ALP levels are often observed in bone disease or liver disease involving the biliary tract. If the source of the isoenzyme is not apparent then estimation of GGT may help differentiate between the two. A raised GGT in the presence of a raised ALP would suggest the liver is the primary source.

Increased ALP (usually normal GGT) is seen in Osteomalacia and Rickets, primary hyperparathyroidism with bone involvement, Pagets disease, secondary carcinoma in bone and some cases of osteogenic sarcoma. Increased levels of ALP (usually with a raised GGT) is seen in cholestasis, hepatitis, cirrhosis, space occupying lesions and malignancy with bone or liver involvement or direct production. Low levels of ALP may be observed in conditions which cause arrested bone growth or in hypophosphatasia.

## PRINCIPLE

Under alkaline conditions, alkaline phosphatase (alkaline phosphohydrolase of ortho-phosphoric monoesters, EC 3.1.3.1. – ALP), splits 4-nitrophenylphosphate into 4-nitrophenol and phosphate.

The measure of the catalytic concentration of the enzyme is the amount of released 4-nitrophenol, which is evaluated photometrically in a kinetic way.

## REAGENT COMPOSITION

### R1

N-methyl-D-glucamine, pH 10.1 (37°C)	625 mmol/l
Sodium chloride	71.7 mmol/l
Magnesium chloride	0.64 mmol/l

### R2

4-Nitrophenylphosphate	102 mmol/l
------------------------	------------

## COMPOSITION OF MIXTURE

N-methyl-D-glucamine, pH 10.1 (37°C)	500 mmol/l
Sodium Chloride	55.8 mmol/l
Magnesium chloride	0.50 mmol/l
4-Nitrophenylphosphate	20 mmol/l

## REAGENT PREPARATION

Reagents are liquid, ready to use.

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

### Two reagents method – substrate start

Reagents are ready to use. After opening, reagents are stable until expiry date at 2–8°C if stored at appropriate conditions, closed carefully and without any contamination.

## Monoreagent method – sample start

Mix 4 portion of reagent R1 with 1 portion of reagent R2.

Stability: 1 day at 15–25°C in dark  
1 week at 2–8°C in dark

Maximum allowable absorbance of the working reagent measured at 420 nm against distilled water is 1.0

## SPECIMEN COLLECTION AND HANDLING

Use serum, plasma (heparin, EDTA).

It is recommended to follow NCCLS procedures (or similar standardized conditions).

**Stability in serum / plasma:** 4 hours at 20–25°C  
3 days at 4–8°C  
2 months at -20°C

Discard contaminated specimens.

## CALIBRATION

Calibration with calibrator Lyonorm CALIBRATOR, Cat. No. BLT00069 is recommended.

## QUALITY CONTROL

For quality control, it is recommended to use following materials:

Lyonorm HUM N, Cat. No. BLT00070 and Lyonorm HUM P, Cat. No. BLT00071.

## UNIT CONVERSION

U/l x 0.017 = µkat/l

## EXPECTED VALUES

fS ALP (µkat/l) 37 °C

Men	0.73–2.60
Women	0.62–2.40
Children (males up 14 years, females up 12 years)	2.35–8.00
Newborns	1.90–7.00

**It is recommended that each laboratory verify this range or derives reference interval for the population it serves.**

## PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

**Limit of quantification:** 0.17 µkat/l

**Linearity:** 25 µkat/l

**Measuring range:** 0.17–25 µkat/l

Intra-assay precision Within run (n=20)	Mean (µkat/l)	SD (µkat/l)	CV (%)
Sample 1	2.29	0.01	0.93
Sample 2	5.78	0.02	0.39

Inter-assay precision Run to run (n=20)	Mean (µkat/l)	SD (µkat/l)	CV (%)
Sample 1	2.27	0.03	1.13
Sample 2	5.76	0.04	0.71

## COMPARISON

A comparison between ALP MEG 500 (y) and a commercially available test (x) using 99 samples gave following results:

y = 1.035 x - 0.049 µkat/l

r = 1.000

## INTERFERENCES

Following substances do not interfere:

Bilirubin up to 1000 µmol/l, triglycerides up to 20 mmol/l, haemoglobin up to 7g/l.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Reagents of the kit are not classified as dangerous.

## FIRST AID

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

## WASTE DISPOSAL

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials.

Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

## PROCEDURE

Wavelength: 420 (405–430) nm

Cuvette: 1 cm

Temperature: 37 °C

Serum/reaction mixture ratio 1/51

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio.

### Two reagents method – substrate start

	Reagent blank	Calibrator	Sample
Reagent 1	0.800 ml	0.800 ml	0.800 ml
Sample	–	–	0.020 ml
Calibrator	–	0.020 ml	–
Distilled water	0.020 ml	–	–
Mix and after 5 min. incubation (at 37 °C ) add:			
Reagent 2	0.200 ml	0.200 ml	0.200 ml

Mix, incubate 1 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA/min).

### Monoreagent method – sample start

	Reagent blank	Calibrator	Sample
Working reagent	1.000 ml	1.000 ml	1.000 ml
Sample	–	–	0.020 ml
Calibrator	–	0.020 ml	–
Distilled water	0.020 ml	–	–

Mix, incubate 1 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA/min).

## CALCULATION

$$\text{Alkaline phosphatase (µkat/l)} = \frac{\Delta A_{\text{sample}}/\text{min}}{\Delta A_{\text{cal}}/\text{min}} \times C_{\text{cal}}$$

C<sub>cal</sub> = calibrator concentration

**Applications for automatic analysers will be supplies on request.**

 Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ  
e-mail: diagnostics@erbamannheim.com, www.erbamannheim.com

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