ALPHA AMYLASE (SINGLE REAGENT)

Cat. No.	Pack Name	Packaging (Content)
BLT00006	AMY SINGLE 100	R1: 5 x 20 ml

INTENDED USE


Diagnostic reagent for quantitative *in vitro* determination of alpha-Amylase in human serum and plasma.

CLINICAL SIGNIFICANCE

α -Amylase is derived mainly from the salivary glands and the exocrine pancreas. α -Amylase catalyses the hydrolysis of α -1-4 glucosidic linkages of starch and other related polysaccharides to produce maltose and other oligosaccharides. The enzyme is a relatively small molecule which is rapidly cleared by the kidneys and excreted in the urine. α -Amylase is most frequently measured in the diagnosis of acute pancreatitis when serum levels may be grossly elevated. In acute pancreatitis α -amylase starts to rise approximately 4 hours after the onset of pain, reaches a peak at 24 hours and remains elevated for 3–7 days. Hyperamylasemia is also associated with other acute abdominal disorders, biliary dysfunction, salivary gland disorders, ruptured ectopic pregnancy and macroamylasemia.

PRINCIPLE

2-Chloro-4-nitrophenol- β -1-4galactopyranosylmaltotrioxide (CNP-G) is a direct substrate for determination of α -amylase activity, which does not require the presence of ancillary enzymes. The rate of 2-chloro-4-nitrophenol formation can be monitored at (400–420) nm and is proportional to the α -amylase activity.

REAGENT COMPOSITION

R1

MES buffer	50 mmol/l
Calcium Chloride	3.81 mmol/l
Sodium Chloride	300 mmol/l
Potassium Thiocyanate	450 mmol/l
Sodium Azide	13.85 mmol/l
CNPG	0.91 mmol/l

REAGENT PREPARATION

Reagent is liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

SPECIMEN COLLECTION AND HANDLING

Use serum, plasma (heparin, EDTA), urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability

in serum/plasma:	7 days	at 20–25°C
	7 days	at 4–8°C
	1 year	at -20°C
in urine:	2 days	at 20–25°C
	10 days	at 4–8°C
	3 weeks	at -20°C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL 4x3, Cat. No. XSYS0034 or XL MULTICAL 10x3, Cat. No. XSYS0122 is recommended.

QUALITY CONTROL

For quality control ERBA NORM 4x5, Cat. No. BLT00080 or ERBA NORM 10x5,

Cat. No. XSYS0123 and ERBA PATH 4x5, Cat. No. BLT00081 or ERBA PATH 10x5, Cat. No. XSYS0124 are recommended.

UNIT CONVERSION

U/l x 0.017 = μ kat/l

EXPECTED VALUES ⁷

at 37°C

Serum: up to 80 U/l

Urine: up to 500 U/l

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 10.8 U/l

Linearity: 1500 U/l

Measuring range: 10.8–1500 U/l

PRECISION

Intra-assay precision Within run (n=20)	Mean (U/l)	SD (U/l)	CV (%)
Sample 1	247.1	2.5	1.0
Sample 2	260.8	2.6	1.0

Inter-assay precision Run to run (n=20)	Mean (U/l)	SD (U/l)	CV (%)
Sample 1	58.3	1.3	2.2
Sample 2	142.9	2.3	1.6

COMPARISON

A comparison between XL-Systems Amylase (y) and a commercially available test (x) using 40 samples gave following results:

y = 0.973 x - 4.80 U/l

r = 0.989

INTERFERENCES

Following substances do not interfere:

haemoglobin up to 2.5 g/l, bilirubin up to 40 mg/dl, triglycerides up to 2000 mg/dl.

Note:

Saliva and skin contain alpha-amylase therefore never pipette reagents by mouth and avoid contamination of samples and reagents. Even trace contamination can affect results.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Hazards identification in accordance with Regulation (EC) No 1272/2008

UFI: 6RYE-AE3K-3H7A-Q7QR

Danger

Contains: potassium thiocyanate

Hazard statement:

H318 Causes serious eye damage.

Precautionary statement:

P280 Wear eye protection.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing. P337 + P313 If eye irritation persists: Get medical advice/attention.

WASTE MANAGEMENT

Please refer to local legal requirements.

ASSAY PROCEDURE

Wavelength: 405 (400–420) nm

Cuvette: 1 cm

Reagent 1	1000 μ l
Sample	20 μ l

Mix, incubate 1 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA /min).

CALCULATION

$$1. \text{ Amylase activity (U/l)} = \frac{\Delta A_{\text{sam}} / \text{min.}}{\Delta A_{\text{cal}} / \text{min.}} \times C_{\text{cal}}$$

C_{cal} = calibrator concentration

2. Using factor:

Amylase activity (U/l) = f x ΔA /min

f = factor

f = 3128 (at 405 nm)

Applications for automatic analysers are available on request.

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	Kinetic
Wavelength 1 (nm)	405
Sample Volume (μ l)	10/20
Working Reagent Volume (μ l)	500/1000
Lag time (sec.)	60
Kinetic interval (sec.)	60
No. of readings	3
Kinetic factor	3128
Reaction temperature (°C)	37
Reaction direction	Increasing
Normal Low (U/l)	0
Normal High (U/l)	80
Linearity Low (U/l)	10.8
Linearity High (U/l)	1500
Blank with	Water
Absorbance limit (max.)	0.14
Units	U/l

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