# HDL CHOLESTEROL DIRECT

Cat. No.	Pack Name	Packing (Content)
BLT00030	HDL C DIRECT 240	R1: 3 x 60 ml, R2: 3 x 20 ml

EN

## INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of HDL Cholesterol in human serum and plasma.

## CLINICAL SIGNIFICANCE

High-density lipoproteins (HDL) compose one of the major classes of plasma lipoproteins. They are synthesized in liver as complexes of apolipoprotein and phospholipid and are capable of picking up cholesterol and carrying it from arteries to the liver, where the cholesterol is converted to bile acids and excreted into the intestine.

An inverse relationship between HDL Cholesterol (HDL-C) levels in serum and the incidence/prevalence of coronary heart disease (CHD) has been demonstrated in a number of epidemiological studies. The importance of HDL-C as a risk factor for CHD is now recognized.

Accurate measurement of HDL-C is of vital importance when assessing patient's risk for CHD.

## PRINCIPLE

Antihuman  $\beta$ -lipoproteins antibody in the Reagent 1 binds to lipoproteins (LDL, VLDL and chylomicrons) other than HDL. The antigen-antibody complexes formed blocked enzyme reaction when R2 is added. Cholesterol esterase (CHES) and cholesterol oxidase (CHOD) in R2 react only with HDL cholesterol. Hydrogen peroxide produced by the enzyme reactions with HDL cholesterol yields a blue color complex upon oxidative condensation of F-DAOS [N-ethyl-N-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxy-4-fluoroaniline, sodium salt] and 4-aminoantipyrine (4-AA) in the presence of peroxidase (POD). The absorbance of the blue color complex (substituted chinondimine) is measured at optimum wavelength 593 nm.

LDL, VLDL, chylomicrons  $\xrightarrow[\beta\text{-lipoprotein antibodies}]{\text{antihuman}}$  antigen-antibody complex

HDL cholesterol + H<sub>2</sub>O + O<sub>2</sub>  $\xrightarrow{\text{CHES\&CHOD}}$   $\Delta^4$  - cholestenone + fatty acid + H<sub>2</sub>O<sub>2</sub>

H<sub>2</sub>O<sub>2</sub> + F-DAOS + 4-aminoantipyrine  $\xrightarrow{\text{POD}}$  blue dye + H<sub>2</sub>O

## REAGENT COMPOSITION

**R1**  
Good's buffer, pH 7.0 30 mmol/l  
4-aminoantipyrine (4-AA) 0.9 mmol/l  
Peroxidase (POD) 40  $\mu$ kat/l  
Ascorbate oxidase (AOD) 45  $\mu$ kat/l  
Antihuman  $\beta$ -lipoproteins antibody

**R2**  
Good's buffer, pH 7.0 30 mmol/l  
Cholesterol esterase (CHES) 66.6  $\mu$ kat/l  
Cholesterol oxidase (CHOD) 333.3  $\mu$ kat/l  
F-DAOS 0.8 mmol/l

## REACTION MIXTURE

Good's buffer, pH 7.0 30 mmol/l  
4-aminoantipyrine (4-AA) 0.67 mmol/l  
Peroxidase (POD) 29.75  $\mu$ kat/l  
Ascorbate oxidase (AOD) 33.47  $\mu$ kat/l  
Cholesterol esterase (CHES) 16.51  $\mu$ kat/l  
Cholesterol oxidase (CHOD) 82.65  $\mu$ kat/l  
F-DAOS 0.198 mmol/l  
Antihuman  $\beta$ -lipoproteins antibody

## REAGENT PREPARATION

Reagents R1 and R2 are liquid, ready to use.

## STORAGE AND STABILITY

If stored at 2–8 °C, reagents are stable until expiry date, that is stated on the package. After first opening, reagents R1 and R2 are stable for 60 days at 2–8 °C in the dark and without contamination.

## SPECIMEN COLLECTION & HANDLING

Use serum or heparin plasma.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability in serum / plasma:

1 day at 20–25 °C

7 days at 4–8 °C

12 weeks at -20 °C

Discard contaminated specimens.

## CALIBRATION

For the calibration, it is recommended to use HDL Cholesterol Calibrator, Cat. No. BLT00072.

## QUALITY CONTROL

For quality control ERBA NORM 4x5, Cat. No. BLT00080 or ERBA NORM 10x5, Cat. No. XSYS0123 and ERBA PATH 4x5, Cat. No. BLT00081 or ERBA PATH 10x5, Cat. No. XSYS0124 are recommended.

## UNIT CONVERSION

mg/dl x 0,026 = mmol/l

## EXPECTED VALUES

fS HDL cholesterol (mmol/l) male 0.91 – 2.05

fS HDL cholesterol (mmol/l) female 1.09 – 2.28

**It is recommended that each laboratory verify this range or derives reference interval for the population it serves.**

## PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

**Limit of quantification:** 0.028 mmol/l

**Linearity:** 5 mmol/l

**Measuring range:** 0.028 – 5 mmol/l

## PRECISION

Intra-assay (n=20)	Mean (mmol/l)	SD (mmol/l)	CV (%)
Sample 1	0.795	0.009	1.12
Sample 2	1.218	0.019	1.54

Inter-assay (n=20)	Mean (mmol/l)	SD (mmol/l)	CV (%)
Sample 1	0.767	0.017	2.25
Sample 2	1.225	0.029	2.39

## COMPARISON

A comparison between HDL CHOLESTEROL DIRECT (y) and a commercially available test (x) using 40 samples gave following results:

N = 40 r = 0,996 y = 1.032 x + 0.004 mmol/l

## INTERFERENCES

Following substances do not interfere:

bilirubin up to 40 mg/dl, triglycerides up to 1200 mg/dl, haemoglobin up to 10 g/l.

Interference by N-acetylcysteine (NAC), acetaminophen and metamizole causes falsely low results. To carry out the test, blood withdrawal should be performed prior to administration of drugs.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

## Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagent R1

Warning

Contains: reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one

**Hazard statement:**

H317 May cause an allergic skin reaction.

**Precautionary statement:**

P280: Wear protective gloves/protective clothing/eye protection.

P303+P361+P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P333+P313: If skin irritation or rash occurs: Get medical advice/attention.

UFI: PAWG-4WA6-3J5S-PM0Q

Reagent R2

EUH208 Contains reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one. May produce an allergic reaction.

## FIRST AID

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

## WASTE DISPOSAL

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials. Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

## PROCEDURE

**Wavelength:** 600/700 nm

**Cuvette:** 1 cm

**Temperature:** 37 °C

Serum / reaction mixture ratio 1/121

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio.

	Reagent blank	Calibrator	Sample
Reagent R1	0.900 ml	0.900 ml	0.900 ml
Sample	-	-	0.010 ml
Calibrator	-	0.010 ml	-
Distilled water	0.010 ml	-	-
Mix and incubate 5 min. at 37 °C. Measure absorbance A <sub>i</sub> of the sample (A <sub>sam</sub> ), calibrator (A <sub>cal</sub> ) and blank (A <sub>bl</sub> ). Then add:			
Reagent R2	0.300 ml	0.300 ml	0.300 ml

Mix and incubate 5 min. at 37 °C. Measure absorbance A<sub>2</sub> of the sample (A<sub>sam</sub>), calibrator (A<sub>cal</sub>) and blank (A<sub>bl</sub>).

Calculate absorbance change  $\Delta A = A_2 - A_1$  for sample, calibrator and blank.

## CALCULATION

HDL Cholesterol (mmol/l) =  $\frac{\Delta A_{\text{sam}} - \Delta A_{\text{bl}}}{\Delta A_{\text{cal}} - \Delta A_{\text{bl}}} \times C_{\text{cal}}$  C<sub>cal</sub> = calibrator concentration

Applications for automatic analysers will be supplied on request.

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ  
e-mail: diagnostics@erba.com, www.erbamannheim.com

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