LACTATE DEHYDROGENASE-P

| Cat. No. | Pack Name | Packaging (Content) |
|----------|-----------|------------------------------|
| BLT00037 | LDH 100 | R1: 4 x 20 ml, R2: 1 x 20 ml |

EN

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of LDH in human serum and plasma (DGKCH method).

CLINICAL SIGNIFICANCE

The enzyme Lactate dehydrogenase (LDH) is concentrated in heart, kidney, liver, muscle and body tissues. Consequently, damage to these results in increased serum levels of LDH. Elevated levels are associated with myocardial infarction, renal damage, hepatitis, anemias, malignancies and muscular disease or damage.

PRINCIPLE

The LDH method is based on the recommendations of DGKCH (from pyruvate). This reagent uses pyruvate and is based on the method of Henry et al.

LDH catalyses the reduction of pyruvate to lactate oxidising reduced nicotinamide adenine dinucleotide (NADH) to NAD. The activity of LDH can be determined by the rate of decrease in absorbance at 340 nm as NAD⁺ is produced.

REAGENT COMPOSITION

| | |
|----------------------|-------------|
| R1 | |
| Tris Buffer (pH 7.5) | 100 mmol/l |
| Pyruvate | 2.0 mmol/l |
| R2 | |
| NADH | 1.66 mmol/l |

REAGENT PREPARATION

Reagents are liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

Two reagents method

Reagents are ready to use.

After opening, reagents are stable until expiry date at 2–8 °C if stored at appropriate conditions, closed carefully and without any contamination.

Monoreagent method

Mix 4 portion of Reagent R1 and 1 portion of Reagent R2.

Stability:

| | |
|----------|---------------------|
| 24 hours | at 15–25 °C at dark |
| 5 days | at 2–8 °C at dark |

SPECIMEN COLLECTION AND HANDLING

Use serum, plasma (heparin, EDTA).

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Loss of activity:

| | | |
|-----------------|-------------|-------|
| within 24 hours | at 15–25 °C | < 2 % |
| within 3 days | at 2–8 °C | < 8 % |

Stability at least 6 weeks at –20 °C

CALIBRATION

Calibration with calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

UNIT CONVERSION

U/l x 0.017 = µkat/l

EXPECTED VALUES ²

At 37 °C: 225 - 450 U/l

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 43.8 U/l

Linearity: 1200 U/l

Measuring range: 43.8–1200 U/l

| Intra-assay precision Within run (n=20) | Mean (U/l) | SD (U/l) | CV (%) |
|--|---------------|-------------|-----------|
| Sample 1 | 767.4 | 3.6 | 0.49 |
| Sample 2 | 760.2 | 7.8 | 0.99 |

| Inter-assay precision Run to run (n=20) | Mean (U/l) | SD (U/l) | CV (%) |
|--|---------------|-------------|-----------|
| Sample 1 | 562.8 | 13.8 | 2.43 |
| Sample 2 | 312.0 | 5.4 | 1.88 |

COMPARISON

A comparison between XL-Systems LDH (y) and a commercially available test (x) using 40 samples gave following results:

y = 1.982 x + 0.06 U/l

r = 0.996

INTERFERENCES

Following substances do not interfere:

bilirubin up to 20 mg/dl, triglycerides up to 500 mg/dl, haemoglobin up to 5.0 g/l.

Significant hemolysis may increase LD concentration because of high levels of LD in the erythrocytes.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagent R2 of the kit is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

ASSAY PROCEDURE

Wavelength 340 nm, Hg 334 nm, Hg 365 nm

Cuvette 1 cm

Two reagents method – substrate start

| | Reagent blank | Calibrator | Sample |
|-----------------|---------------|------------|----------|
| Reagent 1 | 0.800 ml | 0.800 ml | 0.800 ml |
| Sample | – | – | 0.020 ml |
| Calibrator | – | 0.020 ml | – |
| Distilled water | 0.020 ml | – | – |

Mix and after 1 min. incubation (at 37 °C) add:

| | | | |
|-----------|----------|----------|----------|
| Reagent 2 | 0.200 ml | 0.200 ml | 0.200 ml |
|-----------|----------|----------|----------|

Mix, incubate 1 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA/min).

Monoreagent method – sample start

| | Reagent blank | Calibrator | Sample |
|-----------------|---------------|------------|----------|
| Working reagent | 1.000 ml | 1.000 ml | 1.000 ml |
| Sample | – | – | 0.020 ml |
| Calibrator | – | 0.020 ml | – |
| Distilled water | 0.020 ml | – | – |

Mix, incubate 1 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA/min).

CALCULATION

$$1. \text{LDH (U/l)} = \frac{\Delta A_{\text{sam}}/\text{min}}{\Delta A_{\text{cal}}/\text{min}} \times C_{\text{ca}} \quad C_{\text{cal}} = \text{calibrator concentration}$$

2. Using factor:

LDH (U/l) = f x ΔA/min

f = factor

f = 8095 (at 340 nm)

Applications for automatic analysers are available on request.

ASSAY PARAMETERS FOR PHOTOMETERS

| | |
|---------------------------|------------|
| Mode | Kinetic |
| Wavelength 1 (nm) | 340 |
| Sample Volume (µl) | 10/20 |
| Reagent Volume (µl) | 500/1000 |
| Lag Time (sec.) | 60 |
| Kinetic Interval (sec.) | 60 |
| No. of readings | 3 |
| Kinetic factor | 8095 |
| Reaction temperature (°C) | 37 |
| Reaction direction | Decreasing |
| Normal Low (U/l) | 225 |
| Normal High (U/l) | 450 |
| Linearity Low (U/l) | 43.8 |
| Linearity High (U/l) | 1200 |
| Absorbance Limit (Max.) | 0.8 |
| Blank with | Water |
| Units | U/l |

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