LACTATEDEHYDROGENASE-L

Cat. No.	Pack Name	Packaging (Content)
BLT00038	LDH-L 100	R1: 4 x 20 ml, R2: 1 x 20 ml

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of catalytic concentration of lactate dehydrogenase in serum or plasma, IFCC method.

CLINICAL SIGNIFICANCE

The enzyme Lactate dehydrogenase (LDH) is concentrated in heart, kidney, liver, muscle and body tissues. Consequently, damage to these results in increased serum levels of LDH. Elevated levels are associated with myocardial infarction, re-nal damage, hepatitis, anemias, malignancies and muscular disease or damage.

PRINCIPLE

Lactate dehydrogenase catalyses the conversion of L-lactate to pyruvate with the simultaneous reduction of NAD⁺ to NADH. The change in absorbance with the time due to the conversion of NAD to NADH is directly proportional to LD activity.

REAGENT COMPOSITION

R1	
N-Methyl D-glucamin, pH 9.40	406 mmol/l
L-lactate	62.5 mmol/l
R2	
NAD ⁺	50 mmol/l

REACTION MIXTURE

N-Methyl D-glucamin, pH 9.40	325 mmol/l
L-lactate	50 mmol/l
NAD ⁺	10 mmol/l

REAGENT PREPARATION

Reagents are liquid, ready for use.

STABILITY AND STORAGE

Two reagents method

Reagents R1 and R2 are liquid, ready to use.

If stored at 2–8°C, reagents are stable until expiry date, that is stated on the package.

After opening, reagents are stable until expiry date at 2–8°C if stored at appropriate conditions, closed carefully and without any contamination.

Protect from light, particularly reagent 2!

Monoreagent method – sample start

Mix 4 portion of reagent R1 with 1 portion of reagent R2.

Stability: 4 hours at 15–25°C in the dark

24 hours at 2–8°C in the dark

Protect working reagent from light!

The unopened reagents are stable till the expiry date stated on the on the package.

SPECIMEN COLLECTION AND HANDLING

Use serum, plasma (heparin, EDTA).

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Loss of activity:

within 24 hours at 15–25°C < 2%

within 3 days at 2–8°C < 8%

Stability at least 6 weeks at –20 °C

CALIBRATION

Calibration with LYONORM CALIBRATOR, Cat. No. BLT00069 is recommended.

QUALITY CONTROL

For quality control LYONORM HUM N, Cat. No. BLT00070 and LYONORM HUM P, Cat. No. BLT00071 are recommended.

UNIT CONVERSION

U/l x 0.017 = µkat/l

EXPECTED VALUES

fS LDH-L (µkat/l) 37 °C

Male < 4.13

Female < 4.12

The range of reference values is only approximate; it is recommended that each laboratory verify the extent of the reference interval for their particular examined population.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 0.28 µkat/l

Linearity: 20 µkat/l

Measuring range: 0.28 – 20 µkat/l

PRECISION

Intra-assay (n=20)	Mean (µkat/l)	SD (µkat/l)	CV (%)
Sample 1	2.50	0.036	1.43
Sample 2	4.64	0.074	1.60

Inter-assay (n=20)	Mean (µkat/l)	SD (µkat/l)	CV (%)
Sample 1	2.33	0.081	3.46
Sample 2	4.52	0.094	2.09

COMPARISON

A comparison between XL-Systems LDH-L (y) and a commercially available test (x) using 40 samples gave following results

N = 40

r = 0.996

y = 0.0984 x 0,117 µkat/l

INTERFERENCES

Following substances do not interfere:

haemoglobin up to 5 g/l, bilirubin up to 40 mg/dl, triglycerides up to 1750 mg/dl.

Significant hemolysis may increase LD concentration because of high levels of LD in the erythrocytes.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagents of the kit are not classified as dangerous.

FIRST AID

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

WASTE DISPOSAL

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials. Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

ASSAY PROCEDURE

Wavelength: 340, 334, 365 nm

Cuvette: 1 cm

Temperature: 37°C

Serum/reaction mixture ratio 1/51

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio.

Two reagents method – substrate start

	Reagent blank	Calibrator	Sample
Reagent 1	0.800 ml	0.800 ml	0.800 ml
Sample	-	-	0.020 ml
Calibrator	-	0.020 ml	-
Distilled water	0.020 ml	-	-
Mix and incubate for 1-5 minutes at 37 °C. Then add:			
Reagent 2	0.200 ml	0.200 ml	0.200 ml

Mix and incubate for 1 minute at 37 °C and measure the initial absorbance exactly after 1 min and start a stopwatch. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA).

Monoreagent method – sample start

	Reagent blank	Calibrator	Sample
Working reagent	1.00 ml	1.00 ml	1.00 ml
Sample	-	-	0.020 ml
Calibrator	-	0.020 ml	-
0.9 % NaCl	0.020 ml	-	-

Mix and incubate for 1 minute at 37 °C and measure the initial absorbance exactly after 1 min and start a stopwatch. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA).

CALCULATION

$$1. \text{ LDH } (\mu\text{kat/l}) = \frac{\Delta A_{\text{sam}} - \Delta A_{\text{bl}}}{\Delta A_{\text{cal}} - \Delta A_{\text{bl}}} \times C_{\text{cal}} \quad C_{\text{cal}} = \text{calibrator concentration}$$

2. Using factor:

LDH (µkat/l) = f x ΔA/min.

f = factor:

Wavelength	Sample start	Substrate start
334 nm	137.5	171.25
340 nm	134.9	168.00
365 nm	250.0	311.25

NOTES

1. If ΔA/min > 0.150 at 334 (340) nm or 0.080 at 365 nm, dilute the sample in 1+10 ration with saline (0.9 % NaCl). Multiply the results by 11.

2. The reagents contains sodium azide (0.95 g/l) as preservative. Do not swallow! Avoid contamination with skin, and mucous membranes.

Applications for automatic analysers are available on request.

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1. Thomas L. Clinical laboratory diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. 89-94.
2. Moss DW, Henderson AR. Clinical enzymology In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia: W.B. Saunders Company, 1999. 617-721.
3. Schumann G, Bonora R, Ceriotti F, Féraud G et al. IFCC primary reference procedure for the measurement of catalytic concentrations of enzymes at 37 °C. Part 3: Reference procedure for the measurement of catalytic concentration of lactate dehydrogenase. Clin. Chem. Lab Med 2002, 40:643-48.

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