# LDL CHOLESTEROL DIRECT

| Cat. No. | Pack Name        | Packing (Content)            |
|----------|------------------|------------------------------|
| BLT00043 | LDL C DIRECT 80  | R1: 1 x 60 ml, R2: 1 x 20 ml |
| BLT00042 | LDL C DIRECT 240 | R1: 3 x 60 ml, R2: 3 x 20 ml |

EN

## INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of LDL Cholesterol in human serum and plasma.

## CLINICAL SIGNIFICANCE

Low Density Lipoproteins (LDL) are synthesized in the liver by the action of various lipolytic enzymes on triglyceride-rich Very Low Density Lipoproteins (VLDLs). Specific LDL receptors exist to facilitate the elimination of LDL from plasma by liver parenchymal cells. It has been shown that most of the cholesterol stored in atherosclerotic plaques originates from LDL. For this reason the LDL Cholesterol concentration is considered to be the most important clinical predictor, of all single parameters, with respect to coronary atherosclerosis.

Accurate measurement of LDL Cholesterol is of vital importance in therapies which focus on lipid reduction to prevent atherosclerosis or reduce its progress and to avoid plaque rupture. Can be applied on automated analyzers.

## PRINCIPLE

When a sample is mixed with Reagent 1, the protective reagent binds to LDL and protects LDL from enzyme reactions. Cholesterol esterase (CHES) and cholesterol oxidase (CHOD) reacts with non-LDL lipoproteins (chylomicrons, VLDL, HDL). Hydrogen peroxide produced by the enzyme reactions with non-LDL cholesterol is decomposed by catalase in Reagent 1. When Reagent 2 is added, the protective reagent is removed from LDL and catalase is inactivated by sodium azide (NaN<sub>3</sub>). In this second process, CHE and CHOD react only with LDL cholesterol. Hydrogen peroxide produced by the enzyme reactions with LDL cholesterol yields a color complex upon oxidative condensation with HDAOS [N-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline] and 4-aminoantipyrine in the presence of peroxidase (POD). The absorbance of the blue color complex is measured at 600 nm.

## REAGENT COMPOSITION

### R1

|                             |              |
|-----------------------------|--------------|
| Good's buffer, pH 6.8       | 25 mmol/l    |
| Cholesterol esterase (CHES) | 83.3 µkat/l  |
| Cholesterol oxidase (CHOD)  | 83.3 µkat/l  |
| HDAOS                       | 0.64 mmol/l  |
| Catalase                    | 16.66 µkat/l |

### R2

|                          |               |
|--------------------------|---------------|
| Good's buffer, pH 7.0    | 25 mmol/l     |
| 4-aminoantipyrine (4-AA) | 3.4 mmol/l    |
| Peroxidase (POD)         | 333.33 µkat/l |
| NaN <sub>3</sub>         | <0.1 %        |

## REACTION MIXTURE

|                             |              |
|-----------------------------|--------------|
| Good's buffer               | 25 mmol/l    |
| 4-aminoantipyrine (4-AA)    | 0.84 mmol/l  |
| Peroxidase (POD)            | 82.65 µkat/l |
| Cholesterol esterase (CHES) | 61.96 µkat/l |
| Cholesterol oxidase (CHOD)  | 61.96 µkat/l |
| Catalase                    | 12.39 µkat/l |
| HDAOS                       | 0.47 mmol/l  |
| NaN <sub>3</sub>            | 0.025 %      |

## REAGENT PREPARATION

Reagents R1 and R2 are liquid, ready to use.

## STABILITY AND STORAGE

If stored at 2–8 °C, reagents are stable until expiry date, that is stated on the package.

After first opening, reagents R1 and R2 are stable for 30 days at 2–8 °C in the dark and without contamination.

## SPECIMEN COLLECTION & HANDLING

Use serum or heparin plasma.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

|                            |          |             |
|----------------------------|----------|-------------|
| Stability in serum/plasma: | 12 hours | at 20–25 °C |
|                            | 10 days  | at 4–8 °C   |
|                            | 12 weeks | at -20 °C   |

Discard contaminated specimens.

## CALIBRATION

For the calibration, it is recommended to use LDL Cholesterol Calibrator, Cat. No. BLT00073.

## QUALITY CONTROL

For quality control ERBA NORM 4x5, Cat. No. BLT00080 or ERBA NORM 10x5, Cat. No. XSYS0123 and ERBA PATH 4x5, Cat. No. BLT00081 or ERBA PATH 10x5, Cat. No. XSYS0124 are recommended.

## UNIT CONVERSION

mg/dl x 0.026 = mmol/l

## EXPECTED VALUES <sup>2</sup>

|                             |           |
|-----------------------------|-----------|
| fS LDL cholesterol (mmol/l) | < 3.4     |
| Medium risk                 | 3.4 – 4.1 |
| High risk                   | > 4.1     |

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

## PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 0.09 mmol/l

Linearity: 10 mmol/l

Measuring range: 0.09 – 10 mmol/l

| Intra-assay precision<br>Within run (n=20) | Mean<br>(mmol/l) | SD<br>(mmol/l) | CV<br>(%) |
|--|------------------|----------------|-----------|
| Sample 1                                   | 2.53             | 0.023          | 0.90      |
| Sample 2                                   | 3.66             | 0.031          | 0.84      |

| Inter-assay precision<br>Run to run (n=20) | Mean<br>(mmol/l) | SD<br>(mmol/l) | CV<br>(%) |
|--|------------------|----------------|-----------|
| Sample 1                                   | 2.46             | 0.026          | 1.08      |
| Sample 2                                   | 3.56             | 0.045          | 1.25      |

## COMPARISON

A comparison between LDL CHOLESTEROL DIRECT (y) and a commercially available test (x) using 40 samples gave following results:

y = 1.019 x - 0.068 mmol/l

r = 0.994

## INTERFERENCES

Following substances do not interfere:

bilirubin up to 55 mg/dl, triglycerides up to 800 mg/dl, haemoglobin up to 7.5 g/l.

Interference by N-acetylcysteine (NAC), acetaminophen and metamazole causes falsely low results. To carry out the test, blood withdrawal should be performed prior to administration of drugs.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

## Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagent R1:

EUH208 Contains reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one. May produce an allergic reaction.

EUH210 Safety data sheet available on request.

Reagent R2:

EUH208 Contains Alcohols, C11-15-secondary, ethoxylated. May produce an allergic reaction.

EUH210 Safety data sheet available on request.

## FIRST AID

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

## WASTE DISPOSAL

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials.

Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

## PROCEDURE

Wavelength: 600/700 nm

Cuvette: 1 cm

Temperature: 37 °C

Serum/reaction mixture ratio 1/121

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio.

|   | Reagent blank | Calibrator | Sample   |
|---|---------------|------------|----------|
| Reagent R1  | 0.900 ml      | 0.900 ml   | 0.900 ml |
| Sample  | -             | -          | 0.010 ml |
| Calibrator  | -             | 0.010 ml   | -        |
| Distilled water   | 0.010 ml      | -          | -        |
| Mix and incubate 5 min. at 37 °C. Measure absorbance A <sub>1</sub> of the sample (A <sub>sam</sub> ), calibrator (A <sub>cal</sub> ) and blank (A <sub>bl</sub> ). Then add: |               |            |          |
| Reagent R2  | 0.300 ml      | 0.300 ml   | 0.300 ml |

Mix and incubate 5 min. at 37 °C. Measure absorbance A<sub>2</sub> of the sample (A<sub>sam</sub>), calibrator (A<sub>cal</sub>) and blank (A<sub>bl</sub>).

Calculate absorbance change ΔA = A<sub>2</sub> – A<sub>1</sub> for sample, calibrator and blank.

## CALCULATION

$$\text{LDL Cholesterol (mmol/l)} = \frac{\Delta A_{\text{sam}} - \Delta A_{\text{bl}}}{\Delta A_{\text{cal}} - \Delta A_{\text{bl}}} \times C_{\text{cal}}$$

C<sub>cal</sub> = calibrator concentration

## NOTES

- Li-Heparinate plasma values on average are recovered by 3% lower than serum concentrations. For EDTA plasma, cca 9% value decrease against serum is expected.
- For the determination of LDL cholesterol it is recommended to use serum or plasma with blood elements removed within 3 hours after collection.
- If concentration of LDL cholesterol exceeds 10.4 mmol/l, dilute the sample with a saline in ratio 1 + 1 and results multiply with dilution factor.
- Do not use Reagents, which were frozen by mistake.

Applications for automatic analysers will be supplies on request.

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