MICROPROTEIN

Cat. No.	Pack Name	Packaging (Content)
BLT00045	MP 100	R1: 2 x 50 ml, R2 standard: 1 x 1 ml

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Protein in human urine and CSF.

CLINICAL SIGNIFICANCE

The role of the renal system in the conversion of plasma proteins has been recognized for some time. Under normal physiological conditions small molecular weight proteins such as insulin pass through the glomeruli in relatively large amounts. Intermediate size proteins such as Transferrin and Albumin also pass through but only in relatively small amounts. Most of these proteins are reabsorbed in the renal tubules such that normal urine contains less than 150 mg of protein per day. This also includes the protein of non serum origin normally secreted by the distal tubule (muco protein) and collecting ducts. Increased levels of urinary protein, (proteinuria) usually more than 0.15 g per 24 hours (150 mg/24 hours), almost always indicates disease. Proteinuria may be classified as renal proteinurea or proteinuria with normal renal function. Renal proteinuria may be further classified as Glomerular or tubular proteinuria.

Glomerular proteinuria is due to increased glomerular permeability (nephrotic syndrome) and may be seen in glomerular nephritis or secondary to other diseases such as diabetic nephropathy. Albumin is usually the predominant protein in the urine. Tubular proteinuria may be due to renal tubular damage from any cause especially pyelonephritis. Tubular proteinuria results in modest increases in the low molecular weight proteins. Increased protein excretion is seen during normal pregnancy, after strenuous exercise or following prolonged maintenance of an upright posture. Increase in low molecular weight proteins may be due to the production of Bence Jones protein, haemoglobinuria as a result of severe haemolysis and myoglobinuria as a result of severe muscle damage.

PRINCIPLE

Methods employed for the determination of total protein in urine include dye binding, chemical and turbidimetric procedures, the latter being the most commonly employed technique. In this microprotein kit the pyrogallol dye combines with molybdenum acid forming a red complex with maximum absorbance at 470 nm.

When this complex is combined with protein under acidic conditions its maximum absorption is shifted to a longer wavelength $\lambda_{\text{max}} = 600 \text{ nm}$. The concentration of the protein can be obtained by measuring absorbance at 600 nm.

REAGENT COMPOSITION

R1	
Succinate buffer	15 mmol/l
Pyrogallol red	0.060 mmol/l
Ammonium molybdate	0.043 mmol/l
R2 standard	see bottle label

REAGENT PREPARATION

Reagents are liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

SPECIMEN COLLECTION AND HANDLING

Use urine or cerebrospinal fluid (CSF).

It is recommended to follow NCCLS procedures (or similar standardized conditions).

urine:

1 day	at 20–25°C
7 days	at 4–8 °C
1 month	at -20°C

cerebrospinal fluid:

1 day	at 20–25°C
6 days	at 4–8°C
1 year	at -20°C

Discard contaminated specimens.

CALIBRATION

Calibration with the standard included in the kit is recommended.

QUALITY CONTROL

For quality control all control solutions with protein total values determined by this method can be used.

UNIT CONVERSION

mg/dl x 10 = mg/l

EXPECTED VALUES ⁷

Urine: 10–140 mg/day (24 hours sample)

CSF: < 50 mg/dl *

*(the value is an approximate guideline only)

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 1.9 mg/dl

Linearity: 300 mg/dl

Measuring range: 1.9–300 mg/dl

PRECISION

Intra-assay precision Within run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	22.7	0.74	3.28
Sample 2	11.2	0.43	3.86

Inter-assay precision Run to run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	16.8	0.56	3.35
Sample 2	60.6	1.95	3.22

COMPARISON

A comparison between XL-Systems Microprotein (y) and a commercially available test (x) using 40 samples gave following results:

$y = 1.017 x - 0.46 \text{ mg/dl}$

$r = 0.999$

INTERFERENCES

In this method, some kinds of surface active agents may affect the color. Cationic surfactants, in general, give the same color as proteins. Because anions inhibit the color reaction, wash the equipment thoroughly, using distilled water, until no surface active agent remains. Then dry equipment completely before using it.

Haemoglobin shows about one-half the color of albumin. If haematuria is present, a falsely high value will be measured.

The Color Reagent contains components which are unstable in light. Because of this, when the Color Reagent is transferred to another container, it should be a light-resistant one.

Small amounts of protein attached to the cuvette wall after measurement of certain other tests will cause an erroneously high measured value when the test solution is transferred to the cuvette. If this should occur, wash the cuvette completely and measure again. If the absorbed protein cannot be removed completely by washing with water, clean the cuvette with an alkaline solution containing hypochlorite and wash thoroughly with water.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1 UFI: GA8F-XE48-QH77-UD06

Reagent R1 contains methanol.

Warning

Hazard statement(s):

H371 May cause damage to organs (eyes).

Precautionary statement(s):

P260 Do not breathe vapours/spray.

P264 Wash hands thoroughly after handling.

P308+P311: IF exposed or concerned: Call a POISON CENTER or doctor.

Reagent R2 standard is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

ASSAY PROCEDURE

Wavelength: 600 nm

Cuvette: 1 cm

Macro method

Add to tubes marked	Blank	Standard	Test
Sample	–	–	10 µl
Standard	–	10 µl	–
Water	10 µl	–	–
Dye reagent	1000 µl	1000 µl	1000 µl

Mix well. Incubate at 37°C for 10 minutes. Measure absorbance of Standard and test against Blank within 60 minutes.

Micro method

Add to tubes marked	Blank	Standard	Test
Sample	–	–	5 µl
Standard	–	5 µl	–
Water	5 µl	–	–
Dye reagent	500 µl	500 µl	500 µl

Mix well. Incubate at 37°C for 10 minutes. Measure absorbance of Standard and test against the reagent Blank within 60 minutes.

CALCULATION

$$\text{Protein (mg/dl)} = \frac{\Delta A_{\text{sam}}}{\Delta A_{\text{st}}} \times C_{\text{st}} \quad C_{\text{st}} = \text{standard concentration}$$

Applications for automatic analysers are available on request.

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	End point
Wavelength 1 (nm)	600
Sample Volume (µl)	5/10
Reagent Volume (µl)	500/1000
Incubation time (min.)	10
Incubation temp. (°C)	37
Normal Low (mg/dl)	0
Normal High (mg/dl)	15
Linearity Low (mg/dl)	1.9
Linearity High (mg/dl)	300
Concentration of Standard	See bottle label
Blank with	Reagent
Absorbance limit (max.)	0.4
Units	mg/dl

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erba.com, www.erbamannheim.com

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