PHOSPHORUS

Cat. No.	Pack Name	Packaging (Content)
BLT00047	PHOS 100	R1: 2 × 50 ml, R2 standard: 1 × 5 ml
BLT00048	PHOS 250	R1: 1 × 250 ml, R2 standard: 1 × 5 ml

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Phosphorus in human serum, plasma or urine.

CLINICAL SIGNIFICANCE

More than 80 % of the body's phosphate is present in bones as calcium phosphate. The remainder is found intracellularly as organic phosphates such as phospholipids, nucleic acids and ATP or extracellularly as inorganic phosphorus.

There is generally a reciprocal relationship between serum calcium and inorganic phosphorus levels. Increased levels of serum phosphorus is seen in renal diseases, hypoparathyroidism and excessive vitamin D intake.

Decreased levels of phosphorus is seen in rickets, osteomalacia (adult rickets), hyperparathyroidism and in diabetic coma.

PRINCIPLE

Inorganic phosphorus combines with ammonium molybdate in the presence of strong acids to form phosphomolybdate. The formation of reduced phosphomolybdate is measured at 340 nm and is directly proportional to the concentration of inorganic phosphorus present in the sample.

REAGENT COMPOSITION

R1

Ammonium Molybdate 1.00 mmol/l
Sulphuric acid 336 mmol/l

R2 standard See bottle label

REAGENT PREPARATION

Reagents are liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

SPECIMEN COLLECTION AND HANDLING

Use unheamolytic serum or plasma (heparin) or urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability in serum / plasma:

7 days at 4–25°C
3 months at -20°C

Stability in urine:

2 days at 20–25°C at pH < 5

For collection of 24 h urine add 20 to 30 ml of 6 mol/l HCl into the collection bottle to avoid phosphate precipitations.

Dilute urine 1 + 10 with dist. water before determination and multiply the result by 11. Discard contaminated specimens.

CALIBRATION

Calibration with the standard included in the kit or the calibrator XL MULTICAL 4x3 (Cat.No. XSYS0034) or XL MULTICAL 10x3 (Cat.No. XSYS0122) is recommended.

QUALITY CONTROL

For quality control ERBA NORM 4x5 (Cat.No. BLT00080) or ERBA NORM 10x5 (Cat.No. XSYS0123) and ERBA PATH 4x5 (Cat.No. BLT00081) or ERBA PATH 10x5 (Cat.No. XSYS0124) are recommended.

UNIT CONVERSION

mg/dl × 0.32 = mmol/l

EXPECTED VALUES *

Serum

Adult 2.5–4.5 mg/dl
Children 4.0–7.0 mg/dl

Urine, 24 h

Adult 0.4–1.3 g / 24 h

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 0.20 mg/dl

Linearity: 22.9 mg/dl

Measuring range: 0.20–22.9 mg/dl

PRECISION

Intra-assay precision Within run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	5.29	0.05	0.94
Sample 2	3.33	0.04	1.25

Inter-assay precision Run to run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	3.11	0.08	2.65
Sample 2	7.00	0.08	1.13

COMPARISON

A comparison between XL-Systems Phosphorus (y) and a commercially available test (x) using 40 samples gave following results:

y = 1.053x + 0.034 mg/dl

r = 0.996

INTERFERENCES

Following substances do not interfere:

haemoglobin up to 2.5 g/l, bilirubin up to 20 mg/dl, triglycerides up to 500 mg/dl.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1 UFI: 9YDF-8EDQ-8H7H-KTFP

Reagent 1 of the kit contains Triton X-100.

Danger

Hazard statement:

H314 Causes severe skin burns and eye damage.

H412 Harmful to aquatic life with long lasting effects.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P301+P330+P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Reagent 2 of the kit is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

ASSAY PROCEDURE

Wavelength: 340 nm, Hg 334 nm, Hg 365 nm

Cuvette: 1 cm

	Reagent blank	Standard (Cal.)	Sample
Reagent 1	1.00 ml	1.00 ml	1.00 ml
Sample	–	–	0.01 ml
Standard (Cal.)	–	0.01 ml	–
Distilled water	0.01 ml	–	–

Mix and incubate 5 min. at 37°C. Measure absorbance of the sample A_{sam} and standard A_{st} against reagent blank.

CALCULATION

Phosphorus (mmol/l) = $\frac{\Delta A_{sam}}{\Delta A_{st}} \times C_{st}$ C_{st} = standard (calibrator) concentration

Applications for automatic analysers are available on request.

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	End Point
Wavelength 1 (nm)	340
Sample Volume (µl)	5/10
Reagent Volume (µl)	500/1000
Incubation time (min.)	5
Incubation temp. (°C)	37
Normal Low (mg/dl)	2.5
Normal High (mg/dl)	4.5
Linearity Low (mg/dl)	0.20
Linearity High (mg/dl)	22.9
Concentration of standard	See bottle label
Blank with	Reagent blank
Absorbance limit (max.)	2.5
Units	mg/dl

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erba.com, www.erbamannheim.com

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