# TOTAL PROTEIN

| Cat. No. | Pack Name | Packaging (Content)                   |
|----------|-----------|---------------------------------------|
| BLT00054 | TP 250    | R1: 5 x 50 ml, R2 standard: 1 x 3 ml  |
| BLT00055 | TP 500    | R1: 2 x 250 ml, R2 standard: 1 x 3 ml |

## INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Total Protein in human serum, plasma and urine.

## CLINICAL SIGNIFICANCE

Total protein is useful for monitoring gross changes in protein levels caused by various disease states. It is usually performed in conjunction with other tests such as serum albumin, liver function tests or protein electrophoresis. An albumin/globulin ratio is often calculated to obtain additional information.

Increased levels of serum protein are observed in dehydration, multiple myeloma and chronic liver disease.

Decreased levels are encountered in renal diseases and terminal liver failure.

## PRINCIPLE

Biuret method. The peptide bonds of protein react with copper II ions in alkaline solution to form a blue-violet ion complex, (the so called biuret reaction), each copper ion complexing with 5 or 6 peptide bonds. Tartrate is added as a stabiliser whilst iodide is used to prevent auto-reduction of the alkaline copper complex. The colour formed is proportional to the protein concentration and is measured at 546 nm (520-560).

## REAGENT COMPOSITION

### R1

|                           |                  |
|---------------------------|------------------|
| Copper II Sulphate        | 12 mmol/l        |
| Potassium Sodium Tartrate | 31.9 mmol/l      |
| Potassium Iodide          | 30.1 mmol/l      |
| Sodium Hydroxide          | 0.6 mol/l        |
| <b>R2 Standard</b>        | See bottle label |

## REAGENT PREPARATION

Reagents are liquid, ready to use.

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

**R1:** after opening is stable till the expiry date stated on the bottle and kit label when stored at 2–25°C in dark, avoid contamination.

**R2 standard:** after opening is stable till the expiry date stated on the bottle and kit label when stored at 2–8°C in dark, avoid contamination.

## SPECIMEN COLLECTION AND HANDLING

Use unheamolytic serum, plasma (heparin, EDTA), urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

### Stability

|                   |            |
|-------------------|------------|
| 6 days            | at 20–25°C |
| 4 weeks           | at 4–8°C   |
| at least one year | at -20°C   |

Discard contaminated specimens.

## CALIBRATION

Calibration with the standard included in the kit or the calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

## QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

## UNIT CONVERSION

g/dl x 10 = g/l

## EXPECTED VALUES <sup>5</sup>

|               | (g/dl)    |
|---------------|-----------|
| Adults:       | 6.4 – 8.3 |
| Premature     | 3.6 – 6.0 |
| Newborn       | 4.6 – 7.0 |
| 1 week        | 4.4 – 7.6 |
| 7 – 12 months | 5.1 – 7.3 |
| 1 – 2 years   | 5.6 – 7.5 |
| > 2 years     | 6.0 – 8.0 |

**It is recommended that each laboratory verify this range or derives reference interval for the population it serves.**

## PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

**Limit of quantification:** 0.37 g/dl

**Linearity:** 15 g/dl

**Measuring range:** 0.37 – 15 g/dl

## PRECISION

| Intra-assay precision<br>Within run (n=20) | Mean<br>(g/dl) | SD<br>(g/dl) | CV<br>(%) |
|--|----------------|--------------|-----------|
| <b>Sample 1</b>                            | 4.914          | 0.072        | 1.46      |
| <b>Sample 2</b>                            | 7.314          | 0.040        | 0.55      |

| Inter-assay precision<br>Run to run (n=20) | Mean<br>(g/dl) | SD<br>(g/dl) | CV<br>(%) |
|--|----------------|--------------|-----------|
| <b>Sample 1</b>                            | 5.177          | 0.038        | 0.73      |
| <b>Sample 2</b>                            | 6.761          | 0.078        | 1.15      |

## COMPARISON

A comparison between XL-Systems Total Protein (y) and a commercially available test (x) using 40 samples gave following results:

$$y = 0.986x + 0.163 \text{ g/dl} \quad r = 0.997$$

## INTERFERENCES

Following substances do not interfere:  
haemoglobin up to 7.5 g/l, bilirubin up to 40 mg/dl, triglycerides up to 1500 mg/dl.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagent 1 contains 2.4 % sodium hydroxide which is classified as corrosive.

Danger

### Hazard statement:

H314 Causes severe skin burns and eye damage.

H412 Harmful to aquatic life with long lasting effects.

## Precautionary statement:

P260 Do not breathe vapours.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301+P330+P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

R2 Standard is not classified as dangerous. It contains less than 0.1% sodium azide which is classified as very toxic and dangerous substance for the environment.

## WASTE MANAGEMENT

Please refer to local legal requirements.

## ASSAY PROCEDURE

**Wavelength:** (520-560) nm

**Cuvette:** 1 cm

|                       | Reagent blank | Sample  | Standard<br>(Calibrator) |
|-----------------------|---------------|---------|--------------------------|
| Reagent R1            | 1.00 ml       | 1.00 ml | 1.00 ml                  |
| Distilled water       | 0.02 ml       | –       | –                        |
| Sample                | –             | 0.02 ml | –                        |
| Standard (calibrator) | –             | –       | 0.02 ml                  |

Mix and incubate for 10 minutes (in case of automatic procedure incubate for 5 minutes) incubation in the dark. Absorbance of the sample  $A_1$  and the standard (calibrator)  $A_2$  against reagent blank is read in interval 30 minutes.

## CALCULATION

$$\text{Total protein (g/dl)} = \frac{A_1}{A_2} \times C_{st} \quad C_{st} = \text{standard (calibrator) concentration}$$

Applications for automatic analysers are available on request.

## ASSAY PARAMETERS

| Mode                      | End Point        |
|---------------------------|------------------|
| Wavelength 1 (nm)         | 546              |
| Sample Volume (μl)        | 10/20            |
| Reagent Volume (μl)       | 500/1000         |
| Incubation time (min.)    | 10               |
| Incubation temp. (°C)     | 37               |
| Normal Low (g/dl)         | 6.4              |
| Normal High (g/dl)        | 8.3              |
| Linearity Low (g/dl)      | 0.37             |
| Linearity High (g/dl)     | 15               |
| Concentration of Standard | See bottle label |
| Blank with                | Reagent          |
| Absorbance limit (max.)   | 0.4              |
| Units                     | g/dl             |

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ  
e-mail: diagnostics@erbamannheim.com, www.erbamannheim.com

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