TRIGLYCERIDES

Cat. No.	Pack Name	Packaging (Content)
BLT00057	TG 100	R1: 2 x 50 ml, R2 standard: 1 x 3 ml
BLT00058	TG 1000	R1: 1 x 1000 ml
BLT00059	TG 250	R1: 1 x 250 ml, R2 standard: 1 x 3 ml

INTENDED USE


Diagnostic reagent for quantitative *in vitro* determination of Triglycerides in human serum and plasma.

CLINICAL SIGNIFICANCE

Triglycerides are a family of lipids absorbed from the diet and produced endogenously from carbohydrates. Measurement of triglycerides is important in the diagnosis and management of hyperlipidemias. These diseases can be genetic or secondary to other disorders including nephrosis, diabetes mellitus and endocrine disturbances. Elevation of triglycerides has been identified as a risk factor for atherosclerotic disease.

PRINCIPLE

The series of reactions involved in the assay system is as follows:

Triglycerides are enzymatically hydrolyzed by lipase to free acids and glycerol. The glycerol is phosphorylated by adenosine triphosphate (ATP) with glycerol kinase (GK) to produce glycerol-3-phosphate and adenosine diphosphate (ADP). Glycerol-3-phosphate is oxidized to dihydroxy-acetone phosphate (ADP) by glycerol phosphate oxidase producing hydrogen peroxide (H₂O₂). In a Trinder type colour reaction catalyzed by peroxidase, the H₂O₂ reacts with 4-aminoantipyrine (4AAP) and 4-chlorophenol to produce a red coloured dye. The absorbance of this dye is proportional to the concentration present in the sample.

REAGENT COMPOSITION

R1	
Good's buffer (pH 7.2)	50 mmol/l
4-Chlorophenol	4 mmol/l
Mg ²⁺	15 mmol/l
ATP	2 mmol/l
Glycerolkinase	≥ 0.4 KU/l
Peroxidase	≥ 2.0 KU/l
Lipoproteinlipase	≥ 2.0 KU/l
Glycerol-3-phosphate-Oxidase	≥ 0.5 KU/l
4-Aminoantipyrine	0.5 mmol/l
R2 standard	See bottle label

REAGENT PREPARATION

Reagents are liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

SPECIMEN COLLECTION AND HANDLING

Use unheamolytic serum, plasma (EDTA, heparin)
It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability

2 days	at 20–25°C
7 days	at 4–8°C
at least one year	at -20°C

Discard contaminated specimens.

CALIBRATION

Calibration with the standard included in the kit or the calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

UNIT CONVERSION

mg/dl x 0.0113 = mmol/l

EXPECTED VALUES ⁴

Recommended Triglycerides levels for adults:

Normal:	< 150 mg/dl
High:	150–199 mg/dl
Hypertriglyceridemic:	200–499 mg/dl
Very high:	> 499 mg/dl

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification:	9.74 mg/dl
Linearity:	1062 mg/dl
Measuring range:	9.74–1062 mg/dl

PRECISION

Intra-assay precision Within run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	102.7	0.53	0.52
Sample 2	186.7	1.68	0.90

Inter-assay precision Run to run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	108.0	1.68	1.56
Sample 2	184.5	2.83	1.54

COMPARISON

A comparison between XL-Systems Triglycerides (y) and a commercially available test (x) using 40 samples gave following results:

$$y = 1.022x - 1.15 \text{ mg/dl}$$

$$r = 0.999$$

INTERFERENCES

Following substances do not interfere:
haemoglobin up to 10 g/l, bilirubin up to 40 mg/dl.
Interference by N-acetylcysteine (NAC), acetoaminophen and metamazole causes falsely low results. To carry out the test, blood withdrawal should be performed prior to administration of drugs.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagents of the kit are not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

ASSAY PROCEDURE

Wavelength	500 (546) nm
Cuvette	1 cm

	Reagent blank	Standard (calibrator)	Sample
Reagent 1	1.00 ml	1.00 ml	1.00 ml
Sample	–	–	0.01 ml
Standard (calibrator)	–	0.01 ml	–
Distilled water	0.01 ml	–	–

Mix and incubate 10 min. at 37 °C. Measure absorbance of the sample A_{sam} and standard A_{st} against reagent blank. The coloration is stable during one hour.

CALCULATION

$$\text{Triglycerides (mg/dl)} = \frac{\Delta A_{\text{sam}}}{\Delta A_{\text{st}}} \times C_{\text{st}} \quad C = \text{concentration of standard (calibrator)}$$

Applications for automatic analysers are available on request.

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	End Point
Wavelength 1 (nm)	505
Wavelength 2 (nm)	670
Sample Volume (μl)	5/10
Reagent Volume (μl)	500/1000
Incubation time (min.)	10
Incubation temp. (°C)	37
Normal Low (mg/dl)	0
Normal High (mg/dl)	150
Linearity Low (mg/dl)	9.7
Linearity High (mg/dl)	1050
Concentration of Standard	See bottle label
Blank with	Reagent
Absorbance limit (max.)	0.5
Units	mg/dl

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REFERENCES / ЛИТЕРАТУРА / ЛІТЕРАТУРА / LITERATURA / LITERATÚRA / REFERENCIAS

1. Rifai N, Bachorik PS, Alberts JJ. Lipids, lipoproteins and apolipoproteins. In: Burtis CA, Ashwood ER, editors Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 809-61.
2. Cole, TG, Klotzsch SG, McNamara J. Measurement of triglyceride concentration. In: Rifai N, Warnick GR, Dominiczak MH, eds. Handbook of lipoprotein testing. Washington: AACC Press, 1997. p. 155-26
3. Recommendation of the Second Joint Task Force of European and other Societies on Coronary Prevention. Prevention of coronary heart disease in clinical practice. Eur Heart J 1998; 19: 1434-503
4. Tietz Textbook of Clinical Chemistry and Molecular diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Company, 2012.

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