CHOLINESTERASE

Cat. No.	Pack Name	Packaging (Content)
BLT00066	CHE 120	R1: 2 x 50 ml, R2: 1 x 20 ml

EN

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of cholinesterase in human serum and plasma.

CLINICAL SIGNIFICANCE

Two related enzymes have the ability to hydrolyze acetylcholine. One is acetylcholinesterase, which is called true cholinesterase, or choline esterase I. True cholinesterase is found in erythrocytes, in the lungs and spleen, in nerve endings, and in the gray matter of the brain.

The other cholinesterase is acylcholine acylhydrolase, it is usually called pseudocholinesterase, benzoyl cholinesterase, or choline asterase II. Although it is found in the liver, pancreas, heart, white matter of the brain, and serum, its biological role is unknown.

The serum enzyme is the one whose assay is clinically useful.

Cholinesterase levels in serum are useful as a test of liver function, as an indicator of possible insecticide poisoning, or for the detection of patients with atypical forms of the enzyme. Decreased levels of serum enzyme are also found in patients with acute infections, pulmonary embolism, muscular dystrophy, chronic renal disease and in pregnancy.

PRINCIPLE

Cholinesterase (CHE) catalyses the hydrolysis of butyrylthiocholine substrate forming butyrate and thiocholine. Thiocholine reduces hexacyanoferrate (III⁺) to hexacyanoferrate (II⁺). The decrease of absorbance at 405 nm is proportional to the activity of CHE in the sample.

REAGENT COMPOSITION

R1	
Pyrophosphate buffer, pH 7.6 (37 °C)	92 mmol/l
Hexacyanoferrate (III ⁺)	2.5 mmol/l
R2	
Butyrylthiocholine	91 mmol/l

REACTION MIXTURE

Pyrophosphate buffer, pH 7.6 (37 °C)	75.41 mmol/l
Hexacyanoferrate (II ⁺)	2.05 mmol/l
Butyrylthiocholine	14.92 mmol/l

REAGENT PREPARATION

Reagents are liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

Two reagents method

Reagents R1 and R2 are liquid, ready to use.

After first opening, reagents are stable for 90 days at 2–8 °C if stored at appropriate conditions, closed carefully and without any contamination.

SPECIMEN COLLECTION & HANDLING

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Fresh serum, plasma (heparin, EDTA), not haemolysed and promptly separated from the red blood cells.

Do not use sodium fluoride as an anticoagulants because it inhibits cholinesterase.

Stability of cholinesterase in serum and plasma is 7 days at 4–8 °C.

Stability

in serum /plasma: 7 days at 2–8 °C
1 year at -20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with LYONORM CALIBRATOR, Cat. No. BLT00069 is recommended.

QUALITY CONTROL

For quality control LYONORM HUM N, Cat. No. BLT00070 and Lyonorm HUM P, Cat. No. BLT00071 are recommended.

UNIT CONVERSION

U/l x 0.017 = μ kat/l

EXPECTED VALUES

fS CHE (μ kat/l) 37 °C
Male 85 – 195
Female 66.7 – 210

In infants up to 6 months of age, cholinesterase activity is 40% to 50% higher than in adults. In young women, the enzyme activity is approximately 64% to 74% higher than in adult males. The activity decreases in pregnancy.

The range of reference values is only approximate; it is recommended that each laboratory verify the extent of the reference interval for their particular examined population.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 5.7 μ kat/l

Linearity: 330 μ kat/l

Measuring range: 5.7 – 330 μ kat/l

PRECISION

Intra-assay (n=20)	Mean (μ kat/l)	SD (μ kat/l)	CV (%)
Sample 1	50.9	0.59	1.16
Sample 2	73.3	0.87	1.19

Inter-assay (n=20)	Mean (μ kat/l)	SD (μ kat/l)	CV (%)
Sample 1	57.45	2.29	3.99
Sample 2	99.9	2.38	2.38

COMPARISON

A comparison between XL-Systems Cholinesterase (y) and a commercially available test (x) using 40 samples gave following results:

N = 40

r = 0.998

y = 1.032 x – 1.06 μ kat/l

INTERFERENCES

Following substances do not interfere:

haemoglobin up to 5 g/l, bilirubin up to 15 mg/dl, triglycerides up to 500 mg/dl.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Hazards identification in accordance with Regulation (EC) No 127 2/2008

R2 UFI: XFJU-7W7J-SJ5M-36KS

Reagent 1 is not classified as dangerous.

Reagent 2 contains <0.5% maleic acid, which is a skin sensitiser.

Warning

Hazard statement:

H317 May cause an allergic skin reaction.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P302+P352 IF ON SKIN: Wash with plenty of water and soap.

P333+P313 If skin irritation or rash occurs: Get medical advice / attention.

FIRST AID

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

WASTE DISPOSAL

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials. Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

ASSAY PROCEDURE

Wavelength: 405 nm

Cuvette: 1 cm

Temperature: 37 °C

Serum/reaction mixture ratio 1/61

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio.

Two reagents method – substrate start

	Reagent blank	Calibrator	Sample
Reagent 1	1.000 ml	1.000 ml	1.000 ml
Sample	-	-	0.020 ml
Calibrator	-	0.020 ml	-
Distilled water	0.020 ml	-	-
Mix and incubate for 5 minutes at 37 °C. Then add:			
Reagent 2	0.200 ml	0.200 ml	0.200 ml

Mix and incubate for 90 seconds at 37 °C and read initial absorbance of sample, calibrator and reagent blank. Repeat the reading after exactly 30, 60 a 90 seconds. Calculate the mean 30 seconds absorbance change ($\Delta A/30s$).

CALCULATION

$$1. \text{Cholinesterase } (\mu\text{kat/l}) = C_{\text{cal}} \times \frac{\Delta A_{\text{sam}}/30 \text{ s} - \Delta A_{\text{br}}/30 \text{ s}}{\Delta A_{\text{cal}}/30 \text{ s} - \Delta A_{\text{br}}/30 \text{ s}}$$

C_{cal} = calibrator concentration

2. Using factor

Cholinesterase (μ kat/l) = f x $\Delta A/30s$

f = factor:

405 nm	Serum / plasma
Substrate start	2193

Applications for automatic analysers are available on request.

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erba.com, www.erbamannheim.com

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