# BILIRUBIN TOTAL DCA

Cat. No.	Pack Name	Packaging (Content)
BLT00077	BIL T DCA 500	R1: 4 × 100 ml, R2: 1 × 100 ml

## INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Total Bilirubin in human serum and plasma.

## CLINICAL SIGNIFICANCE

Bilirubin is a breakdown product of haemoglobin. Bilirubin formed in the reticulo endothelial system is transported bound by albumin to the liver. This bilirubin is water insoluble and is known as indirect or unconjugated bilirubin. In the liver, bilirubin is conjugated to glucuronic acid to form direct bilirubin. Conjugated bilirubin is excreted via the biliary system into the intestine. Here it is metabolised by bacteria to urobilinogen and stercobilinogen.

TOTAL BILIRUBIN = INDIRECT BILIRUBIN + DIRECT BILIRUBIN

Total Bilirubin is elevated in obstructive conditions of the bile duct, hepatitis, cirrhosis, in haemolytic disorders and several inherited enzyme deficiencies.

Indirect Bilirubin is elevated by pre-hepatic causes such as haemolytic disorders or liver diseases resulting in impaired entry transport or conjugation within the liver. Monitoring of indirect bilirubin in neonates is of special importance as it is the indirect (or free) bilirubin bound to albumin that is able to cross the blood brain barrier more easily increasing the danger of cerebral damage.

## PRINCIPLE

Bilirubin reacts with diazotized 2,4-dichloroaniline to produce an intensely coloured diazo compound. The intensity of color of this dye in solution is proportional to the concentration of total bilirubin.

## REAGENT COMPOSITION

### R1

Phosphate buffer	50 mmol/l
NaCl	150 mmol/l
Surfactant	

### R2

2,4-dichlorophenyl diazonium salt	5 mmol/l
Hydrochloric acid	130 mmol/l
Surfactants	

## REACTION MIXTURE

Phosphate buffer	39 mmol/l
NaCl	118 mmol/l
2,4-dichlorophenyl diazonium salt	1 mmol/l
Hydrochloric acid	25 mmol/l

## REAGENT PREPARATION

Reagents R1 and R2 are liquid, ready to use.

## STABILITY AND STORAGE

If stored at 2–8°C, reagents are stable until expiry date, that is stated on the package. After opening, reagents are stable until expiry date at 2–8°C if stored at appropriate conditions, closed carefully and without any contamination.

Maximum allowable absorbance of the working reagent (mix 4 portion of reagent R1 with 1 portion of reagent R2) measured at 546 nm against distilled water is 0.1.

## SPECIMEN COLLECTION AND HANDLING

Use serum, plasma (heparin, EDTA). Protected from the light!

It is recommended to follow NCCLS procedures (or similar standardized conditions).

<b>Stability:</b>	1 day at 15–25°C
	7 days at 2–8°C
	3 months at -20°C

## FREEZE ONLY ONCE!

Discard contaminated specimens.

## CALIBRATION

Calibration with calibrator XL MULTICAL 4x3 (Cat. No. XSYS0034) or XL MULTICAL 10x3 (Cat. No. XSYS0122) is recommended.

## QUALITY CONTROL

For quality control ERBA NORM 4x5 (Cat. No. BLT00080) or ERBA NORM 10x5 (Cat. No. XSYS0123) and ERBA PATH 4x5 (Cat. No. BLT00081) or ERBA PATH 10x5 (Cat. No. XSYS0124) are recommended.

## UNIT CONVERSION

mg/dl × 17.1 = µmol/l

## EXPECTED VALUES

newborn ≤171  
adults and children (older than 1 month) 3.4–17.1

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

## PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

**Limit of quantification:** 0.73 µmol/l

**Linearity:** 630 µmol/l

**Measuring range:** 0.73–630 µmol/l

Intra-assay precision Within run (n=20)	Mean (µmol/l)	SD (µmol/l)	CV (%)
Sample 1	184.9	1.27	0.69
Sample 2	46.45	0.48	1.03

Inter-assay precision Run to run (n=20)	Mean (µmol/l)	SD (µmol/l)	CV (%)
Sample 1	19.82	0.52	2.63
Sample 2	80.55	1.59	1.97

## COMPARISON

A comparison between BIL T DCA 500 (y) and a commercially available test (x) using 40 samples gave following results:

y = 0.970 x + 0.410 µmol/l

r = 0.998

## INTERFERENCES

Following substances do not interfere:

triglycerides up to 2000 mg/dl, haemoglobin up to 5 g/l.

Eltrombopag medication can lead to falsely low or high results in patient samples.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

## Hazards identification in accordance with Regulation (EC) No 1272/2008

### R1:

UFI: UHWU-HW8S-9J5F-E7S3

**Danger**

Contains: Cetrimonium bromide, Oxytphenoxypolyethoxy-ethanol, Octylphenol-ethoxylate

### Hazard statement:

H290 May be corrosive to metals.

H314 Causes severe skin burns and eye damage.

H412 Harmful to aquatic life with long lasting effects.

### Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 If eye irritation persists: Get medical advice/attention.

### R2:

UFI: 7MWU-0WY5-KJ5X-3KC5

**Danger**

### Hazard statement:

H290 May be corrosive to metals.

H314 Causes severe skin burns and eye damage.

### Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 If eye irritation persists: Get medical advice/attention.

## WASTE MANAGEMENT

Please refer to local legal requirements.

## FIRST AID

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

## WASTE DISPOSAL

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials. Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

## PROCEDURE

**Wavelength:** 546 (540–560) nm

**Cuvette:** 1 cm

**Temperature:** 37 °C

Serum/reaction mixture ratio 1/50

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio.

	Reagent blank	Calibrator	Sample
Reagent R1	1.00 ml	1.00 ml	1.00 ml
Sample	-	-	0.025 ml
Calibrator	-	0.025 ml	-
Distilled water	0.025 ml	-	-

Mix and incubate 5 min at 37 °C or 10 min at (20–25)°C. Measure absorbance A1 of the sample ( $A_{sam}$ ), calibrator ( $A_{cal}$ ) against reagent blank. Then add:

Reagent R2	0.250 ml	0.250 ml	0.250 ml
------------	----------	----------	----------

Mix and incubate 5 min at 37 °C or 10 min at (20–25)°C. Measure absorbance A2 of the sample ( $A_{sam}$ ), calibrator ( $A_{cal}$ ) against reagent blank.

Calculate absorbance change  $\Delta A = A2 - A1$  for sample, calibrator and blank.

## CALCULATION

$$\text{Bilirubin } (\mu\text{mol/l}) = \frac{\Delta A_{sam}}{\Delta A_{cal}} \times C_{cal}$$

$C_{cal}$  = concentration of calibrator

Applications for automatic analysers will be supplied on request.

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ  
e-mail: diagnostics@erba.com, www.erbamannheim.com

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