CONTROL SERUM FOR BIOCHEMISTRY

Assayed Values for ERBA Reagents in ABNORMAL range
SysPack Reagents / Liquid Stable Reagents

Cat. No.	Pack Name	Content
BLT00081	ERBA PATH 4x5	4 x 5 ml
XSYS0124	ERBA PATH 10x5	10 x 5 ml

INTENDED USE

Control serum (assayed) is intended for monitoring the accuracy and precision of commonly used clinical chemistry assays in abnormal range.

PRINCIPLE OF THE PROCEDURE

This human serum control is to be used in the same way as an unknown specimen that would be used in accordance with the direction of the assay procedure.
The results obtained for the control are to be compared with the assigned values given in the "ASSAYED VALUES" section of the insert and an evaluation made by standard statistical techniques to determine if the procedure is within the control limits.

COMPOSITION

The control serum is prepared from human serum with chemical additives and tissue extracts of human and animal origin. Bacteriostatic agents have been added. The control is provided in lyophilized form for increase stability.

RECONSTITUTION

1. Allow the vial to attain room temperature. Carefully open the vial – avoid loss of lyophilisate. Add exactly 5 ml of distilled/deionized water; tightly close the vial and allow to stand for 30 minutes in a light-protected place. Swirl occasionally.
2. Swirl the contents gently to ensure homogeneity before using as sample for testing. Avoid formation of foam.

STORAGE & STABILITY

Prior to reconstitution

The control should be stored at 2–8 °C and is stable till the expiry date printed on the label. Protect from light.

After reconstitution

The constituents when protected from light and contamination are stable:

- 1 day at 18–25 °C
- 1 week at 2–8 °C (Bilirubin: 1 day)
- 1 month at -20 °C

LIMITATIONS

The results obtained using the control are dependent upon several factors. Erroneous results can occur from reconstitution inaccuracy and the technique errors associated with the assay procedure.
Improper storage or handling of the control can also affect the results. If there is a visible evidence of microbial growth in a vial, do not use that vial.

WARNINGS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Control serum contains biological material, which is potentially infectious. It is recommended for this product to be handled with same precautions used for patient specimen.

The blood donations used for production were tested by CE-marked test kits and found to be non-reactive for HBsAg, anti-HIV-1/2 and anti-HCV. In addition HCV and HIV have been tested by PCR. Despite of that the danger of infection for biological material can not be excluded with certainty.

During the work all principles of personal hygiene are to be maintained. Do not eat, drink, smoke. Use the personal protection aids. After finishing the work wash your hands with soap and warm water.

When the vial is broken the contaminated place must be disinfected.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagent is not classified as dangerous.

NOTE

Values for XL analysers are now available on websites:

www.erbalachema.com

- Product Support
→XL-analysers: Calibrator/control values for uploading.

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