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
Intended use
Set for fast staining of blood smears.

Principle
The set contains the solutions for fast staining of blood smears. Result coloration has panoptical character. Staining is performed by dipping smears into colouring solutions. Intensity and tinge of coloration can be changed by number of dips.

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Reagents
R1 Fixing solution
Methanol 24.7 mol/l, Naphtalin green 1 mg/l
R2 Staining solution 1
Eosin Y 1.73 mmol/l, phosphate buffer, pH 6.8, 60 mmol/l
R3 Staining solution 2
Azur II 1.2 g/l, phosphate buffer, pH 6.8, 60 mmol/l
R4 Tablets for rinsing solution
phosphate buffer, pH 7.2, 2.5 mmol/tablet

Storage
Store at (+2 to +25) °C

Preparation of rinsing solution R4
In 250 ml of dist. water dissolve 1 tablet of Reagent 4.
Stability: 2 months at (+2 to +25) °C in the dark.

Procedure
1. On fatless slides make a blood smears in a standard way. Let dry up free on the air.
2. Pour the reagents into convenient staining cuvettes.
3. Fix the smear by dipping 5-times on 1 sec. into Reagent 1. After every dip let the solution flow down and remove the excess by mopping of the drop against the wall of cuvette.
4. Dip the fixed smear 3-times on 1 sec. into Reagent 2. After every dip let the solution flow down and remove the excess by mopping of the drop against the wall of cuvette.
5. Dip the smear 6-times on 1 sec. into Reagent 3. After every dip let the solution flow down and remove the excess by mopping of drop against of wall of cuvette.
6. Rinse the slide by rinsing solution and let dry up, on the air.

Notes
1. The intensity of coloration may be affected by number of dips into staining solutions. The minimal number of dips into each of them is 3.
2. It is possible to increase the intensity of eosinophile coloration by increasing the number of dips into Reagent 2.
3. It is possible to increase the intensity of blue coloration by increasing the number of dips into Reagent 3.
4. All reagents are stable, when protected from direct sunshine. When coloration of standard is different from colouration of smears, change all solutions in staining cuvettes (cca 1-times of week).
5. In case of pathological differential blood picture it is recommended to carry out panoptical coloration according May-Grünwald and Giemsa-Romanowski.
The coloration of the granules of basophiles with LEUKODIF is less remarkable.
This staining is not useful for thrombocytes.

Health protection
For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Danger

Contains: methanol
Hazard statement:
H225 Highly flammable liquid and vapour.
H370 Causes damage to organs.
H301 + H311 + H331 Toxic if swallowed, in contact with skin or if inhaled.
Precautionary statement:
P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P233 Keep container tightly closed.
P260 Do not breathe vapours.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P301 + P310 IF SWALLOWED: Immediately call a POISON CENTER or doctor.

P308+P311 IF exposed or concerned: Call a POISON CENTER or doctor.

UFI: 9VK8-UJVT-8E76-6GPS

First aid


At respiration intoxication transport the patient to fresh air, in contingency start artificial respiration. If swallowed, it is necessary to drink of about 0.5 litre of lukewarm water. On skin or eye contact rinse the place with jet of tap water and exchange the cloths. Attention! If swallowed Reagent 1 drink glass of 40% ethanol and seek medical attention immediately. Mention methanol ingestion! In case of contact of Reagent 1 with eyes, rinse with jet of tap water and transfer casualty to a physician.

Waste disposal

All tested samples should be treated as potentially infectious and with the contingent rest of the reagents should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials.
Put packaging paper waste and rinsed containers to recycling.

Date of revision: 27. 9. 2023

USED SYMBOLS

 Catalogue Number	 Manufacturer	 See Instruction for Use
 Lot Number	 In vitro Diagnostics	 Storage Temperature
 Expiry Date	 Content	

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