## RF CONTROL

Lot: 0000000

Expiry Date: YYYY-MM

Cat. No.	Pack Name	Packaging (Content)
BLT20039	RF CON	1 × 1 ml

## INTENDED USE

Accuracy control for the determination of Rheumatoid Factor (RF) in human serum by turbidimetry and nephelometry

## COMPOSITION

A dilution of human plasma containing a high level of RF with saline. The dilution is liquid stabilised.

Ready to use.

## STORAGE AND STABILITY

After first opening the container, the serum can be used for 6 weeks if stored tightly closed at 2–8 °C after use.

Do not freeze.

## WARNING AND PRECAUTIONS

- For *in Vitro* Diagnostic use.
- Each individual donation intended for use in manufacture of protein control serum was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C virus (anti-HCV) and anti-HIV-1 and HIV-2 by FDA required test. Only donations with negative findings were used for its manufacture. Nevertheless every product obtained from human body fluids should be handled with appropriate care in accordance with recommended procedures for biohazardous materials since absence of infectious agents can never be proven.

**Hazards identification in accordance with Regulation (EC) No 1272/2008**

Reagent is not classified as dangerous.

## ASSIGNED VALUES

Titred value for RF is based on WHO standardisation:

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