ALBUMIN

Cat. No.	Pack Name	Packaging (Content)
XSYS0001	ALB 440	R1: 10 x 44 ml

EN

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Albumin in human serum and plasma.

CLINICAL SIGNIFICANCE

Albumin, a major plasma protein, is synthesised in the liver from amino acids which are absorbed from the ileum. It's functions include regulation of distribution of extra-cellular fluid, transportation of various hormones, vitamins and trace metals.

Increased levels are observed in

- Dehydration due to reduced plasma water content.
- Stasis during venipuncture which causes fluid to escape into the extravascular compartment.

Decreased levels are observed in

- Excessive protein loss (mainly albumin) - from kidney, skin or intestine.
- Decreased synthesis due to dietary, hepatic disease or malabsorption.
- Increased catabolism in fever, untreated diabetes mellitus and hypertension.

PRINCIPLE

Albumin binds with Bromo Cresol Green (BCG) at pH 4.2 causing a shift in absorbance of the yellow BCG dye. The blue-green colour formed is proportional to the concentration of albumin, when measured photometrically between 540–630 nm with maximum absorbance at 625 nm.

REAGENT COMPOSITION

R1

Bromocresol green	0.21 mmol/l
Succinate Buffer	100 mmol/l
Sodium Azide	0.5 g/l

REAGENT PREPARATION

Reagent is liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 30 days if refrigerated (2–10 °C) and not contaminated.

SPECIMEN COLLECTION AND HANDLING

Use unheamolytic serum or plasma (heparin, EDTA).

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability in serum: 1 month at 2–8 °C
1 week at 15–25 °C
at least 3 months at -20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

Traceability:

This calibrator has been standardized with using the ERM-DA470.

QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

CALCULATION

The XL Results are calculated automatically by the instrument.

UNIT CONVERSION

g/dl x 10 = g/l

EXPECTED VALUES ²

Serum:

0–4 d	2.8–4.4 g/dl
4 d – 14 y	3.8–5.4 g/dl
14–18 y	3.2–4.5 g/dl
Adult	
20–60 y	3.5–5.2 g/dl
60–90 y	3.2–4.6 g/dl

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification:	0.1 g/dl
Linearity:	7.2 g/dl
Measuring range:	0.1–7.2 g/dl

Intra-assay precision Within run (n=20)	Mean (g/dl)	SD (g/dl)	CV (%)
Sample 1	3.303	0.017	0.51
Sample 2	4.926	0.044	0.88

Inter-assay precision Run to run (n=20)	Mean (g/dl)	SD (g/dl)	CV (%)
Sample 1	3.345	0.043	1.29
Sample 2	4.612	0.043	0.93

COMPARISON

A comparison between XL-Systems Albumin (y) and a commercially available test (x) using 40 samples gave following results:

y = 1.000 x + 0.065 g/dl

r = 0.997

INTERFERENCES

Following substances do not interfere:

haemoglobin up to 10 g/l, bilirubin up to 40 mg/dl, triglycerides up to 2000 mg/dl.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagent 1 contains < 1.2 % succinic acid.

UFI: Y8RN-RJ05-OE7S-53CM

Warning

Hazard statement:

H319 Causes serious eye irritation.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313 If eye irritation persists: Get medical advice/attention.

WASTE MANAGEMENT

Please refer to local legal requirements.

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