ALPHA AMYLASE (SINGLE REAGENT)

Cat. No.	Pack Name	Packaging (Content)
XSYS0003	AMY 110	R1: 5 x 22 ml

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of alpha-Amylase in human serum, plasma and urine.


CLINICAL SIGNIFICANCE

α -Amylase is derived mainly from the salivary glands and the exocrine pancreas. α -Amylase catalyses the hydrolysis of α -1-4 glucosidic linkages of starch and other related polysaccharides to produce maltose and other oligosaccharides. The enzyme is a relatively small molecule which is rapidly cleared by the kidneys and excreted in the urine.

α -Amylase is most frequently measured in the diagnosis of acute pancreatitis when serum levels may be grossly elevated. In acute pancreatitis α -amylase starts to rise approximately 4 hours after the onset of pain, reaches a peak at 24 hours and remains elevated for 3–7 days. Hyperamylasemia is also associated with other acute abdominal disorders, biliary dysfunction, salivary gland disorders, ruptured ectopic pregnancy and macroamylasemia.

PRINCIPLE

2-Chloro-4-nitrophenol- β -1-4 galactopyranosylmaltotrioxide (CNP-G) is a direct substrate for determination of α -amylase activity, which does not require the presence of ancillary enzymes. The rate of 2-chloro-4-nitrophenol formation can be monitored at (400–420) nm and is proportional to the α -amylase activity.

REAGENT COMPOSITION

MES buffer	50 mmol/l
Calcium Chloride	3.81 mmol/l
Sodium Chloride	300 mmol/l
Potassium Thiocyanate	450 mmol/l
Sodium Azide	13.85 mmol/l
CNPG	0.91 mmol/l

REAGENT PREPARATION

Reagent is liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 30 days if refrigerated (2–10 °C) and not contaminated.

SPECIMEN COLLECTION & HANDLING

Use serum, plasma (heparin, EDTA) or urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability

in serum/plasma:	7 days	at 20–25 °C
	7 days	at 4–8 °C
	1 year	at -20 °C

in urine:	2 days	at 20–25 °C
	10 days	at 4–8 °C
	3 weeks	at -20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL 4x3, Cat. No. XSYS0034 or XL MULTICAL 10x3, Cat. No. XSYS0122 is recommended.

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

Traceability:

This calibrator has been standardized to IFCC formulation.

QUALITY CONTROL

For quality control ERBA NORM 4x5, Cat. No. BLT00080 or ERBA NORM 10x5, Cat. No. XSYS0123 and ERBA PATH 4x5, Cat. No. BLT00081 or ERBA PATH 10x5, Cat. No. XSYS0124 are recommended.

CALCULATION

Results are calculated automatically by the instrument.

UNIT CONVERSION

U/l x 0.017 = μ kat/l

EXPECTED VALUES ⁷

at 37 °C

Serum: up to 80 U/l

Urine: up to 500 U/l

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification:	10.8 U/l
Linearity:	1500 U/l
Measuring range:	10.8–1500 U/l

Intra-assay precision Within run (n=20)	Mean (U/l)	SD (U/l)	CV (%)
Sample 1	247.1	2.5	1.0
Sample 2	260.8	2.6	1.0

Inter-assay precision Run to run (n=20)	Mean (U/l)	SD (U/l)	CV (%)
Sample 1	58.3	1.3	2.2
Sample 2	142.9	2.3	1.6

COMPARISON

A comparison between XL-Systems Amylase (y) and a commercially available test (x) using 40 samples gave following results:

$$y = 0.973 x - 4.80 \text{ U/l}$$

$$r = 0.989$$

INTERFERENCES

Following substances do not interfere:

haemoglobin up to 2.5 g/l, bilirubin up to 40 mg/dl, triglycerides up to 2000 mg/dl.

Note:

Saliva and skin contain alpha-amylase therefore never pipette reagents by mouth and avoid contamination of samples and reagents. Even trace contamination can affect results.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Hazards identification in accordance with Regulation (EC) No 1272/2008

UFI: 6RYE-AE3K-3H7A-Q7QR

Danger

Contains: potassium thiocyanate

Hazard statement:

H318 Causes serious eye damage.

Precautionary statement:

P280 Wear eye protection.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 If eye irritation persists: Get medical advice/attention.

WASTE MANAGEMENT

Please refer to local legal requirements.

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ASSAY PARAMETERS (conventional units)

Instrument	XL-100 EM-100	XL-200 EM-200	XL-300/600 EM-360	XL-640	XL-1000	XL-180
Test Details						
Test	AMY	AMY	AMY	AMY	AMY	AMY
Test Code	6	6	6	6	6	6
Report Name	Amylase	Amylase	Amylase	Amylase	Amylase	Amylase
Unit	U/l	U/l	U/l	U/l	U/l	U/l
Decimal Places	0	0	0	0	0	0
Wavelength-Primary	405	405	415	415	405	405
Wavelength-Secondary	700	700	700	700	700	700
Assay type	Rate-A	Rate-A	Rate-A	Rate-A	Rate-A	Rate-A
Curve type	Linear	Linear	Linear	Linear	Linear	Linear
M1 Start	0	0	0	0	0	0
M1 End	0	0	0	0	0	0
M2 Start	6	6	10	11	6	6
M2 End	12	12	19	21	12	12
Sample replicates	1	1	1	1	1	1
Standard replicates	3	3	3	3	3	3
Control replicates	1	1	1	1	1	1
Control interval	0	0	0	0	0	0
Reaction Direction	Increasing	Increasing	Increasing	Increasing	Increasing	Increasing
React. Abs. Limit	2.1	2.1	2	2	2	2.1
Prozone Limit %	0	0	0	0	0	0
Prozone Check	Lower	Lower	Lower	Lower	Lower	Lower
Linearity Limit %	0	0	0	0	0	0
Delta Abs/Min	0	0	0	0	0	0
Technical Minimum	10.8	10.8	10.8	10.8	10.8	10.8
Technical Maximum	1500	1500	1500	1500	1500	1500
Y=aX+b						
a=	1	1	1	1	1	1
b=	0	0	0	0	0	0
Reagent Abs Min	NA	NA	NA	NA	NA	NA
Reagent Abs Max	0.14	0.14	0.14	0.14	0.14	0.14
Auto Rerun	No	No	No	No	No	No
Total Reagents	1	1	1	1	1	1
Reagent R1	AMY R1	AMY R1	AMY R1	AMY R1	AMY R1	AMY R1
Reagent R2	NA	NA	NA	NA	NA	NA
Reagent R3	NA	NA	NA	NA	NA	NA
Test Volumes						
Test	AMY	AMY	AMY	AMY	AMY	AMY
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM	SERUM
Sample Volumes						
Normal	4	4	4	4	3	4
Dilution Ratio	1	1	1	1	1	1
Increase	8	8	8	8	6	8
Dilution Ratio	1	1	1	1	1	1
Decrease	2	2	2	2	2	2
Dilution Ratio	1	1	1	1	1	1
Standard volume	4	4	4	4	3	4
Reagent Volumes and Stirrer speed						
RGT-1 Volume	200	200	200	200	150	200
R1 Stirrer Speed	High	High	NA	Medium	Medium	High
RGT-2 Volume	0	0	0	0	0	0
R2 Stirrer Speed	NA	NA	NA	NA	NA	NA
RGT-3 Volume	0	0	0	0	0	0
R3 Stirrer Speed	NA	NA	NA	NA	NA	NA
Reference Ranges						
Test	AMY	AMY	AMY	AMY	AMY	AMY
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM	SERUM
Reference Range	Default	Default	Default	Default	Default	Default
Category Male						
Normal-Lower Limit	0	0	0	0	0	0
Normal-Upper Limit	80	80	80	80	80	80
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA
Category Female						
Normal-Lower Limit	0	0	0	0	0	0
Normal-Upper Limit	80	80	80	80	80	80
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA

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