# CHLORIDE

| Cat. No. | Pack Name | Packaging (Content) |
|----------|-----------|---------------------|
| XSYS0008 | CL 120    | R1: 10 x 12 ml      |

## INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Chloride in human serum, plasma and urine.

## CLINICAL SIGNIFICANCE

Chloride is estimated in conjunction with the other electrolytes and used to support interpretation of the other electrolytes. The anion gap ( $\text{Na} + \text{K} - (\text{Cl} + \text{HCO}_3)$ ) between these electrolytes is frequently affected in disease states. Increased levels can occur in nephritis, prostatic obstruction, eclampsia and dehydration. Decreased levels can occur with impaired gastrointestinal or renal function.

## PRINCIPLE

When chloride is mixed with a solution of undissociated mercuric thiocyanate, the chloride preferentially combines with mercury forming mercuric chloride. The thiocyanate that is released then combines with ferric ions present in the solution forming strongly coloured ferric thiocyanate with an absorption maxima at 480 nm.

## REAGENT COMPOSITION

|                      |              |
|----------------------|--------------|
| <b>R1</b>            |              |
| Mercuric Thiocyanate | 1.01 mmol/l  |
| Ferric Nitrate       | 37.13 mmol/l |
| Nitric Acid          | 62.43 mmol/l |
| Mercuric Nitrate     | 105 µmol/l   |

## REAGENT PREPARATION

Reagent R1 is liquid, ready to use.

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.  
On board stability: min. 30 days if refrigerated (2–10 °C) and not contaminated.

## SPECIMEN COLLECTION & HANDLING

Use serum or plasma (heparin) or urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

**Stability:** 7 days at 20–25 °C  
7 days at 4–8 °C  
at least one year at -20 °C

Discard contaminated specimens.

## CALIBRATION

Calibration with calibrator XL MULTICAL 4x3, Cat. No. XSYS0034 or XL MULTICAL 10x3, Cat. No. XSYS0122 is recommended.

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

## Traceability:

This calibrator has been standardized to coulometry.

## QUALITY CONTROL

For quality control ERBA NORM 4x5, Cat. No. BLT00080 or ERBA NORM 10x5, Cat. No. XSYS0123 and ERBA PATH 4x5, Cat. No. BLT00081 or ERBA PATH 10x5, Cat. No. XSYS0124 are recommended.

## CALCULATION

Results are calculated automatically by the instrument.

## UNIT CONVERSION

mmol/l = mEq/l

## EXPECTED VALUES <sup>1</sup>

### Serum:

|           |               |
|-----------|---------------|
| Adult     | 98–107 mmol/l |
| Cord      | 96–104 mmol/l |
| Premature | 95–110 mmol/l |
| 0–30 d    | 98–113 mmol/l |

**It is recommended that each laboratory verify this range or derives reference interval for the population it serves.**

## PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

|                                 |                |
|---------------------------------|----------------|
| <b>Limit of quantification:</b> | 3.7 mmol/l     |
| <b>Linearity:</b>               | 160 mmol/l     |
| <b>Measuring range:</b>         | 3.7–160 mmol/l |

| Intra-assay precision<br>Within run (n=20) | Mean<br>(mmol/l) | SD<br>(mmol/l) | CV<br>(%) |
|--|------------------|----------------|-----------|
| Sample 1                                   | 116.1            | 1.25           | 1.08      |
| Sample 2                                   | 108.3            | 1.64           | 1.52      |

| Inter-assay precision<br>Run to run (n=20) | Mean<br>(mmol/l) | SD<br>(mmol/l) | CV<br>(%) |
|--|------------------|----------------|-----------|
| Sample 1                                   | 117              | 1.47           | 2.75      |
| Sample 2                                   | 94.6             | 1.43           | 1.52      |

## COMPARISON

A comparison between XL-Systems Chloride (y) and a commercially available test (x) using 40 samples gave following results:

$$y = 1.000 x + 1.0 \text{ mmol/l}$$

$$r = 0.927$$

## INTERFERENCES

Following substances do not interfere:

haemoglobin up to 10 g/l, bilirubin up to 40 mg/dl, triglycerides up to 500 mg/dl.

## Note:

Grossly haemolysed, lipaemic and jaundiced samples give falsely elevated results and should not be used. High levels of immunoglobulins, as found in case of multiple myeloma will interfere due to the development of turbidity.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

**Hazards identification in accordance with Regulation (EC) No 1272/2008**

UFI: 058F-WERG-3H7T-HPU2

Contains: methanol, Mercury dithiocyanate, Mercury dinitrate.

**Danger**

## Hazard statement:

H302 Harmful if swallowed.  
H314 Causes severe skin burns and eye damage.  
H371 May cause damage to organs.

## Precautionary statement:

P260 Do not breathe vapours.  
P280 Wear protective gloves/protective clothing/eye protection.  
P301+P330+P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.  
P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.  
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P308+P313 IF exposed or concerned: Get medical advice/attention.

## WASTE MANAGEMENT

Please refer to local legal requirements.

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