ALT/GPT

Cat. No.	Pack Name	Packaging (Content)
XSYS0017	ALT/GPT 330	R1: 6 x 44 ml, R2: 6 x 11 ml
XSYS0074	ALT/GPT 564 XL-1000	R1: 6 x 72 ml, R2: 6 x 22 ml

EN

INTENDED USE


Diagnostic reagent for quantitative *in vitro* determination of ALT/GPT (Alanine Amino transferase) in human serum and plasma.

CLINICAL SIGNIFICANCE

ALT/GPT is present in high concentration in liver and to a lesser extent in kidney, heart, skeletal muscle, pancreas, spleen and lung. Increased levels of ALT/GPT however is generally a result of liver disease associated with some degree of hepatic necrosis such as cirrhosis, viral or toxic hepatitis and obstructive jaundice. Characteristically ALT/GPT is generally higher than AST/GPT in acute viral or toxic hepatitis, whereas for most patients with chronic hepatic disease, ALT/GPT levels are generally lower than AST/GPT levels. Elevated ALT/GPT levels have also been found in extensive trauma and muscle disease, circulatory failure with shock, hypoxia, myocardial infarction and haemolytic disease.

PRINCIPLE

This ALT/GPT reagent is based on the recommendations of the IFCC without pyridoxal phosphate. ALT/GPT catalyzes the reaction between L-alanine and 2-oxoglutarate. The pyruvate formed is reduced by NADH in a reaction catalyzed by lactate dehydrogenase (LDH) to form L-lactate and NAD⁺.

The reaction is monitored by measuring the rate of decrease in absorbance at 340 nm due the oxidation of NADH. The rate of oxidation of NADH is proportional to the activity of ALT/GPT present in the sample. Endogenous sample pyruvate is reduced by LDH during initial incubation period.

REAGENT COMPOSITION

R1

Tris buffer (pH 7.5)	137.5 mmol/l
L-Alanine	709 mmol/l
LDH (microbial)	≥ 2000 U/l

R2

CAPSO	20 mmol/l
2-oxoglutarate	85 mmol/l
NADH	1.05 mmol/l

REAGENT PREPARATION

Reagents are liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 30 days if refrigerated (2–10 °C) and not contaminated.

SPECIMEN COLLECTION AND HANDLING

Use unheamolytic serum or plasma (heparin, EDTA).

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Loss of activity:

within 3 days at 2–8 °C < 10 %

within 3 days at 15–25 °C < 17 %

Stability at least 3 months at -20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL 4x3 (Cat. No. XSYS0034) or XL MULTI-CAL 10x3 (Cat. No. XSYS0122) is recommended.

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

Traceability:

This calibrator has been standardized to IFCC formulation (without PDP).

QUALITY CONTROL

For quality control ERBA NORM 4x5 (Cat. No. BLT00080) or ERBA NORM 10x5 (Cat. No. XSYS0123) and ERBA PATH 4x5 (Cat. No. BLT00081) or ERBA PATH 10x5 (Cat. No. XSYS0124) are recommended.

CALCULATION

The XL Results are calculated automatically by the instrument.

UNIT CONVERSION

U/l x 0.017 = µkat/l

EXPECTED VALUES ⁴

At 37 °C

Men up to 45 U/l

Women up to 34 U/l

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification:	4.4 U/l
Linearity:	360 U/l
Measuring range:	4.4–360 U/l

PRECISION

Intra-assay precision Within run (n=20)	Mean (U/l)	SD (U/l)	CV (%)
Sample 1	58.32	2.52	4.31
Sample 2	114.72	1.38	1.20

Inter-assay precision Run to run (n=20)	Mean (U/l)	SD (U/l)	CV (%)
Sample 1	40.20	1.20	3.12
Sample 2	112.8	2.40	1.95

COMPARISON

A comparison between XL-Systems ALT/GPT (y) and a commercially available test (x) using 40 samples gave following results:

y = 0.979 x - 1.80 U/l

r = 0.996

INTERFERENCES

Following substances do not interfere:

haemoglobin up to 2.5 g/l, bilirubin up to 30 mg/dl, triglycerides up to 2000 mg/dl.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagent is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erba.com, www.erbamannheim.com

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1. Thomas L. Alanine aminotransferase (ALT), Aspartate aminotransferase (AST). In: Thomas L, editor. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 55-65.
2. Moss DW, Henderson AR. Clinical enzymology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 617-721.
3. Schumann G, Bonora R, Ceriotti F, Férard G et al. IFCC primary reference procedure for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 5: Reference procedure for the measurement of catalytic concentration of alanine aminotransferase. Clin Chem Lab Med 2002;40:718-24.
4. Tietz Textbook of Clinical Chemistry. Burtis CA and Ashwood ER, Fifth Edition, 2012.

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