# ANTISTREPTOLYSIN (O)

Cat. No.	Pack Name	Packaging (Content)
XSYS0046	ASO	R1: 2 x 40 ml (Buffer), R2: 2 x 10 ml (Latex)

## INTENDED USE

Quantitative determination of Antistreptolysin (O) in human serum by turbidimetric immunoassay.

## CLINICAL SIGNIFICANCE

The group A  $\beta$ -haemolytic streptococci produces various toxins that can act as antigens, one of these exotoxins is streptolysin O. The affected organism produces specific antibodies against streptolysin O. The concentration of ASO in the patient's serum will enable to establish the degree of infection due to  $\beta$ -haemolytic streptococci.

## METHODOLOGY

Measurement of antigen-antibody reaction by the end-point method.

## REAGENT COMPOSITION

### R1 (Buffer)

Phosphate buffered saline (pH 7.43)

Polyethylene glycol 40 g/l

Sodium azide (< 0.1 %)

### R2 (Latex)

Glycine Buffer (pH 8.2)

ASO sensitized Latex 0.17 %

Sodium azide (< 0.1 %)

## REAGENT PREPARATION

Liquid reagents are ready to use.

## STABILITY AND STORAGE

The reagents are stable until expiry date when kept at 2–8 °C. Stability in the instrument is at least 4 weeks if contamination is avoided. Do not freeze.

## SAMPLE COLLECTION

Use fresh serum. If the test can not be carried out on the same day, the serum may be stored at 2–8 °C for 48 hours. If stored for a longer period, the sample should be frozen.

## MATERIALS REQUIRED BUT NOT PROVIDED

- Any instrument with temperature control of  $37 \pm 0.5$  °C that is capable of reading absorbance accurately at 600 nm may be used.
- Analyser specific consumables such as sample cups.
- Controls.
- Saline (9 g/l NaCl)

## ASSAY PROCEDURE

Refer to the assay parameters for details.

## CALIBRATION

Blank: Saline

Cat. No.:	Product name	Pack name	Content
XSYS0051	ASO CALIBRATOR	ASO CAL SH	1 x 1 ml

Calibration curve: generate a 5 point calibration curve by diluting the calibrator 1:16, 1:8, 1:4, 1:2 and undiluted in saline.

## Calibration frequency

Calibration verification: Not necessary.

## Quality control

For quality control use

Cat. No.:	Product name	Pack name	Content
BLT20004	ASO CONTROL	ASO CON	1 x 1 ml
BLT20034	MULTICONTROL LEVEL 1	MULTICON L1	1 x 1 ml
BLT20035	MULTICONTROL LEVEL 2	MULTICON L2	1 x 1 ml

## CALCULATION

Results are calculated automatically by the instrument.

## EXPECTED VALUES

### Reference Values

Normal values: 0–200 IU/ml (WHO)

**This range is given for orientation only. Each Laboratory should establish its own reference values.**

## PERFORMANCE DATA

Data contained within this section are representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

## COMPARISON

A comparison between XL-Systems ASO (y) and a commercially available test (x) using 40 samples gave following results:

$y = 0.930 x - 5.80$  IU/ml

$r = 0.997$

Measuring Range: 0–480 IU/ml

Detection Limit: 7.2 IU/ml

Hook Effect: No

## PRECISION

Intra-assay precision Within run (n =20)	Mean (IU/ml)	SD (IU/ml)	CV (%)
Sample 1	256	3.54	1.38
Sample 2	329	3.07	0.94

Inter-assay precision Run to run (n =20)	Mean (IU/ml)	SD (IU/ml)	CV (%)
Sample 1	172	5.53	3.22
Sample 2	307	11.6	3.77

## Specificity:

Monospecific

## Interferences:

No interference for: Hemoglobin (1 g/dl), Bilirubin (40 mg/dl), Triacylglycerols (1 g/dl).

## Limitations:

None

## WARNING AND PRECAUTIONS

1. For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.
2. Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion. Flush drains with water thoroughly after disposing of fluids containing sodium azide.
3. Polyethylene glycol is non biohazardous.
4. Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV1 and HIV2 antibodies, as well as for the hepatitis B surface antigen and anti-hepatitis C antibodies, using a method approved by the FDA.

## Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagents R1 and R2 are not classified as dangerous.

## WASTE MANAGEMENT

Please refer to local legal requirements.

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