# HbA1c CONTROL

Lot: XXXXXXXX; Expiry Date: mm/yyyy

Cat. No.	Pack Name	Packaging (Content)
XSYS0056	HBA1C CON H	1 x 0.25 ml

EN

## INTENDED USE

HbA1c liquid control is a human-based control. The HbA1c concentration in HbA1c Control High is pathological.

## STORAGE AND STABILITY

The controls both unopened and reconstituted must be stored at 2–8 °C, protected from light and heat.

Unopened and opened: until the end of the indicated month of expiry if contamination and evaporation are avoided after having opened the bottles.

Proper storage and handling of this product must be observed.

## WARNINGS AND PRECAUTIONS

Reagent of the kit is not classified like dangerous but contains less than 0.1 % sodium azide classified as very toxic and dangerous substance for the environment.

Each individual blood donation used for production of HbA1c control was found to be non-reactive when tested with approved methods for HBsAg, anti-HIV-1/2 and anti-HCV. As there is no possibility to exclude definitely that products derived from human blood transmit infectious agents, it is recommended to handle the control with the same precautions used for patient specimens.

## PREPARATION

HbA1c control is ready to use. Controls must be treated the same way as patient samples. Please refer to the package insert of the reagent HbA1c, Cat. No. XSYS0054

### Control preparation:

Hemolyzing Solution (R3) 500 µl  
Control 10 µl

Mix and allow to stand for 5 minutes or until complete lysis is apparent.

## ASSIGNED VALUES

HbA1c Control High	Targer value (mmol/mol)	Range (mmol/mol)	1 SD (mmol/mol)
Target values according to IFCC (mmol/mol):	xx.x	xx.x – xxx	x.xx

HbA1c Control High	Targer value (%)	Range (%)	1 SD (%)
Target values according to DCCT/NGSP (%):	x.xx	x.xx – xx.x	x.xx

## PROCEDURE

Please refer to the reagent package insert for instructions for use.

## ASSAY VALUES

HbA1c control values according to DCCT/NGSP in % and according to IFCC in mmol/mol have been derived from percentage values according to IFCC by calculation.

## WASTE MANAGEMENT

Please refer to local legal requirements.

## LITERATURE

- The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes in the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med.1993; 329:977-86.
- Little RR, Rohlfing CI, Wiedmeyer HM, Myers GL et al. The national Glycohemoglobin Standardization Program: A five-Year Progress Report. Clin Chem 2001; 47:1985-92.
- Jeppsson JO, Kobold U, Barr J, Finke A et al. Approved IFCC Reference method for the measurement of HbA1c in human blood. Clin Chem Lab Med 2002; 40:78-89.
- Hoelzel W, Weykamp C et al. IFCC Reference system for measurement of HbA1c in Human Blood and the national standardization Schemes in the United States, Japan and Sweden: A Method Comparison Study. Clin Chem 2004; 50:1:1666-74.
- Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 1393-1401.
- Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).

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- The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes in the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med.1993; 329:977-86.
- Little RR, Rohlfing CI, Wiedmeyer HM, Myers GL et al. The national Glycohemoglobin Standardization Program: A five-Year Progress Report. Clin Chem 2001; 47:1985-92.
- Jeppsson JO, Kobold U, Barr J, Finke A et al. Approved IFCC Reference method for the measurement of HbA1c in human blood. Clin Chem Lab Med 2002; 40:78-89.
- Hoelzel W, Weykamp C et al. IFCC Reference system for measurement of HbA1c in Human Blood and the national standardization Schemes in the United States, Japan and Sweden: A Method Comparison Study. Clin Chem 2004; 50:1:1666-74.
- Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 1393- 1401.
- Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).

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1. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes in the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med. 1993; 329:977-86.
2. Little RR, Rohlfing CI, Wiedmeyer HM, Myers GL et al. The national Glycohemoglobin Standardization Program: A five-Year Progress Report. Clin Chem 2001; 47:1985-92.
3. Jeppsson JO, Kobold U, Barr J, Finke A et al. Approved IFCC Reference method for the measurement of HbA1c in human blood. Clin Chem Lab Med 2002; 40:78-89.
4. Hoelzel W, Weykamp C et al. IFCC Reference system for measurement of HbA1c in Human Blood and the national standardization Schemes in the United States, Japan and Sweden: A Method Comparison Study. Clin Chem 2004; 50:1:1666-74.
5. Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 1393-1401.
6. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).

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