## XL WASH

Cat. No.	Pack Name	Packaging (Content)
XSYS0066	XL WASH	4 × 100 mL, Instruction for use

### INTENDED USE

XL WASH is washing solution contains a surfactant which solubilises proteins and cleans the microcuvettes of analysers ERBA XL intended for *in vitro* diagnostics. For professional use in clinical laboratory only.

### PRINCIPLE

The performance of any Erba XL analyzer depends on the cleanliness of the microcuvettes. The use of serum or plasma or urine samples results in protein deposition. Over a period of time this may interfere with the assay resulting in erroneous results. XL WASH contains a surfactant which solubilises proteins and cleans the microcuvettes.

### REAGENT DESCRIPTION AND COMPOSITION

Sodium hydroxide	0.5 – <1 %
Polyethylene glycol dodecyl ether with >5–20 mol ethylene oxide	0.25 – <1 %
Emulsifier K30	0.25 – <1 %

### REAGENT PREPARATION

The contents of each bottle (100 mL) of XL WASH is to be diluted to 10 L with deionized water (100 mL XL WASH + 9900 mL deionized water) before using for cleaning (Cleaning solution Can) in Erba XL analyzer.

### MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Erba XL – analyzer.

### STABILITY AND STORAGE

The reagent is stable till the expiry date mentioned on the label when stored at room temperature 15–25 °C.

### WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and / or the patient is established.

### Hazards identification in accordance with Regulation (EC) No 1272/2008

UFI: TAHV-CWEQ-CJ5U-90D7

Danger

### Hazard statement:

- H290 May be corrosive to metals.
- H314 Causes severe skin burns and eye damage.
- H412 Harmful to aquatic life with long lasting effects.

### Precautionary statement:

- P280 Wear protective gloves/protective clothing/eye protection/face protection.
- P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.
- P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

### WASTE MANAGEMENT

Please refer to local legal requirements.

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