XL Urine Diluent

Cat. No.	Pack Name	Packaging (Content)
XSYS0067	XL URINE DIL	1 × 52 mL, Instruction for use

EN

INTENDED USE

Urine diluent is a buffered reagent used as a sample diluent for analysing electrolytes in urine samples using ISE module of automatic analysers ERBA XL intended for *in vitro* diagnostics. For professional use in clinical laboratory only.

PRINCIPLE

Urine diluent dilutes urine samples during analyses using ISE module of the Clinical Chemistry Auto-analysers.

REAGENT DESCRIPTION AND COMPOSITION

4-Morpholinepropanesulfonic acid	<0.2%
Polyvinylpyrrolidone	<0.02%
Preservatives	<0.01%

REAGENT PREPARATION

Ready to use.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Erba XL – analyzer with ISE.

STABILITY AND STORAGE

The reagent is stable till the expiry date mentioned on the label when stored at room temperature 15–25 °C.

ASSAY PROCEDURE

Urine samples is to be diluted by a factor of 10 (Urine sample to urine diluent ratio of 1 : 9). Mix the diluent and sample thoroughly before placing the diluted sample on the sample tray. Urine 140 µL (10 times diluted with urine diluent) would be required for the electrolyte measurement.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and / or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

The reagent is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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