ERBA XL AutoWash AC/AL

Cat. No.	Pack Name	Packaging (Content)
XSYS0082	XL AUTOWASH AC/AL	AC 5 × 44 mL, AL: 5 × 44 mL, instruction for use

INTENDED USE

ERBA XL AutoWash AC/AL is a washing solution for cleaning of microcuvettes of automatic analysers ERBA XL intended for *in vitro* diagnostics. For professional use in clinical laboratory only.

PRINCIPLE

The performance of any Auto analyzer depends on the cleanliness of the micro cuvettes. The use of serum or plasma or urine samples results in protein deposition. Over a period of time this may interfere with the assay resulting in erroneous results.

REAGENT DESCRIPTION AND COMPOSITION

AC		AL	
HCl	0.3 %	NaOH	0.2 %
H ₃ PO ₄	<0.75	detergent	
detergent		pH	13.1 ± 0.5
pH	1.20 ± 0.5		

REAGENT PREPARATION

Ready to use.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Erba XL – analyzer.

STABILITY AND STORAGE

The solution is stable till the expiry date mentioned on the label when stored at room temperature 15–30 °C.

ASSAY PROCEDURE

Place the bottles into the reagent tray. Further info about positions is available in User manual for XL analyzers.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and / or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

AC UFI: K3WS-H9Q3-KM7Q-6QFJ

AL UFI: E5WS-19DG-WM77-V21M

Danger

Hazard statement:

H290 May be corrosive to metals.

H314 Causes severe skin burns and eye damage.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

WASTE MANAGEMENT

Please refer to local legal requirements.

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